

# User Guide

## High Sensitivity Quantitative Cartridge Kit (C105105-Q/C105205-Q)

### A. Specifications

| Specification               | Description |
|-----------------------------|-------------|
| DNA Sizing Range            | 20-5,000 bp |
| DNA Sizing Precision*       | 2% CV       |
| DNA Quantitative Precision* | 4% CV       |
| Sample Number               | 100 runs    |
| Shelf Life                  | 4 months    |

\* Precision is determined by the 15-622 DNA Size Marker (C109200).

### B. Kit Components and Storage Conditions

| Item                                                             | Storage Condition                                                |
|------------------------------------------------------------------|------------------------------------------------------------------|
| High Sensitivity Quantitative Cartridge<br>(C105105-Q/C105205-Q) | 4-30°C (Do Not Freeze)                                           |
| 20-1,500 bp Quantitative Marker<br>(C109109-60Q, 60 µL)          | Short-Term (≤ 3 months): 4-30°C<br>Long-Term (> 3 months): -20°C |
| 15-622 bp Size Marker<br>(C109200-50, 50 µL)                     | Short-Term (≤ 3 months): 4-30°C<br>Long-Term (> 3 months): -20°C |
| Separation Buffer<br>(C104406, 50 mL)                            | 4-30°C                                                           |
| Dilution Buffer<br>(C104405, 15 mL)                              | 4-30°C                                                           |
| Mineral Oil<br>(C104404, 8 mL)                                   | 4-30°C                                                           |

⚠ Please always store cartridges in a light-proof bag, and then store in the cartridge box after analysis.

⚠ To maintain marker stability, avoid repeated freeze-thaw cycles.

### C. Cartridge Unpacking Preparation


The new cartridge must undergo **HV check** and **calibration**. Please follow the instructions provided in the unpacking guide and calibrate using C109109 Quantitative Marker.

- Calibration mix for cartridge unpacking:
  - 20 bp-1,500 bp Quantitative Marker (C109109): 5 µL
  - Dilution Buffer (C104405): 15 µL

⚠ **CAUTION:** To prevent cartridge lockout after 10 failed calibration attempts, please contact your local distributor or BiOptic if the issue persists after following the troubleshooting steps in the cartridge unpacking guide.

### D. Buffer, Marker, and Sample Preparation

MA-2: C109200 Size marker mix




*Qsep<sub>100</sub>*


S: Separation Buffer  
P/W/C: d<sub>2</sub>H<sub>2</sub>O  
\*Add to the groove

MB-2: C109301 Size marker mix

⚠ Hold the holder steady and press the alignment marker tube firmly into the well with your thumb.



*Qsep<sub>1</sub>*



*Qsep<sub>1</sub>-Plus*

⚠ Ensure the buffer tray is pushed to the end until the color bar aligns with the edge of the holder.

### 1. Compatible Sample Tubes

|                                                         | Name                       | Cat. No. | Volume | Image |
|---------------------------------------------------------|----------------------------|----------|--------|-------|
| <i>Qsep<sub>100</sub></i><br>/ 8-well<br>Sample<br>Tray | 0.1 mL<br>PCR Tube         | -        | 20 µL  |       |
|                                                         | 0.2 mL<br>PCR Tube         | -        |        |       |
| 12-well<br>Sample<br>Tray                               | 0.1 mL<br>Strip Tube       | C104252  |        |       |
| 16-strip<br>Sample<br>Tray                              | 16-strip<br>Sample<br>Tube | C104254  |        |       |
| 18-strip<br>Sample<br>Tray                              | 18-strip<br>Sample<br>Tube | C104257  |        |       |

### 2. Size Marker Mix for Different Size Ranges

- For Sample Size Range from 20 bp to 1,500 bp:
  - 10X diluted** 20 bp-1,500 bp QM (C109109): 5 µL
  - 10X diluted** 15-622 bp SM (C109200): 10 µL
  - Dilution Buffer (C104405): 5 µL
- For Sample Size Range from 20 bp to 5,000 bp:
  - 10X diluted** 20 bp-5,000 bp QM (C109102): 5 µL
  - 5X diluted** 50 bp-3,000 bp SM (C109301): 10 µL
  - Dilution Buffer (C104405): 5 µL

### 3. Sample Mix Preparation

- 10X diluted** QM (C109109/C109102): 5 µL
- Sample: X (2~15) µL
- Dilution Buffer (C104405): 15-X (0~13) µL

### 4. Recommended Sample Concentration

Fragmented sample: 10-500 pg/µL [Best: 10-200 pg/µL]

Smear sample: 50-1,000 pg/µL [Best: 50-500 pg/µL]

**i** The recommended sample concentration indicates the final concentration for loading onto the instrument.

# User Guide

## High Sensitivity Quantitative Cartridge Kit (C105105-Q/C105205-Q)

### I. Software Operation Guide

Qsep<sub>100</sub>



- ① Sample Position: Place the sample and select the corresponding position. Input sample information if necessary.

- ② Method:

- (a) Analysis Type to **Quantitative** and input the sample volume for the sample mix.
- (b) Select Quantitative Marker.
- (c) Select Analytic Method in the Method Selector.

- Adjust injection conditions based on sample concentration.

| Sample Concentration | High (2kV, 10s) | Regular (4kV, 10s) | Low (8kV, 10s) |
|----------------------|-----------------|--------------------|----------------|
| Fragmented DNA       | > 200 pg/µL     | 10-200 pg/µL       | < 10 pg/µL     |
| Smear DNA            | > 1000 ng/µL    | 50-1000 pg/µL      | < 50 pg/µL     |

- ③ Sample Duration: Adjust the sample injection time to increase/decrease injection amount.

- Modify injection conditions based on sample concentration.

- ④ Runs: Set the repetition time.

- ⑤ Separation Duration: Adjust the duration to extend/reduce the separation time.

(Optional)

- ⑥ Result Name: Input the result name for the result file.

- ⑦ Para: Choose between (a) Reference Marker Table and (b) Create Size Marker for calculation.

- When using "Create Size Marker" function, select the appropriate size marker you use. For example, "20-1.5K" is paired with C109200, "20-5k" is paired with C109301.


8. Click "Run" to start the sequence analysis.

# User Guide

## High Sensitivity Quantitative Cartridge Kit (C105105-Q/C105205-Q)



Qsep<sub>1</sub>




- Click "Open Sample Door" and insert the desired sample tray
- Click "Park" to move the sample tray into position.

| Qsep1                                                                                             | Qsep1-Plus                                                                                                                      |
|---------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------|
| <ul style="list-style-type: none"> <li>8-well Sample Tray</li> <li>12-well Sample Tray</li> </ul> | <ul style="list-style-type: none"> <li>8-well Sample Tray</li> <li>12-well Sample Tray</li> <li>18-strip Sample Tray</li> </ul> |


### ① Sample Position:

- Choose Alignment Marker type (20-1.5K or 20-5K).
- Place the paired Size Marker Mix at S4 position if required, and input the sample information if necessary.



### ② Method:

- Analysis Type to **Quantitative** and input the sample volume for the sample mix.
- Check Quantitative Marker.
- Select Analytic Method in the Method Selector.



- Adjust injection conditions based on sample concentration.

| Sample Concentration | High (2kV, 10s) | Regular (4kV, 10s) | Low (8kV, 10s) |
|----------------------|-----------------|--------------------|----------------|
| Fragmented DNA       | > 200 pg/μL     | 10-200 pg/μL       | < 10 pg/μL     |
| Smear DNA            | > 1000 pg/μL    | 50-1000 pg/μL      | < 50 pg/μL     |

### ③ Sample Duration: Adjust the sample injection time to increase/decrease injection amount.

- Modify injection conditions based on sample concentration.


### ④ Runs: Set the repetition time.

### ⑤ Separation Duration: Adjust the duration to extend/reduce the separation time.

(Optional)

### ⑥ Result Name: Input the result name for the result file.

### ⑦ Para: Choose between (a) Reference Marker Table and (b) Create Size Marker for calculation.



- When using "Create Size Marker" function, select the appropriate size marker you use. For example, "20-1.5K" is paired with C109200, "20-5K" is paired with C109301.

### 8. Click "Run" to start the sequence analysis.


# User Guide

## High Sensitivity Quantitative Cartridge Kit (C105105-Q/C105205-Q)


### F. Result & Application

#### Quantitative Marker & Size Marker Result



20-1.5K & C109200 15-622 bp Size Marker



20-5K & C109301 50-3000 bp Size Marker



#### Cell-Free DNA Validation



### G. Troubleshooting

Before attempting any troubleshooting, ensure that the system is functioning properly and that all operations are following the instructions.

If encountering unstable current during sample injection or separation steps, which may be caused by unknown substances in PCR reagent buffer or other kit buffers, consider the following solutions:

1. Dilute the sample using dilution buffer.
2. Centrifuge the sample for a period to allow residues to accumulate at the bottom of the tube.
3. Insert a "T-purge-120" method between several sample runs. For example, insert one run of "T-Purge-120" every 5-10 sample runs.

| Sequence |                 | Assay       |                 |      |                     |                     |      |        |
|----------|-----------------|-------------|-----------------|------|---------------------|---------------------|------|--------|
| Open     |                 | Save        | Save As         |      |                     |                     |      |        |
| SN       | Sample Position | Method      | Sample Duration | Runs | Separation Duration | Result Name         | Para | Add    |
| 1        | A-01.A...       | cDNA        | 20              | 1    | 200                 | Test 1<br>Sample ID |      | Insert |
| 2        |                 | T-Purge-120 | 0               | 1    | 0                   | Sample ID           |      | Delete |
| 3        | B-01.B...       | cDNA        | 20              | 1    | 200                 | Test 2<br>Sample ID |      | Up     |
|          |                 |             |                 |      |                     |                     |      | Down   |

### H. Cartridge Disposal

Please wear gloves before discarding the cartridge.

#### Gel reservoir



#### Cartridge tip

1. Bend the cartridge tip.
2. Open the cap on the gel reservoir and remove the inner cap.
3. Pour the gel into the chemical waste container.
4. Dispose of the cartridge in the trash bin.

#### Contact Information:

Company Name: BiOptic Inc.

Address: 5F., No.108, Minquan Rd., Xindian District, New Taipei City 23141, Taiwan

Tel: +886-2-2218-8726, Fax: +886-2-2218-8727, E-mail: service@bioptic.com.tw