HC2A-AW & HC2A-AW-USB

Short Instruction Manual

1 GENERAL DESCRIPTION

The Rotronic probes type HC2A-AW and HC2A-AW-USB measures the humidity and the temperature of hygroscopic products. The digital probe signal is processed by an indicator such as the HygroLab or from a PC (HC2A-AW-USB). The desired scale, 0...100 %rh equilibrium humidity or 0...1 Aw water activity, can be selected through the indicator.

For more information, please visit www.processsensing.com or www.rotronic.com or scan the QR code (also on the measurement head).

2 GENERAL GUIDELINES

To make accurate and reproducible measurement, follow the following general guidelines to measure water activity:

- Measure water activity only in rooms with constant temperature; ideal are rooms with small daily and yearly temperature variations. Avoid heating of the measuring equipment by sun radiation or radiators, draught, body radiation etc.
- Put the samples into the disposable sample cup and put the lid on it. The probe and the sample must stay in the same room for acclimation (normally room temperature). An error, which occurs very often in the field is putting the sample a night before the measurement will be done into the fridge without acclimation. Another error is to measure a sample which comes directly from production. Use a climate chamber to measure water activity on another temperature, then the ambient temperature. Set the desired temperature and put the probe and sample into it.
- Avoiding warming up the probe and sample through body heat. Do not hold it too long in your hands.
- Check the proper connection between the probe and the indicator. The green LED light ring illuminates when both, the power and the communication are correct.


3 MEASUREMENT ARRANGEMENT

For each measurement, a sample cup is necessary.

Exception: For the measurement of paper, cardboard and other similar products simply put the probe on the sample.

The whole measuring system, probe + sample cup + sample, adapts to the ambient temperature.

4 STARTING UP THE HC2A-AW


Connect the probe cable with the HygroLab. The HygroLab powers the HC2A-AW probe, press the power button to switch the probe on. The LED light ring starts to light up in green as soon as the HC2A-AW is powered up and has made connection with the HygroLab.

5 STARTING UP THE HC2A-AW-USB

The HygroSoft water activity version with AW-Quick function must be installed properly. For detailed information regarding the installation and handling of the HygroSoft please visit the Rotronic Online Manual. Connect the HC2A-AW-USB via USB interface to the PC. The PC powers the HC2A-AW-USB probe, press the power button to switch the probe on. The LED light ring starts to light up as soon as the HC2A-AW-USB is powered up and has made connection with the PC.

6 GENERAL GUIDELINES

Put the sample into the sample cup. To avoid cross-contamination do not fill more than up to the fill line.

Put the probe immediately on the sample cup after having filled in the sample to avoid humidity exchange with the ambient. The measurement time may vary depends on the probe. Humidity and temperature value are ready to read when the values are stable for at least 1–2 minutes.


7 PERIODIC CALIBRATIONS

Both the Pt100 RTD temperature sensor and the corresponding electronics are very stable and do not normally need to be changed or calibrated after factory calibration. The long-term stability of the Rotronic HYGROMER humidity probes is typically better than 1 %rh per year. For maximum accuracy we recommend calibration of the probe about every 6 to 12 months. More frequent calibration can be necessary in applications where the sensor is exposed to pollutants. The calibration can be performed by the user himself on site or in the laboratory / workshop. For routine calibrations the probe should be checked at least at two points, better at three points. The electronics of the HC2A-AW-(USB) do not require calibration in the field. The electronics cannot be repaired in the field and should be returned to the manufacturer in the case of problems.

8 CONTAMINATIONS

Contaminated filters may cause measuring errors and prolong the measuring time. Clean the filter periodically or, if necessary, replace it depending upon the degree of contamination. In order not to damage the sensors, remove the filter for cleaning. Clean the filter with soap water, alcohol or a cleaning agent suitable for removing the contamination and rinse it thoroughly at the end with water. Do not install the filter until it is completely dry. Replace filters which can no longer be cleaned by new ones. Should the sensors be severely contaminated despite of the protecting filter, then we recommend having these replaced by our service department.

9 DIMENSIONS

10 TECHNICAL DATA

	HC2A-AW	HC2A-AW-USB
Sensor	HYGROMER WA-1	
Accuracy @ 23 ±5 °C	± 0.008 aw / 0.8 %rh / ±0.1 K	
Operating conditions	-40...85°C / 0...100%rh	
Storage conditions	-40...85°C / 0...100%rh	
Long-term stability	< 1 %rh / year (normal conditions)	
Supply voltage	3.3...5 V	
LED light ring	Yes, only green light	
AW Quick Mode	Depending on the indicator	
Display	No	
Button	ON / OFF touch button	
Filter	Stainless steel wire mesh	
Cable length	1m	3m
Connector	HC2	USB-A
HygroSoft compatibility	Yes (with AC3001)	Yes
HygroLab compatibility	Yes	No
Dimensions	Ø68 x 65 mm	
Housing material	Aluminum, stainless steel & PC	
Weight	540 g	
CE conform	CE 2014/30/EU	
IP protection rating	IP65	

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