

FlashGel® System with FlashGel® FLX Cassettes Quick Start Guide

For FlashGel® FLX Cassettes (Cat. No. 317258, 317259, 317260, 317281)

Guidelines and instructions for use

For complete safety information, please refer to the FlashGel® System Manual – Your guide to setup and usage of the FlashGel® System.

Important points

- Do not exceed 5 µL sample volume per lane for 12 + 1 well and 16 + 1 well cassettes.
- Optimal sample concentrations are approximately 1/5 the typical per-band concentration of an ethidium bromide gel.
- Flood sample wells with water prior to loading, then load samples into wells filled with water.
- FlashGel® FLX Cassettes must use the appropriate **DNA or RNA Prestain Loading Buffer** for sample detection.
- Use the FlashGel® Mask when running double-tier cassettes.
- For best results with restriction enzyme digests: optimal loading levels are in the range of 20 – 40 ng of DNA. Minimize salt concentrations by using <1.5 µL of digest in the total sample (5 µL) for the FlashGel® FLX Cassette.

Instructions

1. Refer to Table 1 for recommended sample preparation and run conditions.
2. Remove white well seals and tear-off vent seals from cassette. Do not remove the black tips of vent seals remaining on the cassette.
3. Flood sample wells with distilled or deionized water. Tilt cassette to move excess fluid to the edge and blot off with a lint-free wipe. Do not blot wells directly. Leave wells filled with water for sample loading.
4. Insert cassette into dock. Insert FlashGel® Mask under the central tier of sample wells if using double-tier cassettes.
5. Load samples.
6. Plug in high voltage cables, turn on power supply and set to recommended voltage.
7. Plug in low voltage power supply and turn on viewing light as needed during run.
8. Run for recommended time or until separation of desired fragments is complete.
9. Photograph using FlashGel® Camera or other standard camera and transilluminator.

Table 1. Recommended sample preparation and run conditions for FlashGel® FLX Cassettes

Sample volume	5 µL
Optimal sample concentrations	Optimal DNA load levels are 5 – 20 ng/band in a 5 µL load. For best results, do not exceed 20 ng/band. Optimal RNA load levels will vary depending upon RNA sample. For best results, do not exceed 200 ng/band in a 5 µL load. For FlashGel® RNA Markers, 2.5 µL of a 10-fold dilution is suggested.
Sample preparation	DNA Samples: for best results, dilute 4 µL DNA sample in 1 µL FlashGel® FLX DNA Prestain Loading Buffer. You may need to pre-dilute DNA samples prior to this step to obtain final optimal DNA load levels as described above. Denatured RNA samples: dilute RNA samples in RNase-free water, then combine 2.5 µL diluted RNA sample and 2.5 µL FlashGel® FLX RNA Prestain Loading Buffer. The final 5 µL load should not exceed 200 ng RNA/band. Denature for 5 minutes at 70°C.
Voltage and run time	DNA Samples: Single-tier: 275 V for 2 – 7 minutes Double-tier: 275 V for 2 – 5 minutes RNA Samples: Single-tier: 200 V for 10 minutes Double-tier: 200 V for 3 – 5 minutes
Separation range	1.2%: 100 bp – 10 kb DNA 0.5 – 9 kb RNA 2.2%: 20 bp – 1 kb DNA Separation of fragments > 4 kb will be improved by running longer at lower voltage. Detection of small fragments will be improved by running longer at lower voltage.
Recommended markers	DNA Samples: 1.2%: FlashGel® DNA Marker 100 bp – 4 kb 2.2%: FlashGel® DNA Marker 50 bp – 1.5 kb Double-tier cassettes: FlashGel® DNA Marker 100 bp – 3 kb RNA Samples: FlashGel® RNA Marker 0.5 kb – 9 kb

Ordering information

FlashGel® System – FlashGel® FLX Format

Cat. no.	Product name	Size/format
317258	FlashGel® FLX Cassettes	1.2% agarose, 12+1 well single-tier format, 9 pk
317259	FlashGel® FLX Cassettes	1.2% agarose, 16+1 well double-tier format, 9 pk
317260	FlashGel® FLX Cassettes	2.2% agarose, 12+1 well single-tier format, 9 pk
317281	FlashGel® FLX Cassettes	2.2% agarose, 16+1 well double-tier format, 9 pk
316992	FlashGel® FLX DNA Prestain Loading Buffer	5x loading buffer with DNA stain, 1 mL
316993	FlashGel® FLX RNA Prestain Loading Buffer	2x loading buffer with RNA stain, 1 mL
50473	FlashGel® DNA Marker, 100 bp – 4 kb	500 µL ready-to-load, recommended for 1.2% cassettes Band sizes: 100/200/300/500/800/1250/2000/4000 bp
57033	FlashGel® DNA Marker, 50 bp – 1.5 kb	500 µL ready-to-load, recommended for 2.2% cassettes Band sizes: 50/100/150/200/300/500/800/1500 bp
57034	FlashGel® DNA Marker, 100 bp – 3 kb	500 µL ready-to-load, recommended for double-tier cassettes Band sizes: 100/300/500/800/1500/3000 bp
50577	FlashGel® RNA Marker	0.5 bp – 9 kb, 50 µg (1 µg/ mL) Band sizes: 100/300/500/800/1500/3000 bp

Lonza Rockland, Inc. – Rockland, ME 04841

For research use only. Not for use in diagnostic procedures. All trademarks belong to Lonza Group Ltd. and its affiliates (collectively and individually, "Lonza"), and are registered in the USA, EU and/or CH, or used in common law, or belong to third-party owners and are used for only informational purposes. All third-party copyrights have been reproduced with permission from their owners. User assumes all risks of product use and handling. Lonza makes efforts to include accurate and up-to-date information. However, Lonza makes no representations or warranties, express or implied, including as to the accuracy or completeness of information, or the use or handling of these products. The user is responsible for determining if the products supplied by Lonza and the information and recommendations given by Lonza are (i) suitable for intended process or purpose, (ii) are in compliance with environmental, health and safety regulations, and other regulations in the regions and countries where they are purchased, offered for sale, marketed, sold, and used and (iii) will not infringe any third party's intellectual property rights. The user bears the sole responsibility for determining the existence of any such third-party rights, as well as obtaining any necessary licenses and approvals. For more information, including regarding legal disclaimers, Lonza's intellectual property rights, and how Lonza collects, uses and protects personal information: <http://www.lonza.com/legal>, <https://www.lonza.com/about-us/strategy/intellectual-property> and <https://www.lonza.com/privacy>

©2025 Lonza. All rights reserved.

MB-OP003 08/25

For more technical assistance, contact our Scientific Support Team:

North America

Phone: +1 800 521 0390 (toll free)

E-mail: scientific.support@lonza.com

Europe

Phone: +32 87 321 611

E-mail: scientific.support.eu@lonza.com

International

Contact your local Lonza distributor

Customer Service: +1 301 898 7025

scientific.support@lonza.com

bioscience.lonza.com

lonza.com/flashgel