SMARTLine™

NEO NITRILE™ 290

CATEGORY III

AQL 1.0, G1

User's Instructions

Catalogue number	Size	Scan to documents
10 831	Extra Small (XS/6)	
10 832	Small (S/7)	
10 833	Medium (M/8)	
10 834	Large (L/9)	
10 835	Extra Large (XL/10)	
10 836	Extra Extra Large (XXL/11)	

	<p>Prior to use, please read the instructions</p> <p>This glove complies with EN ISO 21420:2020 + A1:2024 - General test requirements and test methods for protective gloves</p>
--	--

	<p>The product complies with and is certified in accordance with the requirements of EU Regulation 2016/425. This product is EU conformity to type based on internal production control plus supervised product checks at random intervals (Module C2) by BSI Group, The Netherlands B.V. (Notified Body No: 2797), Say Building - John M. Keynesplein 9, Unit 4.2, 1066 EP Amsterdam - The Netherlands</p>
---	---

	<p>The product complies with and is certified in accordance with the requirements of Regulation 2016/425 on personal protective equipment, as amended to apply in Great Britain. This product is UKCA conformity to type based on internal production control plus supervised product checks at random intervals (Module C2) by BSI Assurance UK Limited (Approved Body No: 0086) Kitemark Court - Davy Avenue - Knowhill - Milton Keynes - MK5 8PP - United Kingdom</p>
---	--

Determination of material resistance to permeation by chemicals (EN ISO 374-1:2016 + A1:2018) and degradation (EN ISO 374-4:2019)				
	Code letter and chemical	CAS number	Permeation performance and breakthrough time	Degradation (DR) Mean in %
	K: Sodium hydroxide 40%	1310-73-2	Level 6, > 480 min.	-2.5
	P: Hydrogen peroxide 30%	7722-84-1	Level 6, > 480 min.	31.1
	T: Formaldehyde 37%	50-00-0	Level 6, > 480 min.	26.4
	O: Ammonium hydroxide 25%	1336-21-6	Level 4, > 120 min.	11.0
	M: Nitric acid 65%	7697-37-2	Level 3, > 60 min.	88.3
	S: Hydrofluoric acid 40%	7664-39-3	Level 4, > 120 min.	57.0

**Marking of gloves protecting against viruses, bacteria and fungi, as defined in EN ISO 374-5:2016
(Protective gloves against dangerous chemicals and micro-organisms
- Part 5: Terminology and performance requirements for micro-organisms risks). Protection
against viruses: tested according to ISO 16604:2004 (Procedure B) Phi-X, Herpes, HIV.**

 	Performance level	Acceptable quality level unit	Inspection level
	Level 3	< 0.65	G1
	Level 2	< 1.5	G1
	Level 1	< 4.0	S4

EN 388:2016+A1:2018 (ISO 23388:2018) Protective Glove against Mechanical Risk

 	Test	Performance Level
	6.1 - Abrasion resistance	Level 2
	6.2 - Blade cut Resistance	Test unsuitable for glove design
	6.4 - Tear Resistance	Test unsuitable for glove design
	6.5 - Puncture resistance	Test unsuitable for glove design
	6.3 - Cut Resistance	Test not performed
<p>"X": This indicates that the test was either not performed or was deemed unsuitable for the glove design</p>		

(en) Storage recommendations: Store in a cool and dry place. Storage temperature 5°C (41°F) to 40°C (104°F). Gloves should be protected from exposure to direct sunlight, intense artificial light, X-ray machines and other sources of ozone. Recommendation for packaging whilst being transported: Gloves should be transported in their original packaging to maintain the protective properties and usability characteristics of the gloves. Disposal: Dispose of in accordance with local authority regulations. This is an MDR Class 1 examination glove according to Regulation (EU) 2017/745, a PPE Category III protective glove according to Regulation (EU) 2016/425 and Regulation 2016/425 on Personal Protective Equipment, as amended to apply in Great Britain. Gloves provide protection against dangerous chemicals, micro-organisms and mechanical risks. For multiple layer gloves, the EN 388 performance rating applies to the entire glove and may not reflect the protection of the outermost layer. Gloves shall not be worn when there is a risk of entanglement by moving parts of machines. Caution: These gloves should not be used for thermal protection. This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400 mm - where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical is used in a mixture. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. Before usage, inspect the gloves for any defect or imperfections. The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen. DONNING: 1. Clean and prepare the work surface before opening the packaging. 2. Remove first glove from its packaging and check for any possible mechanical damage like holes or tears before use. 3. Touching the inner cuff area, pick up the glove and don pulling it all the way up to the wrist or forearm. Avoid any contact with the outer area of the glove. 4. Ensure that the glove fits properly. 5. Repeat steps 1-4 for second glove. DOFFING: 1. Hold the outer edge of the gloves near the cuff/ wrist or forearm area. 2. Carefully pull the glove away from the hand, turning it inside out. Hold it in the opposite hand. 3. Slip two fingers into the cuff of the remaining gloved hand without touching the outer area of the glove. 4. Peel off the glove, inside out, until it is completely removed ensuring that the glove removed encapsulates the first glove in its palm. 5. Discard both gloves appropriately. Declaration of conformity accessible from: <http://shieldscientific.com/resource-center/glove-declaration-of-conformity>.

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SHIELD Scientific B.V. Costerweg 1-B, 6702 AA Wageningen, The Netherlands

EU	UK
NOTIFIED BODY NO: 2797	APPROVED BODY NO: 0086
BSI Group The Netherlands B.V. Say Building, John M. Keynesplein 9, Unit 4.2, 1066 EP Amsterdam, The Netherlands	BSI Assurance UK Ltd Kitemark Court, Davy Avenue, Knowlhill, Milton Keynes, MK 5 8PP, United Kingdom
MODULE B EU TYPE EXAMINATION UNDERTAKEN BY NOTIFIED BODY NO: 2797	MODULE B UKCA TYPE EXAMINATION UNDERTAKEN BY APPROVED BODY NO: 0086
BSI Group The Netherlands B.V. Say Building, John M. Keynesplein 9, Unit 4.2, 1066 EP Amsterdam, The Netherlands	BSI Assurance UK Ltd Kitemark Court, Davy Avenue, Knowlhill, Milton Keynes, MK 5 8PP, United Kingdom

“SMARTLine™, Your First Line of Defence”

www.shieldscientific.com