


Thermo Scientific
Nunc EasyFill Cell Factory System



Versatile
and easy to use

Thermo
SCIENTIFIC

Committed to perfecting cell culture at every scale

The Thermo Scientific™ Nunc™ EasyFill™ Cell Factory™ system is an out-of-the-box, multi-layer solution, as easy to use as a T-flask. Scale-up with confidence. Featuring our certified Thermo Scientific™ Nunclon™ Delta cell culture surface, which is available on most Thermo Scientific Nunc cell culture products to assure consistent performance from lot to lot and from format to format.



One EasyFill Cell Factory 10 compared
to 36 standard T175 flasks

- **Save valuable space** — each EasyFill Cell Factory 10 is equivalent to 36 T175 flasks
- **Enhance productivity** — fill and empty 5 times faster vs. T175 flasks
- **Achieve fast start-up** — to get started, simply pour your media and cells directly in the large opening
- **Versatile port** designed to facilitate pouring and aseptic filling techniques

More confidence — less risk

Certified non-cytotoxic, non-pyrogenic and adhesive-free


EasyFill Cell Factory multi-layer systems are designed for a variety of biological and commercial applications. With a certified Sterility Assurance Level of 10^{-6} achieved following ISO 11137-2 guidelines and Nunclon Delta surface treatment, you can be assured of consistent performance from layer to layer and lot to lot. Our cell culture certification includes verification of monolayer formation and cloning efficiency using three cell lines and one primary cell.

Instructions for use

Nunc EasyFill Cell Factory System



1 Pour media and cells directly into the Cell Factory system.



2 **Equilibrate** by placing the Cell Factory system on the side with the small blue port cap.
Note: Contamination during leveling is prevented as there is no risk of the filter cap from becoming wet.



3 Turn the Cell Factory system **90°**, so that the filling inlet is up. Medium will be separated with equal amounts in each chamber.



With filling inlet up, **carefully tilt** the Cell Factory system to a horizontal incubation position. **Handle by holding on to the sides - do not manipulate by holding on to the upper rim.**



5 **Incubate** following appropriate protocol.



6 **Loosen filter cap** and remove blue port cap. Empty directly into reservoir.
Note: Controlled pouring from the traditional port prevents the loss of cells or supernatant resulting from spills

For a video guide of our EasyFill Cell Factory system please visit www.thermoscientific.com/easyfill

Nunc EasyFill Cell Factory Systems Nunclon Delta Surface. Polystyrene. Sterile					
Cat. No.	No. of Trays	Cell Culture Area, cm ²	Nominal Dimensions LxWxH mm (in)	Qty per pack/case	
140000	1	632	335x205x35 (13.2 x 8.1 x 1.4)	1/6	
140250	2	1,264	335x205x53 (13.2 x 8.1 x 2.1)	1/6	
140360	4	2,528	335x205x87 (13.2 x 8.1 x 3.4)	1/4	
140400	10	6,320	335x205x190 (13.2 x 8.1 x 7.5)	1/2	
140410	10	6,320	335x205x190 (13.2 x 8.1 x 7.5)	1/6	
140440	40	25,280	335x205x698 (13.2 x 8.1 x 27.5)	1/2	

For a complete list of our Thermo Scientific Cell Factory Solutions, please visit www.thermoscientific.com/cellfactory

Each EasyFill Cell Factory system includes:



Port cover cap, blue
Cat. No. 170615



Replacement Vent Cap
Cat. No. 146008




Thermo Scientific CO₂ Incubators

Discovery thrives in a culture of confidence

From the largest capacity incubators to a small personal-sized model, there is a Thermo Scientific CO₂ incubator fit for your laboratory needs.

- Choice of volume capacities ranging from 1.4 cu. ft. (40 L) to 29 cu. ft. (821 L)
- Conveniently stackable models for space-constrained labs
- Space to accommodate shakers, stirrers, culture devices or large sample throughput
- Delivering long-term performance, optimal growth conditions and proven contamination prevention

www.thermoscientific.com/co2



www.thermoscientific.com/EasyFill

© 2014 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

ANZ: Australia: 1300 735 292, New Zealand: 0800 933 966

Asia: China Toll-free: 800-810-5118 or 400-650-5118; India: +91 22 6716 2200, India Toll-free: 1 800 22 8374; Japan: +81 3 5826 1616; Other Asian countries: 65 68729717

Europe: Austria: +43 1 801 40 0; Belgium: +32 53 73 42 41; Denmark: +45 4631 2000; France: +33 2 2803 2180; Germany: +49 6184 90 6940, Germany Toll-free: 08001-536 376; Italy: +39 02 02 95059 or 434-254-375; Netherlands: +31 76 571 4440; Nordic/Baltic countries: +358 9 329 100; Russia/CIS: +7 (812) 703 42 15; Spain/Portugal: +34 93 223 09 18; Switzerland: +41 44 454 12 12; UK/Ireland: +44 870 609 9203

North America: USA/Canada +1 585 586 8800; USA Toll-free: 800 625 4327

South America: USA sales support: +1 585 586 8800

Countries not listed: +49 6184 90 6940 or +33 2 2803 2180

BRLLSPNUNCEFCF 0814

Thermo
SCIENTIFIC

A Thermo Fisher Scientific Brand