

Biochemical identification

Giving back bench time to microbiologists

Thermo Scientific™ RapID™ Systems




No oil. No pipetting. No 24-hr incubation.

Waiting for an ID panel to do its job can be frustrating when the results are needed to make a critical, informed decision. Traditionally, manual, confirmatory ID panels provide results in 24–48 hours, because they utilize carbohydrate fermentation reactions that are dependent on organism growth.

RapID Systems work differently. They use enzyme technology to reduce the time-to-result to 4 hours. The advantage – it provides a same-day definitive answer to a significant bacterial identification faster.

RapID™ Systems are ideal for stand-alone testing or for complementing automated systems when testing for anaerobes, Enterobacteriaceae, staphylococci, yeasts, Neisseria-Haemophilus, streptococci, Corynebacteria, non-fermenters, and urinary tract isolates.



No oil required

Reactions are based on the detection of preformed bacterial enzymes, which removes the requirement for anaerobic incubation.

No pipetting

RapID panels allow for the simultaneous inoculation of each cavity, decreasing set-up time and the amount of labor involved in rapidly identifying over 400 organisms.

No large air bubbles

Surfactant is placed on the back of the panel for a smooth, fluid inoculation to minimize air bubbles in reaction cavities.

No 24 hour incubation

Panels require incubation for up to 4 hours, which facilitates same-day reporting of results



Broad range of RapID Systems

RapID ANA II

For clinically important anaerobic bacteria. Panel includes 18 substrates for the identification of more than 90 species.

RapID CB PLUS

For *Corynebacterium* spp., and other irregular, Gram-positive coryneform bacilli. Panel includes 18 substrates for the identification of more than 50 species.

RapID NH

For *Neisseria*, *Haemophilus*, *Moraxella*, and other related microorganisms. Panel includes 13 substrates for the identification of more than 30 species.

RapID NF PLUS

For glucose fermenting and nonfermenting, Gram-negative bacteria not belonging to the Enterobacteriaceae family. Panel includes 17 substrates for the identification of more than 70 species.

RapID ONE

For Enterobacteriaceae and other Gram-negative, oxidase-negative bacteria. Panel includes 19 substrates for the identification of more than 70 species.

RapID STAPH PLUS

For clinically important *Staphylococcus* species and related organisms. Panel includes 18 substrates for the identification of more than 40 species.

RapID SS/u

For commonly isolated urinary tract microorganisms from common urinary tract pathogens grown on agar. Panel includes 11 substrates for the identification of 13 species in only 2 hours.

RapID STR

For *Streptococcus* species and other related Gram-positive bacteria. Panel includes 14 substrates for the identification of more than 30 species.

RapID YEAST PLUS

For clinically important yeast, yeast-like and related microorganisms. Panel includes 18 substrates for the identification of more than 40 species.



Globally available

Our goal is to provide the highest quality products, outstanding service, fast delivery, and excellent technical support to exceed your expectations.

ThermoFisher
SCIENTIFIC

RapID Systems and product guide

Color guides are included in every kit

Product	Format	Product Code
RapID ANA II	20 panels/kit	R8311002
RapID Inoculation Fluid 1 mL [†]	20 tubes/pack	R8325102
McFarland Standard No. 3	Each	R20413
RapID Spot Indole 15 mL [†]	Each	R8309002
RapID ANA II QC Set/3	5 loops/pack	R4653050
RapID CB Plus	20 panels/kit	R8311008
RapID Inoculation Fluid 2 mL [†]	20 tubes/pack	R8325106
McFarland Standard No. 4	Each	R20414
RapID Nitrate A 15 mL [†]	Each	R8309003
RapID Nitrate B 15 mL [†]	Each	R8309004
RapID CB Plus QC Set/3	5 loops/pack	R4653048
RapID NH Plus	20 panels/kit	R8311001
RapID Inoculation Fluid 1 mL [†]	20 tubes/pack	R8325102
McFarland Standard No. 3	Each	R20413
RapID Spot Indole 15 mL [†]	Each	R8309002
RapID Nitrate A 15 mL [†]	Each	R8309003
RapID Nitrate B 15 mL [†]	Each	R8309004
Thermo Scientific™ Remel™ BactiDrop™ Oxidase 0.75 mL	50 ampules/pack	R21540
RapID NH Plus QC Set/4	5 loops/pack	R4653051
RapID NF Plus	20 panels/kit	R8311005
RapID Inoculation Fluid 1 mL [†]	20 tubes/pack	R8325102
McFarland Standard No. 1	Each	R20411
McFarland Standard No. 3	Each	R20413
RapID Spot Indole 15 mL [†]	Each	R8309002
RapID Nitrate A 15 mL [†]	Each	R8309003
BactiDrop Oxidase 0.75 mL	50 ampules/pack	R21540
RapID NF Plus QC Set/4	5 loops/pack	R4653054

[†] Accessory product required for use with RapID Systems.

Product	Format	Product Code
RapID ONE	20 panels/kit	R8311006
RapID Inoculation Fluid 2 mL [†]	20 tubes/pack	R8325106
McFarland Standard No. 2	Each	R20412
RapID Spot Indole 15 mL [†]	Each	R8309002
BactiDrop Oxidase 0.75 mL	50 ampules/pack	R21540
RapID ONE QC Set/4	5 loops/pack	R4653056
RapID Staph Plus	20 panels/kit	R8311009
RapID Inoculation Fluid 2 mL [†]	20 tubes/pack	R8325106
McFarland Standard No. 3	Each	R20413
RapID Nitrate A 15 mL [†]	Each	R8309003
RapID Nitrate B 15 mL [†]	Each	R8309004
RapID Staph Plus QC Set/4	5 loops/pack	R4653047
RapID SS/u	20 panels/kit	R8311004
RapID Inoculation Fluid 1 mL [†]	20 tubes/pack	R8325102
McFarland Standard No. 1	Each	R20411
RapID Spot Indole 15 mL [†]	Each	R8309002
RapID SS/u QC Set/5	5 loops/pack	R4653053
RapID STR	20 panels/kit	R8311003
RapID Inoculation Fluid 1 mL [†]	20 tubes/pack	R8325102
McFarland Standard No. 1	Each	R20411
RapID STR QC Set/5	5 loops/pack	R4653052
RapID Yeast Plus	20 panels/kit	R8311007
RapID Inoculation Fluid 2 mL [†]	20 tubes/pack	R8325106
RapID Yeast Plus QC Set/5	5 loops/pack	R4653060



Designed exclusively for use with RapID Systems, the latest version of ERIC™ Software features updated taxonomy for extensive coverage and confidence. All laboratories can take advantage of this easy-to-use database which is available free online or as the ERIC RapID Compendium CD.

 To learn more about our complete range of RapID products and the free online ERIC software visit: thermofisher.com/RapID

Thermo Fisher Scientific products are distributed globally so uses, applications, and availability of product in each country depend on local regulatory marketing authorization status. © 2023 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. **LT2926A August 2023**

thermo scientific