Key Code TSMX5241E
www.thermofisher.com

Europe +800 135 79 135 US 1 855 2360 190
CA 1 855 805 8539 ROW +31 20 794 7071

DrySpot Staphytest Plus

REF DR0100M.....120 Tests

EN

1. INTENDED USE

DrySpot™ Staphytest™ Plus kit is a qualitative latex slide agglutination assay for the differentiation of *Staphylococcus aureus* from other *Staphylococcus* species isolates grown on agar by the detection of clumping factor, Protein A and certain capsular polysaccharides found in methicillin-resistant *Staphylococcus aureus* (MRSA). Used in a diagnostic workflow to aid clinicians in treatment options for patients suspected of having staphylococcal infection or colonization. The device is not automated, is for professional use only and is not a companion diagnostic.

2. PRINCIPLE OF THE TEST

Traditionally, differentiation between coagulase-positive and coagulase-negative staphylococci has been performed with the tube coagulase test that detects extracellular staphylocoagulase or the slide coagulase test that detects the clumping factor (bound coagulase) present on the bacterial cell surface. Several other differentiation tests are also available including the passive haemagglutination test (Oxoid Staphylase – DR595) and the DNase test.

It has been reported that approximately 97% of human strains of *Staphylococcus aureus* possess both bound coagulase and extracellular staphylocoagulase.

Protein A is found on the cell surface of about 95% of human strains of *S. aureus* and has the ability to bind the Fc portion of immunoglobulin G (IgG)².

It has been observed that certain methicillin-resistant strains of *S. aureus* (MRSA) may express undetectable levels of clumping factor and Protein A^{3,4,5}. It has been shown, however, that these strains all possess capsular polysaccharide⁶. The capsule can mask both Protein A and the clumping factor thereby preventing agglutination.

DrySpot Staphytest Plus uses blue latex particles coated with both porcine fibrinogen and rabbit IgG including specific polyclonal antibodies raised against capsular polysaccharides of *S. aureus*^{7,8}.

The reagent is dried onto the reaction card. When the reagent is mixed on the card with colonies of *S. aureus* emulsified in saline rapid agglutination occurs through the reaction between (i) fibrinogen and clumping factor, (ii) Fc portion of IgG and Protein A, (iii) specific IgG and capsular polysaccharide. Agglutination may also occur with other species which possess clumping factor or Protein A such as *Staphylococcus hyicus* and *Staphylococcus intermedius*. If neither clumping factor, Protein A nor specific capsular polysaccharides are present, agglutination will not occur and the result will be regarded as negative. The most frequent

coagulase and Protein A negative isolates of staphylococci are *Staphylococcus epidermidis*.

3. COMPONENTS OF THE KIT

DrySpot Staphytest Plus Reagent Cards.

Blue latex particles coated with both porcine fibrinogen and rabbit IgG together with specific polyclonal antibodies raised against capsular polysaccharide of *S. aureus* (Test Reaction Area).

Blue latex particles sensitised with non-reactive globulin (Control Reaction Area).

4 pouches each containing 10 cards and a moisture absorbent sachet. There are 3 test and 3 control reaction areas on each card. 120 tests in total.

Plastic pouch clip for storage of opened pouches.

Instructions for Use (IFU).

4. MATERIALS REQUIRED BUT NOT PROVIDED:

Saline (0-85%)

Timer

Pipette or dropper (50 µl) Sterile Microbiological Loops

A suitable laboratory disinfectant

Positive Control: *S. aureus* strain such as ATCC® 25923

Negative Control: *S. epidermidis* strain such as ATCC® 12228.

5. PRECAUTIONS

This product is for *in vitro* diagnostic use only.

Specimen material may contain pathogenic organisms. Handle with the appropriate precautions.

Directions should be read and followed closely.

Please refer to the Safety Data Sheet (SDS) on company website and product labelling for information on potentially hazardous components.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

In the event of malfunction do not use device.

6. STORAGE AND OPENING

This kit must be stored between 2 and 25°C. If stored in a cold environment, allow the pouches to reach room temperature before opening to prevent condensation of moisture on the cards. The DrySpot reagents will deteriorate and may give false results if they are allowed to absorb moisture.

Open the pouches by cutting with scissors just below the seal.

Once opened remove the number of cards required for immediate testing (testing within the next 10 minutes) and reseal the pouch immediately by clamping the open end of the bag between the two halves of the plastic clip.

If fewer tests are required cut the reaction cards along the indicated lines and return the unused portions to the pouch. Do not return used portions to the pouch because they will cause contamination of remaining cards in the pouch.

Under these conditions the reagents will retain their activity until the expiry date shown on the kit box.

7. CONTROL PROCEDURES

On each occasion the kit is used, the following control procedures must be performed:

1. Positive Control: Use a known *S. aureus* strain such as ATCC® 25923 (Thermo Scientific Culti-Loops™ R4607010). Follow the method given in Test Procedure. Ensure that agglutination occurs within 20 seconds.
2. Negative Control: Use a known *S. epidermidis* strain such as ATCC® 12228 (Thermo Scientific Culti-Loops™ R4606500). Follow the method given in Test Procedure. Ensure that the reagent remains smooth and non-agglutinated for the entire 20 seconds of the test.

Do not use the test if reactions with the control organisms are incorrect.

8. IMPORTANT PROCEDURE NOTES

Do not touch the circles on the reaction cards as this may cause contamination and affect the reaction.

In a high humidity environment do not leave pouches open for more than 2 minutes. If there is evidence of moisture in the spots do not use.

Do not place the drop of saline directly onto the dry latex spots. Pouch clips can be retained for future use to allow multiple packs to be opened. Although suitable for room temperature storage the kit or pouches must not be stored near a heat source or where exposure to sunlight may cause increased temperatures.

Ensure that sufficient growth is removed from the culture plate; an insufficient quantity may give rise to false-negative reactions.

9. SPECIMEN COLLECTION AND PREPARATION

For details of specimen collection and treatment a standard test book should be consulted⁹.

Gram-positive, catalase-positive colonies may be tested from any of the following culture media:

Baird-Parker Agar
Blood Agar
Brilliance MRSA 2 Agar
Columbia Blood Agar
Columbia CNA Agar
Mannitol Salt Agar
Mueller-Hinton Agar with 5% blood
Nutrient Agar
Oxacillin Resistance Screening Agar (ORSA)
Tryptone Soya Agar
Tryptone Soya Agar with 5% blood

The use of fresh cultures grown overnight is recommended (18 to 36 hours incubation). The tendency of colonies to cause autoagglutinating reactions increases with incubation beyond the recommended 36 hour period.

10. STANDARD TEST METHOD

1. Add 1 drop (50 µl) of saline (0-85%) to the small rings (at the bottom of each oval) in both the test and control reaction areas ensuring that the liquid does not mix with the dried latex reagents.
2. Use a sterile loop to pick-up the equivalent of 5 average-sized suspect staphylococcal colonies (equivalent to 2–3 mm diameter of growth) from a culture media plate and carefully emulsify in the saline drop. Ensure that the resulting suspension is smooth.

3. Mix the suspension into the dry Control Latex spots until completely suspended and spread to cover the reaction area using the loop. Discard the loop appropriately.
4. Using a separate loop proceed in the same way with the Test Latex.
5. Pick up and rock the card for up to 20 seconds and observe for agglutination under normal lighting conditions. Do not use a magnifying glass.
6. When the test is completed dispose of reaction cards safely into disinfectant.

11. TEST METHOD FOR OXACILLIN RESISTANCE SCREENING AGAR

1. Use a sterile loop to pick-up the equivalent of 5 average-sized suspect staphylococcal colonies (equivalent to 2–3 mm diameter of growth) from a culture media plate and smear onto the Control Reaction Area. Using the loop, spread the colony material into a thin film keeping the colony material clear of the dried latex reagents.
2. Dispense 1 drop of saline (0-85%) directly onto the thin film to form a smooth suspension and mix IMMEDIATELY.

3–6. As Standard Test Method.

12. READING AND INTERPRETATION OF RESULTS

Positive Result

A result is positive if agglutination of the blue test latex particles occurs within 20 seconds. This presumptively identifies the strain as a *S. aureus*.

Negative Result

A negative result is obtained if no agglutination occurs and a smooth blue suspension remains after 20 seconds in the test circle. This presumptively identifies the strain as a non *S. aureus*.

Equivocal Result

Slight graininess of the test latex accompanied by no change in the appearance of the control latex should be interpreted as an equivocal result. Strains should be re-tested following subculture onto non-selective media.

Uninterpretable Result

The test is uninterpretable if the Control Reagent shows agglutination. This indicates that the culture causes autoagglutination.

Granular or Stringy Reactions

Occasionally granular or stringy reactions may be seen due to the particulate nature of the test material. When such reactions are seen to occur they should be interpreted using the following criteria:

The result is **positive** when, using the Test Reagent, greater clearing of the blue background is observed compared with the reaction of the Control Reagent.

The result is **negative** when there are no differences between clearing of the blue background using the Test and Control Reagents.

Reactions occurring after 20 seconds should be ignored.

13. LIMITATIONS OF THE PROCEDURE

1. The tendency of isolated colonies to cause autoagglutination increases with incubation times beyond the recommended 36 hour period.
2. The antibody used in DrySpot Staphytest Plus has been optimised to avoid potential cross-reactions with shared antigens from coagulase negative staphylococci. It should be noted that this has been shown to reduce sensitivity to some type 18 MRSA strains¹⁰.

3. Some species of staphylococci other than *S. aureus* notably *S. hyicus*, *S. intermedius*, *S. lugdunensis*, *S. xylosus*, *S. schleiferi* and *S. haemolyticus*^{11,12,13,14}, may give positive results in coagulase tests and/or rapid latex procedures. If necessary the species may be identified by biochemical test procedures e.g. using a test for PYRase activity (Oxoid O.B.I.S. PYR ID0580M), *S. aureus* and *S. hyicus* will be PYRase negative and all the other strains named above will be positive^{15,16}. *S. hyicus* and *S. intermedius* are rarely encountered in the clinical laboratory.
4. Staphylococci isolated from urine specimens which give a weak positive¹⁷ result with Staphytest Plus may be *Staphylococcus saprophyticus*. Further identification of such isolates may be conducted using biochemical tests and novobiocin sensitivity (*S. saprophyticus* is resistant to novobiocin).
5. Some streptococci and possibly other organisms possessing immunoglobulin or plasma binding factors may react in the latex test and some species such as *Escherichia coli* are able to agglutinate latex particles^{18,19} non-specifically. To overcome these non-specific results a Gram-stain should be performed so only typical staphylococci are tested.

14. PERFORMANCE CHARACTERISTICS (DATA ON FILE)

The performance characteristics of Oxoid DrySpot Staphytest Plus have been determined using data from the studies detailed below. It is important to note, however, that *S. aureus* is known to show considerable antigenic variation with respect to different geographical locations.

Clinical Study

Oxoid DrySpot Staphytest Plus was evaluated at a large Australian teaching hospital. A total of 300 isolates was tested using tube coagulase as the gold standard method. In the following data analysis, results from strains known to be cross reactors^{11,12,13,14} and results from autoagglutinating strains have been omitted (n=284). The relative sensitivity was 100% and the relative specificity was 96·3%.

Industrial Study

Oxoid DrySpot Staphytest Plus was evaluated in food laboratories in a multi-centre study in the United Kingdom. A total of 621 samples was evaluated, these were colonies taken from Baird-Parker Agar, which had been inoculated with food or environmental material. The gold standard method was tube coagulase. In the following data analysis, results from strains known to be cross reactors^{11,12,13,14} and results from autoagglutinating strains have been omitted (n=603). The relative sensitivity was 98·4% and the relative specificity was 96·9%.

15. REFERENCES:

1. Essers, L. and Radebold, K. (1980). "Rapid and Reliable Identification of *Staphylococcus aureus* by a Latex Agglutination Test". J.Clin.Microbiol. 12: 641–643.

2. Taussig, M. J. (1984). Processes in Pathology and Microbiology. 2nd Edn. 520–530. Blackwell, Oxford.

3. Ruane, P. J., Morgan, M. A., Citron, D. M. and Mulligan, M. E. (1986). "Failure of Rapid Agglutination Methods to Detect Oxacillin-Resistant *Staphylococcus aureus*". J.Clin.Microbiol. 24: 490–492.

4. Roberts, J. I. S. and Gaston, M. A. (1987). "Protein A and coagulase expression in epidemic and non-epidemic *Staphylococcus aureus*". J.Clin.Pathol. 40: 837–840.

5. Wanger, A. R., Morris, S. L., Ericsson, C., Singh, K. V. and LaRocco,

M. T. (1992). "Latex Agglutination-Negative Methicillin-Resistant *Staphylococcus aureus* Recovered from Neonates: Epidemiologic Features and Comparison of Typing Methods". J.Clin.Microbiol. 30: 2583–2588.

6. Fournier, J. M., Boutonnier, A. and Bouvet, A. (1989). "*Staphylococcus aureus* Strains Which Are Not Identified by 7 Rapid Agglutination Methods Are of Capsular Serotype 5". J.Clin.Microbiol. 27: 1372–1374.

7. Fournier, J. M., Bouvet, A., Boutonnier, A., Audurier, A., Goldstein, F., Pierre, J., Bure, A., Lebrun, L. and Hochkeppel, H. K. (1987). Predominance of Capsular Polysaccharide Type 5 among Oxacillin-Resistant *Staphylococcus aureus*". J.Clin.Microbiol. 25: 1932–1933.

8. Karakawa, W. W., Fournier, J. M., Vann, W. F., Arbeit, R., Schneerson, R. S. and Robbins, J. B. (1985). "Method for the Serological Typing of the Capsular Polysaccharides of *Staphylococcus aureus*". J.Clin.Microbiol. 22: 445–447.

9. Kloos, W. E. and Jorgensen, J. H. (1988). Staphylococci. pp. 143–153. In Manual of Clinical Microbiology. 4th Edn. (Eds) Lennette, E. H., Balows, A., Hausser, W. J. and Shadomy, H. J.: Assoc. Amer. Microbiol. Washington.

10. Data on file at Oxoid Ltd.

11. Jean-Pierre, H., Darbas, H., Jean-Rousseng, A. and Boyer, G. (1989). "Pathogenicity in Two Cases of *Staphylococcus schleiferi*, a Recently Described Species". J.Clin.Microbiol. 27: 2110–2111.

12. Freney, J., Brun, Y., Bes, M., Heugnier, H., Grimont, F., Grimont, P. A. D., Nervie, C. and Fleurette, J. (1988). "*Staphylococcus lugdunensis* sp. nov. and *Staphylococcus schleiferi* sp. nov., Two species from Human Clinical Specimens". Int.J.Sup.Bacteriol. 38: 168–172.

13. Phillips, W. E. and Kloos, W. E. (1981). "Identification of Coagulase-Positive *Staphylococcus intermedius* and *Staphylococcus hyicus* subsp. *hyicus* isolates from Veterinary Clinical Specimens". J.Clin.Microbiol. 14: 671–673.

14. van Griethuysen, A., Bes, M., Etienne, J., Zbinden, R. and Kluytmans, J. (2000). "An International Multicenter Evaluation of a new Latex Agglutination Test for Identification of *Staphylococcus aureus*". J.Clin. Microbiol. 39: 86–89.

15. Schnitzler, N., Rainer, M., Conrads, G., Frank, D. and Haase, G. (1988). "*Staphylococcus lugdunensis*: Report of a case of Peritonitis and an Easy-To-Perform Screening Strategy". J.Clin.Microbiol. 26: 1939–1949.

16. Ann-Herbert, A., Crowder, C. G., Hancock, G. A., Jarvis, W. R. and Thornsberry, C. (1998). "Characteristics of Coagulase- Negative-Staphylococci That Help Differentiate These Species of the Family *Micrococcaceae*". J.Clin.Microbiol. 36: 812–813.

17. Gregson, D. B., Low, D. E., Skulnick, M. and Simor, A. E. (1988). "Problems with Rapid Agglutination Methods for 8 Identification of *Staphylococcus aureus* When *Staphylococcus saprophyticus* Is Being Tested". J.Clin.Microbiol. 26: 1398–1399.

18. Myhre, E. B. and Kuusela, P. (1983). "Binding of Human Fibronectin to Group A, C and G Streptococci". Infect.Immun. 40: 29–34.

19. Runeheggen, A., Schonbeck, C., Hedner, U., Hessel, B. and Kronvall, G. (1981). "Binding of Fibrinogen Degradation Products to *S. aureus* and to beta-hemolytic Streptococci Group A, C and G". Acta Pathologica Microbiologica Scandinavica Section B Microbiology (now called APMS), Sect B. 89: 49–55.

16. SYMBOL DEFINITIONS

	Catalogue Number
	In Vitro Diagnostic Medical Device
	Consult Instructions for Use (IFU)
	Temperature Limitations (Storage temp.)
	Contains sufficient for <N> tests
	Do not use if package is damaged
	Do not re-use
	Batch Code (Lot Number)
	Use By (Expiration Date)
	Importer
	Unique Device Identifier
	Authorized representative in the European Community
	UK Conformity Assessed
	European Conformity Assessment
	Manufacturer

For technical assistance www.thermofisher.com

 Oxoid Limited, Wade Road, Hampshire
Basingstoke, RG24 8PW, UK
+44 (0) 1256 841144 thermofisher.com

Version	Date of modifications introduced
X5241E	April 2024 Updated to conform to IVDR

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6. Fournier, J. M., Boutonnier, A. and Bouvet, A. (1989). “*Staphylococcus aureus* Strains Which Are Not Identified by 7 Rapid Agglutination Methods Are of Capsular Serotype 5”. J.Clin.Microbiol. 27: 1372–1374.

7. Fournier, J. M., Bouvet, A., Boutonnier, A., Audurier, A., Goldstein, F., Pierre, J., Bure, A., Lebrun, L. and Hochkeppel, H. K. (1987). Predominance of Capsular Polysaccharide Type 5 among Oxacillin-Resistant *Staphylococcus aureus*. J.Clin.Microbiol. 25: 1932–1933.

8. Karakawa, W. W., Fournier, J. M., Vann, W. F., Arbeit, R., Schneerson, R. S. and Robbins, J. B. (1985). “Method for the Serological Typing of the

Capsular Polysaccharides of *Staphylococcus aureus*”. J.Clin.Microbiol. 22: 445–447.

9. Kloos, W. E. and Jorgensen, J. H. (1988). Staphylococci. pp. 143–153. In Manual of Clinical Microbiology. 4th Edn. (Eds) Lennette, E. H., Balows, A., Hausser, W. J. and Shadomy, H. J.: Assoc. Amer. Microbiol. Washington.

10. Data on file at Oxoid Ltd.

11. Jean-Pierre, H., Darbas, H., Jean-Rousseny, A. and Boyer, G. (1989). “Pathogenicity in Two Cases of *Staphylococcus schleiferi*, a Recently Described Species”. J.Clin.Microbiol. 27: 2110–2111.

12. Freney, J., Brun, Y., Bes, M., Heugnier, H., Grimont, F., Grimont, P. A. D., Nervie, C. and Fleurette, J. (1988). “*Staphylococcus lugdunensis* sp. nov. and *Staphylococcus schleiferi* sp. nov., Two species from Human Clinical Specimens”. Int.J.Sup.Bacteriol. 38: 168–172.

13. Phillips, W. E. and Kloos, W. E. (1981). “Identification of Coagulase-Positive *Staphylococcus intermedius* and *Staphylococcus hyicus* subsp. *hyicus* isolates from Veterinary Clinical Specimens”. J.Clin.Microbiol. 14: 671–673.

14. van Griethuysen, A., Bes, M., Etienne, J., Zbinden, R. and Kluytmans, J. (2000). “An International Multicenter Evaluation of a new Latex Agglutination Test for Identification of *Staphylococcus aureus*”. J.Clin. Microbiol. 39: 86–89.

15. Schnitzler, N., Rainer, M., Conrads, G., Frank, D. and Haase, G. (1988). “*Staphylococcus lugdunensis*: Report of a case of Peritonitis and an Easy-To-Perform Screening Strategy”. J.Clin.Microbiol. 26: 1939–1949.

16. Ann-Herbert, A., Crowder, C. G., Hancock, G. A., Jarvis, W. R. and Thornsberrry, C. (1998). “Characteristics of Coagulase- Negative-Staphylococci That Help Differentiate These Species of the Family *Micrococcaceae*”. J.Clin.Microbiol. 36: 812–813.

17. Gregson, D. B., Low, D. E., Skulnick, M. and Simor, A. E. (1988). “Problems with Rapid Agglutination Methods for 8 Identification of *Staphylococcus aureus* When *Staphylococcus saprophyticus* Is Being Tested”. J.Clin.Microbiol. 26: 1398–1399.

18. Myhre, E. B. and Kuusela, P. (1983). “Binding of Human Fibronectin to Group A, C and G Streptococci”. Infect.Immun. 40: 29–34.

19. Runehegen, A., Schonbeck, C., Hedner, U., Hessel, B. and Kronvall, G. (1981). “Binding of Fibrinogen Degradation Products to *S. aureus* and to -Hemolytic Streptococci Group A, C and G”. Acta.path.microbiol. Scand., Sect B. 89: 49–55.

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6. Fournier, J. M., Boutonnier, A. and Bouvet, A. (1989). “*Staphylococcus aureus* Strains Which Are Not Identified by 7 Rapid Agglutination Methods Are of Capsular Serotype 5”. J.Clin.Microbiol. 27: 1372–1374.
7. Fournier, J. M., Bouvet, A., Boutonnier, A., Audurier, A., Goldstein, F., Pierre, J., Bure, A., Lebrun, L. and Hochkeppel, H. K. (1987). Predominance of Capsular Polysaccharide Type 5 among Oxacillin-Resistant *Staphylococcus aureus*”. J.Clin.Microbiol. 25: 1932–1933.
8. Karakawa, W. W., Fournier, J. M., Vann, W. F., Arbeit, R., Schneerson, R. S. and Robbins, J. B. (1985). “Method for the Serological Typing of the Capsular Polysaccharides of *Staphylococcus aureus*”. J.Clin. Microbiol. 22: 445–447.
9. Kloos, W. E. and Jorgensen, J. H. (1988). Staphylococci. pp. 143–153. In Manual of Clinical Microbiology. 4th Edn. (Eds) Lennette, E. H., Balows, A., Hausser, W. J. and Shadomy, H. J.: Assoc. Amer. Microbiol. Washington.
10. Reģistrētie dati uzņēmumā Oxoid Ltd.
11. Jean-Pierre, H., Darbas, H., Jean-Rousseng, A. and Boyer, G. (1989). “Pathogenicity in Two Cases of *Staphylococcus schleiferi*, a Recently Described Species”. J.Clin.Microbiol. 27: 2110–2111.
12. Freney, J., Brun, Y., Bes, M., Heugnier, H., Grimont, F., Grimont, P. A. D., Nervie, C. and Fleurette, J. (1988). “*Staphylococcus lugdunensis* sp. nov. and *Staphylococcus schleiferi* sp. nov., Two species from Human Clinical Specimens”. Int.J.Sup.Bacteriol. 38: 168–172.
13. Phillips, W. E. and Kloos, W. E. (1981). “Identification of Coagulase-Positive *Staphylococcus intermedius* and *Staphylococcus hyicus* subsp. *hyicus* isolates from Veterinary Clinical Specimens”. J.Clin. Microbiol. 14: 671–673.
14. van Griethuysen, A., Bes, M., Etienne, J., Zbinden, R. and Kluytmans, J. (2000). “An International Multicenter Evaluation of a new Latex Agglutination Test for Identification of *Staphylococcus aureus*”. J.Clin. Microbiol. 39: 86–89.
15. Schnitzler, N., Rainer, M., Conrads, G., Frank, D. and Haase, G. (1988). “*Staphylococcus lugdunensis*: Report of a case of Peritonitis and an Easy-To-Perform Screening Strategy”. J.Clin.Microbiol. 26: 1939–1949.
16. Ann-Herbert, A., Crowder, C. G., Hancock, G. A., Jarvis, W. R. and Thornsberrry, C. (1998). “Characteristics of Coagulase- Negative-Staphylococci That Help Differentiate These Species of the Family *Micrococcaceae*”. J.Clin.Microbiol. 36: 812–813.
17. Gregson, D. B., Low, D. E., Skulnick, M. and Simor, A. E. (1988). “Problems with Rapid Agglutination Methods for 8 Identification of *Staphylococcus aureus* When *Staphylococcus saprophyticus* Is Being Tested”. J.Clin.Microbiol. 26: 1398–1399.
18. Myhre, E. B. and Kuusela, P. (1983). “Binding of Human Fibronectin to Group A, C and G Streptococci”. Infect.Immun. 40: 29–34.
19. Runehegen, A., Schonbeck, C., Hedner, U., Hessel, B. and Kronvall, G. (1981). “Binding of Fibrinogen Degradation Products to *S. aureus* and to beta-hemolytic Streptococci Group A, C and G”. Acta Pathologica Microbiologica Scandinavica Section B Microbiology (now called APMIS)., Sect B. 89: 49–55.

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- Jean-Pierre, H., Darbas, H., Jean-Roussenq, A. and Boyer, G. (1989). "Pathogenicity in Two Cases of *Staphylococcus schleiferi*, a Recently Described Species". J.Clin.Microbiol. 27: 2110–2111.
- Freney, J., Brun, Y., Bes, M., Heugnier, H., Grimont, F., Grimont, P. A. D., Nerville, C. and Fleurette, J. (1988). "*Staphylococcus lugdunensis* sp. nov. and *Staphylococcus schleiferi* sp. nov., Two species from Human Clinical Specimens". Int.J.Sup.Bacteriol. 38: 168–172.
- Phillips, W. E. and Kloos, W. E. (1981). "Identification of Coagulase-Positive *Staphylococcus intermedius* and *Staphylococcus hyicus* subsp. *hyicus* isolates from Veterinary Clinical Specimens". J.Clin.Microbiol. 14: 671–673.
- van Griethuysen, A., Bes, M., Etienne, J., Zbinden, R. and Kluytmans, J. (2000). "An International Multicenter Evaluation of a new Latex Agglutination Test for Identification of *Staphylococcus aureus*". J.Clin. Microbiol. 39: 86–89.
- Schnitzler, N., Rainer, M., Conrads, G., Frank, D. and Haase, G. (1988). "*Staphylococcus lugdunensis*: Report of a case of Peritonitis and an Easy-To-Perform Screening Strategy". J.Clin.Microbiol. 26: 1939–1949.
- Ann-Herbert, A., Crowder, C. G., Hancock, G. A., Jarvis, W. R. and Thornsberry, C. (1998). "Characteristics of Coagulase- Negative-Staphylococci That Help Differentiate These Species of the Family *Micrococcaceae*". J.Clin.Microbiol. 36: 812–813.
- Gregson, D. B., Low, D. E., Skulnick, M. and Simor, A. E. (1988). "Problems with Rapid Agglutination Methods for 8 Identification of *Staphylococcus aureus* When *Staphylococcus saprophyticus* Is Being Tested". J.Clin.Microbiol. 26: 1398–1399.
- Myhre, E. B. and Kuusela, P. (1983). "Binding of Human Fibrinectin to Group A, C and G Streptococci". Infect.Immun. 40: 29–34.
- Runeheggen, A., Schonbeck, C., Hedner, U., Hessel, B. and Kronvall, G. (1981). "Binding of Fibrinogen Degradation Products to *S. aureus* and to beta-hemolytic Streptococci Group A, C and G". Acta Pathologica Microbiologica Scandinavica Section B Microbiology (now called APMIS), Sect B. 89: 49–55.

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6. Fournier, J. M., Boutonnier, A. and Bouvet, A. (1989). "*Staphylococcus aureus* Strains Which Are Not Identified by 7 Rapid Agglutination Methods Are of Capsular Serotype 5". J.Clin.Microbiol. 27: 1372-1374.

7. Fournier, J. M., Bouvet, A., Boutonnier, A., Audurier, A., Goldstein, F., Pierre, J., Bure, A., Lebrun, L. and Hochkeppel, H. K. (1987). Predominance of Capsular Polysaccharide Type 5 among Oxacillin-Resistant *Staphylococcus aureus*". J.Clin.Microbiol. 25: 1932-1933.

8. Karakawa, W. W., Fournier, J. M., Vann, W. F., Arbeit, R., Schneerson, R. S. and Robbins, J. B. (1985). "Method for the Serological Typing of the Capsular Polysaccharides of *Staphylococcus aureus*". J.Clin. Microbiol. 22: 445-447.

9. Kloos, W. E. and Jorgensen, J. H. (1988). Staphylococci. pp. 143–153. In Manual of Clinical Microbiology. 4th Edn. (Eds) Lennette, E. H., Balows, A., Hausser, W. J. and Shadomy, H. J.: Assoc. Amer. Microbiol. Washington.

10. Údaje v evidencii spoločnosti Oxoid Ltd.

11. Jean-Pierre, H., Darbas, H., Jean-Rousseng, A. and Boyer, G. (1989). "Pathogenicity in Two Cases of *Staphylococcus schleiferi*, a Recently Described Species". J.Clin.Microbiol. 27: 2110-2111.

12. Freney, J., Brun, Y., Bes, M., Heugnier, H., Grimont, F., Grimont, P. A. D., Nervie, C. and Fleurette, J. (1988). "*Staphylococcus lugdunensis* sp. nov. and *Staphylococcus schleiferi* sp. nov., Two species from Human Clinical Specimens". Int.J.Sup.Bacteriol. 38: 168-172.

13. Phillips, W. E. and Kloos, W. E. (1981). "Identification of Coagulase-Positive *Staphylococcus intermedius* and *Staphylococcus hyicus* subsp. *hyicus* isolates from Veterinary Clinical Specimens". J.Clin.Microbiol. 14: 671-673.

14. van Griethuysen, A., Bes, M., Etienne, J., Zbinden, R. and Kluytmans, J. (2000). "An International Multicenter Evaluation of a new Latex Agglutination Test for Identification of *Staphylococcus aureus*". J.Clin. Microbiol. 39: 86-89.

15. Schnitzler, N., Rainer, M., Conrads, G., Frank, D. and Haase, G. (1988). "*Staphylococcus lugdunensis*: Report of a case of Peritonitis and an Easy-To-Perform Screening Strategy". J.Clin.Microbiol. 26: 1939-1949.

16. Ann-Herbert, A., Crowder, C. G., Hancock, G. A., Jarvis, W. R. and Thornsberry, C. (1998). "Characteristics of Coagulase- Negative-Staphylococci That Help Differentiate These Species of the Family Micrococcaceae". J.Clin.Microbiol. 36: 812-813.

17. Gregson, D. B., Low, D. E., Skulnick, M. and Simor, A. E. (1988). "Problems with Rapid Agglutination Methods for 8 Identification of *Staphylococcus aureus* When *Staphylococcus saprophyticus* Is Being Tested". J.Clin.Microbiol. 26: 1398-1399.

18. Myhre, E. B. and Kuusela, P. (1983). "Binding of Human Fibronectin to Group A, C and G Streptococci". Infect.Immun. 40: 29-34.

19. Runeheggen, A., Schonbeck, C., Hedner, U., Hessel, B. and Kronvall, G. (1981). "Binding of Fibrinogen Degradation Products to *S. aureus* and to beta-hemolytic Streptococci Group A, C and G". Acta Pathologica Microbiologica Scandinavica Section B Microbiology (now called APMIS), Sect B. 89: 49-55.

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