Key Code TSMX5326D
www.thermofisher.com

Europe 00800 135 79 135 US 1 855 2360 910
CA 1 855 805 8539 ROW +31 20 794 7071

DrySpot E.coli O157 Latex Test

REF DR0120M..... 120

1. INTENDED USE

A qualitative latex slide agglutination test for the confirmatory identification of *E. coli* serogroup O157 isolates grown on agar. The device is used in a diagnostic workflow to aid clinicians in treatment options for patients suspected of having bacterial infections. The device is not automated, is for professional use only and is not a companion diagnostic.

2. PRINCIPLE OF THE TEST

Enterohaemorrhagic *Escherichia coli* belong to a number of O serogroups with O157 being most significant in human disease. The potentially verocytotoxin (VT)-producing strains are associated with a range of symptoms from non-bloody diarrhoea, fever and vomiting to cases of haemorrhagic colitis (HC) and haemolytic uraemic syndrome (HUS).^{1,2,3,4,5,6,7,8}

The DrySpot *E. coli* O157 latex test will demonstrate by slide agglutination *E. coli* strains possessing the O157 serogroup antigen. The test is best used in conjunction with Sorbitol MacConkey Agar (Oxoid CM0813) or Sorbitol MacConkey Agar with Cefixime Tellurite Supplement (SR0172).

The majority of *E. coli* O157:H7 strains do not ferment sorbitol and therefore give colourless colonies on this medium. The majority of non-O157 *E. coli* isolates do ferment sorbitol and give characteristic pink colonies.

Sorbitol MacConkey Agar should be used as the primary screen. Non-sorbitol fermenting colonies can then be tested with the latex reagents, to determine whether the isolate belongs to the O157 serogroup and therefore is a potential VT-producing strain.

Production of verotoxin may be confirmed with the Oxoid VTEC-RPLA kit (TD960A).

3. COMPONENTS OF THE KIT (DR0120M)

DR0121M DrySpot E. coli O157 Test Reagent Cards
A suspension of synthetic blue latex particles coated with antibody specifically reactive with the *E. coli* O157 serogroup antigen (Test reaction area). A suspension of synthetic blue latex coated with pre-immune rabbit globulins (Control reaction area). Four pouches each containing 10 cards and a moisture absorbent sachet. There are 3 test and 3 control reaction areas on each card – 120 tests in total.

DR0122M Positive control strips (10 sticks – pink spots)
Pink-dyed inactivated antigenic extract of *E. coli* O157.

DR0123M Negative control strips (10 sticks – green spots)
Green-dyed inactivated antigenic extract of *E. coli* O116.

Mixing paddles (2 bags)
Pouch Clips (2)
Instructions for Use (IFU)

Materials required but not provided
Saline (0.9% NaCl prepared using distilled or deionised water)
Microbiological loop and bunsen burner
A suitable laboratory disinfectant
A timer.

4. PRECAUTIONS

IVD

This product is for *in vitro* diagnostic use only.

Specimen materials may contain pathogenic organisms. Handle with the appropriate precautions.


For professional use only.

In the event of malfunction do not use device.

Please refer to the Material Safety Data Sheet (MSDS), available on ThermoFisher website, and product labeling for information on potentially hazardous components.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

5. STORAGE AND OPENING

This kit must be stored between 2°C and 25°C. If stored in a cold environment, allow pouches to equilibrate to room temperature before opening to prevent condensation of moisture on the cards. The dry reagents will deteriorate and give false results if they are allowed to absorb moisture.

Open the pouches by cutting with scissors just below the seal. Once opened, remove the number of cards required for immediate testing (testing within the next 10 minutes) and re-seal the pouch immediately by clamping the open end of the bag between the two halves of the plastic clip.

If fewer tests are required cut the reaction cards along the indicated lines and return the unused portions to the pouch. Do not return used portions because they will cause contamination of remaining cards in the pouch.

The control sticks are also provided in a moisture impermeable pouch. Ensure that the same techniques are used to avoid moisture damage.

Under these conditions the reagents will retain their activity until the expiry date shown on the kit box.

6. CONTROL PROCEDURES

The dried control sticks provided should be used in the following way to check the correct working of the latex reagent each day before routine tests are performed.

Add a 50 µl drop of saline to the small circle at the base of the test oval reaction area. Remove a positive control card stick from the pouch by tearing one off from the others on the strip, take care to avoid touching the flexible end where the dry spots are located. Re-seal the inner bag and pouch. Turn the stick over so that the coloured spots are at the bottom and place the stick on the card with the spots touching the liquid. Push down so that the end bends at the hinge and mix in a circular manner for 10 seconds to rehydrate the dried reagent. Continue to use the stick to mix the liquid into the DrySpot O157 test reagent until all the reagents are fully rehydrated and homogeneous. Rock the card and look for agglutination. This procedure should be repeated using a negative control stick.

The positive control DR0122M must show agglutination with the dried test reagent within 1 minute.

The negative control DR0123M must show no agglutination within 1 minute.

Do not use the test if reactions with the control reagents are incorrect.

7. IMPORTANT PROCEDURE NOTES

Do not touch the circles on the reaction cards as this may cause contamination and thereby affect the reaction.

In a high humidity environment do not leave pouches open for more than 2 minutes. If there is evidence of moisture in the spots do not use.

Do not place the drop of saline directly onto the dry latex spots.

Pouch clips should be retained for future use to allow multiple packs to be opened. Although suitable for room temperature storage the kit or pouches must not be stored near heat sources. They should not be stored where exposure to sunlight may cause increased temperatures.

The controls provided should be used to check the correct working of the latex reagents each day before routine tests are performed.

The positive control must cause visible agglutination with the latex reagent within 1 minute.

The negative control must not cause agglutination within 1 minute.

Do not use the test if reactions with the control suspensions are incorrect.

8. CULTURE MATERIAL

Non-sorbitol fermenting (NSF) colonies may be taken from Sorbitol MacConkey Agar or Sorbitol MacConkey Agar with Cefixime Tellurite Supplement.

Alternatively NSF isolates may be inoculated onto non-selective media such as Nutrient Agar for testing. It is necessary to test up to 10 separate NSF colonies to ensure a high probability of detecting any O157 strains which may be in mixed culture with NSF *E. coli* of other serotypes. The use of the control latex will ensure that the isolate is not an autoagglutinating strain.

9. TEST METHOD

1. Add 1 drop (50 µl) of saline (0.9%) to the small ring (at the bottom of each oval) in both the test and control reaction areas ensuring that the liquid does not mix with the dried latex reagents.
2. Use a sterile loop (or one of the paddles provided) to pick-up a portion of the suspect colony from a culture media plate and carefully emulsify in the saline drop. Ensure that the resulting suspension is smooth.

3. Using the provided paddle mix the suspension into the dry latex spots until completely suspended and spread to cover the reaction area. Discard the loop-paddle appropriately.
4. Using a separate loop or paddle, proceed in the same way with the Control Latex.
5. Pick up and rock the card for up to 60 seconds and look for agglutination under normal lighting conditions. Do not use a magnifying glass.
6. When the test is completed dispose of reaction cards safely into disinfectant.

10. READING AND INTERPRETATION OF RESULTS

Positive Result
A result is positive if agglutination of the latex particles occurs within 1 minute.

This indicates the presence of *E. coli* serogroup O157.

Negative Result
A negative result is obtained if no agglutination occurs and a smooth blue suspension remains after 60 seconds in the test area.

Reactions occurring after 60 seconds should be ignored.

Uninterpretable Result
The test is uninterpretable if the control reagent shows agglutination. This indicates that the culture causes autoagglutination.

Granular or Stringy Reactions
Occasionally granular or stringy reactions may be seen due to the particulate nature of the test material. When such reactions are seen to occur they should be interpreted using the following criteria:

The result is positive when, using the test reagent, greater clearing of the blue background occurs compared with the reaction of the control reagent. The result is negative, when there are no significant changes in the clearing of the blue background using the test and control reagents.

If stringiness is found to be too severe for a correct judgement to be made then another colony should be suspended in 0.3 ml of saline. Allow the lumps to settle and re-test this smooth supernatant.

11. LIMITATIONS OF THE TEST

Mixing paddles provided are not sterile. These may be sterilised locally if required.

If a positive result is obtained on a colony of unknown species, then biochemical tests should be performed to confirm that the organism is an *E. coli* strain.

Neither the Sorbitol MacConkey Agar nor the *E. coli* O157 latex test will directly confirm the isolate as a toxin-producing strain.

Other serotypes have been found which produce the verocytotoxin.

Strains of *Escherichia hermanii* cross react with *E. coli* O157 sera and the latex test due to a shared antigen.¹

In addition other bacterial strains have been identified with epitopes that mimic those of *E. coli* O157 lipopolysaccharide.¹⁰

E. hermanii may be differentiated from *E. coli* by the former's fermentation of cellobiose growth in the presence of potassium cyanide (KCN), negative reaction in the 4-Methylumbelliferone glucuronide (MUG) assay and yellow pigmentation, which may be delayed.¹¹

12. PERFORMANCE CHARACTERISTICS

The performance of DrySpot *E. coli* O157 kit was evaluated at one clinical microbiology laboratory in Great Britain. Previously characterised strains of *E. coli* O157 and other non-sorbitol fermenting enterobacteriaceae were tested with this kit and with two other commercially available latex tests.

The sensitivity of DrySpot *E. coli* O157 kit was calculated to be 98.3%.

The specificity of DrySpot *E. coli* O157 kit was calculated to be 100%.

The results of the evaluation are summarised below:

Organism	Total tested	DrySpot E. coli O157 test result	
		Positive	Negative
<i>E. coli</i> O157	60	59	1*
Other non-sorbitol fermenting enterobacteriaceae	30	0	30

* A false negative reaction was recorded for an *E. coli* O157:H9 (VT-) which gave a negative result with all latex kits tested.

13. REFERENCES

1. Borczyk, A., Lior, H. and Ciebin, B. (1987). False positive identifications of *Escherichia coli* O157 in foods. *Int. J. Food Microbiol.* 4, 347–349.

2. Konowalchuk, J., Speirs, J. and Stavric, S. (1977). Vero Response to a Cytotoxin of *Escherichia coli*. *Vero Infect. Immun.* 18, 775–779.

3. Scotland, S., Day, N. and Rowe, B. (1980). Production of a Cytotoxin Affecting Vero Cells by strains of *Escherichia coli* belonging to Traditional Enteropathogenic Serogroups. *FEMS Microbiol. Lett.* 7, 15–17.

4. Centers for Disease Control. (1982). Isolation of *E. coli* O157:H7 from Sporadic Cases of Hemorrhagic colitis United States. *Morbid. Mortal. Wkly* 31, 580–585.

5. Karmali, M., Still, B. et al. (1983). Sporadic cases of Haemolytic-Uraemic Syndrome Associated with Faecal Cytotoxin and Cytotoxin-Producing *Escherichia coli* in stools. *Lancet* i, 619–620.

6. Johnson, W., Lior, H. and Bezanson, G. S. (1983). Cytotoxic *Escherichia coli* O157:H7 Associated with Haemorrhagic colitis in Canada. *Lancet* i, 76.

7. March, S. and Ratnam (1986). Sorbitol-MacConkey Medium for Detection of *Escherichia coli* O157:H7 Associated with Hemorrhagic Colitis. *J. Clin. Microbiol.* 23, 869–872.

8. Krishnan, C., Fitzgerald, et al. (1987). Laboratory Investigation of Outbreak of Haemorrhagic Colitis caused by *Escherichia coli* O157:H7. *J. Clin. Microbiol.* 25, 1043–1047.

9. Chapman, P. A. (1989). Evaluation of Commercial Latex Slide test for Identifying *Escherichia coli* O157. *J. Clin. Pathol.* 42, 1109–1110.

10. Westerman, R. B., Yongsheng, H. et al. (1997). Production and Characterisation of Monoclonal Antibodies Specific for the Lipopolysaccharide of *Escherichia coli* O157. *J. Clin. Microbiol.* 35, 679–684.

11. Brenner et al. (1982) Atypical biogroups of *Escherichia coli* found in clinical specimens and description of *Escherichia*

14. SYMBOL LEGEND

	Catalogue Number
	In Vitro Diagnostic Medical Device
	Consult Instructions for Use (IFU)
	Temperature Limitations (Storage temp.)
	Contains sufficient for <N> tests
	Batch Code (Lot Number)
	Use By (Expiration Date)
	Importer
	Unique Device Identifier
	Authorized representative in the European Community
	UK Conformity Assessed
	European Conformity Assessment
	Manufacturer

Oxoid Limited, Wade Road, Basingstoke,
Hampshire, RG24 8PW, UK

For technical assistance www.thermofisher.com

Version	Date of modifications introduced
X5326D	April 2024
	Updated to meet IVDR requirements

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13. REFERENCE

1. Borczyk, A., Lior, H. and Ciebin, B. (1987). False positive identifications of *Escherichia coli* O157 in foods. Int. J. Food Microbiol. 4, 347–349.
2. Konowalchuk, J., Speirs, J. and Stavric, S. (1977). Vero Response to a Cytotoxin of *Escherichia coli*. Vero Infect. Immun. 18, 775–779.
3. Scotland, S., Day, N. and Rowe, B. (1980). Production of a Cytotoxin Affecting Vero Cells by strains of *Escherichia coli* belonging to Traditional Enteropathogenic Serogroups. FEMS Microbiol. Lett. 7, 15–17.
4. Centers for Disease Control. (1982). Isolation of *E. coli* O157:H7 from Sporadic Cases of Hemorrhagic colitis United States. Morbid. Mortal. Wkly 31, 580–585.
5. Karmali, M., Still, B. et al. (1983). Sporadic cases of Haemolytic-Uraemic Syndrome Associated with Faecal Cytotoxin and Cytotoxin-Producing *Escherichia coli* in stools. Lancet i, 619–620.
6. Johnson, W., Lior, H. and Bezanson, G. S. (1983). Cytotoxic *Escherichia coli* O157:H7 Associated with Haemorrhagic colitis in Canada. Lancet i, 76.
7. March, S. and Ratnam (1986). Sorbitol-MacConkey Medium for Detection of *Escherichia coli* O157:H7 Associated with Hemorrhagic Colitis. J. Clin. Microbiol. 23, 869–872.
8. Krishnan, C., Fitzgerald, et al. (1987). Laboratory Investigation of Outbreak of Haemorrhagic Colitis caused by *Escherichia coli* O157:H7. J. Clin. Microbiol. 25, 1043–1047.

9. Chapman, P. A. (1989). Evaluation of Commercial Latex Slide test for Identifying *Escherichia coli* O157. J. Clin. Pathol. 42, 1109–1110.
10. Westerman, R. B., Yongsheng, H. et al. (1997). Production and Characterisation of Monoclonal Antibodies Specific for the Lipopolysaccharide of *Escherichia coli* O157. J. Clin. Microbiol. 35, 679–684.
11. Brenner et al. (1982) Atypical biogroups of *Escherichia coli* found in clinical specimens and description of *Escherichia hermannii* sp. nov Journal of Clinical Microbiology sv. 15 č.4 str.703

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13. VIITED

- Borczyk, A., Lior, H. ja Ciebin, B. (1987). False positive identifications of *Escherichia coli* O157 in foods. *Int. J. Food Microbiol.* 4, 347–349.
- Konowalchuk, J., Speirs, J. ja Stavric, S. (1977). Vero Response to a Cytotoxin of *Escherichia coli*. *Vero Infect. Immun.* 18, 775–779.
- Scotland, S., Day, N. ja Rowe, B. (1980). Production of a Cytotoxin Affecting Vero Cells by strains of *Escherichia coli* belonging to Traditional Enteropathogenic Serogroups. *FEMS Microbiol. Lett.* 7, 15–17.
- Centers for Disease Control. (1982). Isolation of *E. coli* O157:H7 from Sporadic Cases of Hemorrhagic colitis United States. *Morbid. Mortal. Wkly* 31, 580–585.
- Karmali, M., Still, B. et al. (1983). Sporadic cases of Haemolytic-Uraemic Syndrome Associated with Faecal Cytotoxin and Cytotoxin-Producing *Escherichia coli* in stools. *Lancet* i, 619–620.
- Johnson, W., Lior, H. ja Bezanson, G. S. (1983). Cytotoxic *Escherichia coli* O157:H7 Associated with Haemorrhagic colitis in Canada. *Lancet* i, 76.
- March, S. ja Ratnam (1986). Sorbitol-MacConkey Medium for Detection of *Escherichia coli* O157:H7 Associated with Hemorrhagic Colitis. *J. Clin. Microbiol.* 23, 869–872.
- Krishnan, C., Fitzgerald, et al. (1987). Laboratory Investigation of Outbreak of Haemorrhagic Colitis caused by *Escherichia coli* O157:H7. *J. Clin. Microbiol.* 25, 1043–1047.
- Chapman, P. A. (1989). Evaluation of Commercial Latex Slide test for Identifying *Escherichia coli* O157. *J. Clin. Pathol.* 42, 1109–1110.
- Westerman, R. B., Yongsheng, H. et al. (1997). Production and Characterisation of Monoclonal Antibodies Specific for the Lipopolysaccharide of *Escherichia coli* O157. *J. Clin. Microbiol.* 35, 679–684.
- Brenner et al. (1982) Atypical biogroups of *Escherichia coli* found in clinical specimens and description of *Escherichia hermannii* sp. nov *Journal of Clinical Microbiology*, Vol. 15 No.4, lk 703

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1. Borczyk, A., Lior, H. and Ciebin, B. (1987). False positive identifications of *Escherichia coli* O157 in foods. Int. J. Food Microbiol. 4, 347–349.
2. Konowalchuk, J., Speirs, J. and Stavric, S. (1977). Vero Response to a Cytotoxin of *Escherichia coli*. Vero Infect. Immun. 18, 775–779.
3. Scotland, S., Day, N. and Rowe, B. (1980). Production of a Cytotoxin Affecting Vero Cells by strains of *Escherichia coli* belonging to Traditional Enteropathogenic Serogroups. FEMS Microbiol. Lett. 7, 15–17.
4. Centers for Disease Control. (1982). Isolation of *E. coli* O157:H7 from Sporadic Cases of Hemorrhagic colitis United States. Morbid. Mortal. Wkly 31, 580–585.
5. Karmali, M., Still, B. et al. (1983). Sporadic cases of Haemolytic-Uraemic Syndrome Associated with Faecal Cytotoxin and Cytotoxin-Producing *Escherichia coli* in stools. Lancet i, 619–620.
6. Johnson, W., Lior, H. and Bezanson, G. S. (1983). Cytotoxic *Escherichia coli* O157:H7 Associated with Haemorrhagic colitis in Canada. Lancet i, 76.
7. March, S. and Ratnam (1986). Sorbitol-MacConkey Medium for Detection of *Escherichia coli* O157:H7 Associated with Hemorrhagic Colitis. J. Clin. Microbiol. 23, 869–872.
8. Krishnan, C., Fitzgerald, et al. (1987). Laboratory Investigation of Outbreak of Haemorrhagic Colitis caused by *Escherichia coli* O157:H7. J. Clin. Microbiol. 25, 1043–1047.
9. Chapman, P. A. (1989). Evaluation of Commercial Latex Slide test for Identifying *Escherichia coli* O157. J. Clin. Pathol. 42, 1109–1110.

10. Westerman, R. B., Yongsheng, H. et al. (1997). Production and Characterisation of Monoclonal Antibodies Specific for the Lipopolysaccharide of *Escherichia coli* O157. J. Clin. Microbiol. 35, 679–684.

11. Brenner et al. (1982) Atypical biogroups of *Escherichia coli* found in clinical specimens and description of *Escherichia hermannii* sp. nov Journal of Clinical Microbiology Vol. 15 No.4 p703

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- Borczyk, A., Lior, H. and Ciebin, B. (1987). False positive identifications of *Escherichia coli* O157 in foods. Int. J. Food Microbiol. 4, 347–349.
- Konowalchuk, J., Speirs, J. and Stavric, S. (1977). Vero Response to a Cytotoxin of *Escherichia coli*. Vero Infect. Immun. 18, 775–779.
- Scotland, S., Day, N. and Rowe, B. (1980). Production of a Cytotoxin Affecting Vero Cells by strains of *Escherichia coli* belonging to Traditional Enteropathogenic Serogroups. FEMS Microbiol. Lett. 7, 15–17.
- Centers for Disease Control. (1982). Isolation of *E. coli* O157:H7 from Sporadic Cases of Hemorrhagic colitis United States. Morbid. Mortal. Wkly 31, 580–585.
- Karmali, M., Still, B. et al. (1983). Sporadic cases of Haemolytic-Uraemic Syndrome Associated with Faecal Cytotoxin and Cytotoxin-Producing *Escherichia coli* in stools. Lancet i, 619–620.
- Johnson, W., Lior, H. and Bezanson, G. S. (1983). Cytotoxic *Escherichia coli* O157:H7 Associated with Haemorrhagic colitis in Canada. Lancet i, 76.
- March, S. and Ratnam (1986). Sorbitol-MacConkey Medium for Detection of *Escherichia coli* O157:H7 Associated with Hemorrhagic Colitis. J. Clin. Microbiol. 23, 869–872.
- Krishnan, C., Fitzgerald, et al. (1987). Laboratory Investigation of Outbreak of Haemorrhagic Colitis caused by *Escherichia coli* O157:H7. J. Clin. Microbiol. 25, 1043–1047.
- Chapman, P. A. (1989). Evaluation of Commercial Latex Slide test for Identifying *Escherichia coli* O157. J. Clin. Pathol. 42, 1109–1110.
- Westerman, R. B., Yongsheng, H. et al. (1997). Production and Characterisation of Monoclonal Antibodies Specific for the Lipopolysaccharide of *Escherichia coli* O157. J. Clin. Microbiol. 35, 679–684.
- Brenner et al. (1982) Atypical biogroups of *Escherichia coli* found in clinical specimens and description of *Escherichia hermannii* sp. nov Journal of Clinical Microbiology Vol. 15 No.4 p703

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13. REFERENCES

1. Borczyk, A., Lior, H. and Ciebin, B. (1987). False positive identifications of Escherichia coli O157 in foods. Int. J. Food Microbiol. 4, 347–349.
2. Konowalchuk, J., Speirs, J. and Stavric, S. (1977). Vero Response to a Cytotoxin of Escherichia coli. Vero Infect. Immun. 18, 775–779.
3. Scotland, S., Day, N. and Rowe, B. (1980). Production of a Cytotoxin Affecting Vero Cells by strains of Escherichia coli belonging to Traditional Enteropathogenic Serogroups. FEMS Microbiol. Lett. 7, 15–17.
4. Centers for Disease Control. (1982). Isolation of *E. coli* O157:H7 from Sporadic Cases of Hemorrhagic colitis United States. Morbid. Mortal. Wkly 31, 580–585.
5. Karmali, M., Still, B. et al. (1983). Sporadic cases of Haemolytic-Uraemic Syndrome Associated with Faecal Cytotoxin and Cytotoxin-Producing Escherichia coli in stools. Lancet i, 619–620.
6. Johnson, W., Lior, H. and Bezanson, G. S. (1983). Cytotoxic Escherichia coli O157:H7 Associated with Haemorrhagic colitis in Canada. Lancet i, 76.
7. March, S. and Ratnam (1986). Sorbitol-MacConkey Medium for Detection of Escherichia coli O157:H7 Associated with Hemorrhagic Colitis. J. Clin. Microbiol. 23, 869–872.
8. Krishnan, C., Fitzgerald, et al. (1987). Laboratory Investigation of Outbreak of Haemorrhagic Colitis caused by Escherichia coli O157:H7. J. Clin. Microbiol. 25, 1043–1047.
9. Chapman, P. A. (1989). Evaluation of Commercial Latex Slide test for Identifying Escherichia coli O157. J. Clin. Pathol. 42, 1109–1110.
10. Westerman, R. B., Yongsheng, H. et al. (1997). Production and Characterisation of Monoclonal Antibodies Specific for the Lipopolysaccharide of Escherichia coli O157. J. Clin. Microbiol. 35, 679–684.
11. Brenner et al (1982) Journal of Clinical Microbiology Vol. 15 No. 4 p.703

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2. Konowalchuk, J., Speirs, J. and Stavric, S. (1977). Vero Response to a Cytotoxin of Escherichia coli. Vero Infect. Immun. 18, 775–779.

3. Scotland, S., Day, N. and Rowe, B. (1980). Production of a Cytotoxin Affecting Vero Cells by strains of Escherichia coli belonging to Traditional Enteropathogenic Serogroups. FEMS Microbiol. Lett. 7, 15–17.

4. Centers for Disease Control. (1982). Isolation of *E. coli* O157:H7 from Sporadic Cases of Hemorrhagic colitis United States. Morbid. Mortal. Wkly 31, 580–585.

5. Karmali, M., Still, B. et al. (1983). Sporadic cases of Haemolytic-Uraemic Syndrome Associated with Faecal Cytotoxin and Cytotoxin-Producing Escherichia coli in stools. Lancet i, 619–620.

6. Johnson, W., Lior, H. and Bezanson, G. S. (1983). Cytotoxic Escherichia coli O157:H7 Associated with Haemorrhagic colitis in Canada. Lancet i, 76.

7. March, S. and Ratnam (1986). Sorbitol-MacConkey Medium for Detection of Escherichia coli O157:H7 Associated with Hemorrhagic Colitis. J. Clin. Microbiol. 23, 869–872.

8. Krishnan, C., Fitzgerald, et al. (1987). Laboratory Investigation of Outbreak of Haemorrhagic Colitis caused by Escherichia coli O157:H7. J. Clin. Microbiol. 25, 1043–1047.

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Oxoid Ltd Wade Road,
Basingstoke, Hampshire,
RG24 8PW, Regno Unito

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3. Scotland, S., Day, N. and Rowe, B. (1980). Production of a Cytotoxin Affecting Vero Cells by strains of *Escherichia coli* belonging to Traditional Enteropathogenic Serogroups. *FEMS*

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8. Krishnan, C., Fitzgerald, et al. (1987). Laboratory Investigation of Outbreak of Hemorrhagic Colitis caused by *Escherichia coli* O157:H7. *J. Clin. Microbiol.* 25, 1043–1047.

9. Chapman, P. A. (1989). Evaluation of Commercial Latex Slide test for Identifying *Escherichia coli* O157. *J. Clin. Pathol.* 42, 1109–1110.

10. Westerman, R. B., Yongsheng, H. et al. (1997). Production and Characterisation of Monoclonal Antibodies Specific for the Lipopolysaccharide of *Escherichia coli* O157. *J. Clin. Microbiol.* 35, 679–684.

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Oxoid Ltd Wade Road,
Basingstoke, Hampshire,
RG24 8PW, Reino Unido

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