Key Code TSMX6501E
www.thermofisher.com

Europe + 800 135 79 135
CA 1 855 805 8539

US 1 855 236 0910
ROW +31 20 794 7071

Oxoid Biochemical Identification System (O.B.I.S.) mono

REF ID0600M.....▽60

INTENDED USE

The Oxoid Biochemical Identification System (O.B.I.S.) mono is a rapid, qualitative colorimetric test for the detection of D-alanyl aminopeptidase (DALAase). It is used the presumptive differentiation of *Listeria monocytogenes* from other *Listeria* species grown on plated media. The device is used in a diagnostic workflow to aid clinicians in the treatment options for patients suspected of having bacterial infections.

The device is not automated, is for professional use only and is not a companion diagnostic.

PRINCIPLE OF THE TEST

The O.B.I.S. mono test offers a rapid screening method for differentiation of *Listeria monocytogenes* from other *Listeria* species. This reduces the need for full biochemical identification of all suspect colonies.

Listeria species, with the exception of *Listeria monocytogenes*, possess the enzyme D-alanyl aminopeptidase^{1,2,3}. D-alanyl-7-amido-4-methylcoumarin (DALA) is provided as a suspension. An acidic solution of dimethylaminocinnamaldehyde is used as a colour developer. If the substrate is hydrolysed by DALAase, free 7-amino-4-methylcoumarin (7AMC) combines with the developer to produce a purple Schiff's base⁴.

COMPONENTS OF THE O.B.I.S. MONO KIT (ID0600M)

Each O.B.I.S. mono kit contains the following reagents with enough material for 60 tests:

ID090M O.B.I.S. Reaction Sleeves. One pouch containing 30 plastic bags (76mm x 102mm).

ID0121M O.B.I.S. - Buffer (DALA). One white capped bottle containing 4.5ml of a 0.5% w/v D-alanyl-7-amido-4-methylcoumarin suspension in distilled water.

ID0220M O.B.I.S. - Developer (DALA). One purple capped bottle containing 4.5ml of a 0.5% w/v dimethylaminocinnamaldehyde in 1M hydrochloric acid.

ID0601M O.B.I.S. mono Test Cards. One pouch containing 10 cards. There are six reaction areas per card, labelled 'DALA'.

Materials required but not included

Sterile plastic disposable inoculating loops.

Positive and negative quality control organisms.

37°C ± 2°C incubator.

PRECAUTIONS

This product is for *in vitro* diagnostic use only.

Do not use the O.B.I.S. mono reagents beyond stated expiry date.

Specimen material may contain pathogenic organisms, handle with appropriate precautions.

The O.B.I.S. mono Developing Solution contains acid. Wear suitable protective equipment. If the reagents come into contact with the skin, mucous membranes or eyes, immediately flush the area with water.

Used O.B.I.S. Test Cards, reaction sleeves and inoculating loops should be disposed of as biohazardous waste. This should be incinerated, or autoclaved at 121°C for at least 15 minutes.

Please refer to the Safety Data Sheet (SDS) on company website and product labelling for information on potentially hazardous components.

Directions should be read and followed carefully.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

In the event of malfunction do not use device.

STORAGE AND OPENING

The O.B.I.S. mono kit must be stored at 2°C to 8°C. Allow the pouches to equilibrate to room temperature before use to prevent the formation of condensation on the Test Cards.

Open the pouches by cutting at the notch between the end seal and the clip-lock opening.

Remove the number of Test Cards required and reseal the pouch.

If fewer tests are required than the number on the Test Card, cut the card and return the unused portion to the pouch. Do not return used Test Cards to the pouch.

When stored as described, O.B.I.S. mono reagents will retain their activity until the expiry date shown on the box.

QUALITY CONTROL PROCEDURE

Each day the kit is used the following procedure should be performed:

1. **Positive Control** – Use a known DALAase positive strain such as *Listeria innocua* ATCC™ 33090 (Thermo Scientific™ Culti Loops™ R4609005). Follow the method given in the test procedure.

Ensure that a purple colour forms within 20 seconds.

2. **Negative Control** – Use a known DALAase negative strain such as *Listeria monocytogenes* ATCC™ 7644 (Thermo Scientific™ Culti-Loops™ R4603970). Follow the method given in the test procedure.

Ensure that no purple colour forms within 20 seconds.

SPECIMENS

The test is designed for use from purity plates, not from primary isolation media, as the colonies on primary isolation media are too small to carry out an effective test.

Pick colonies which have typical *Listeria* morphology from selective *Listeria* isolation media such as Oxford Agar (CM0856), PALCAM Agar (CM0877) and chromogenic *Listeria* Media, and streak onto a purity plate.

O.B.I.S. mono tests must be performed from purity plates recommended by international standards, such as Tryptone Soya Agar (CM0131), Tryptone Soya Yeast Extract Agar (TSA-YE)5,6 or Oxoid Chromogenic *Listeria* Media.

TEST PROCEDURE AND INTERPRETATION OF RESULTS

Important notes:

Before use shake the O.B.I.S. Buffer (DALA) (white capped bottle) well to suspend the reagent.

Ensure that the culture to be tested is a Gram-positive, catalase-positive, oxidase-negative bacterium.

- Remove one of the O.B.I.S. mono Test Cards from the pack.
- Using a plastic loop, take colonial material equivalent to

five colonies (each greater than 1mm in diameter) from the purity plate.

- Transfer the colonial material on to the circular test area, ensuring the material is smeared thinly and evenly within the circular test area.
- Shake the O.B.I.S. - Buffer (DALA) well and dispense one drop onto the inoculated reaction zone.
- Place the Test Card in a plastic Reaction Sleeve, so that the reaction zones are within the sleeve and the portion of the Test Card with the Oxoid logo is outside the sleeve.
- Place the Test Card and Reaction Sleeve into a 37°C ± 2°C incubator for 10 minutes.
- Remove from incubator. Dispense one drop of O.B.I.S. - Developer (DALA) onto each of the inoculated reaction zones.
- The appearance of a purple colour within 20 seconds, is a positive D-alanyl aminopeptidase reaction. A positive reaction indicates the organism is NOT a *Listeria monocytogenes*.
- If no colour develops within 20 seconds this is a negative reaction, and indicates the organism is a presumptive *Listeria monocytogenes*.

INTERPRETATION CHART

Typical reactions:

Organism	DALAase
<i>Listeria monocytogenes</i>	–
Other <i>Listeria</i> species	+

Trial Results

Results from pure culture:

Organism	No. of Strains Tested	Positive DALA Reaction
<i>Listeria monocytogenes</i>	40	–
<i>Listeria seeligeri</i>	9	+*
<i>Listeria grayii</i>	2	+
<i>Listeria innocua</i>	18	+
<i>Listeria ivanovii</i>	10	+
<i>Listeria murrayi</i>	1	+
<i>Listeria welshimeri</i>	1	+

*8/9

In an independent trial of 34 food samples O.B.I.S. mono gave the following results:

Sensitivity	100%
Specificity*	99%

LIMITATIONS OF THE TEST

O.B.I.S. mono is intended for the detection of DALAase in Gram positive, catalase-positive, oxidase-negative short rod shaped bacteria, capable of growing on selective *Listeria* primary isolation media. It can be used as a screen to differentiate *Listeria monocytogenes* from other *Listeria*.

Occasionally aesculin positive *Bacillus* species may grow on *Listeria* isolation media. Bacilli are DALAase positive, but the colonies are different from *Listeria* and bacilli appear as large rods upon Gram staining.

Due to the small size of *Listeria* colonies on primary isolation media it is not possible to carry out the test directly from these plates. Use of multiple colonies from the primary isolation plate is not recommended as this may lead to a mixed culture and an

incorrect result. International standards recommend subculturing presumptive *Listeria* species onto purity plates TSA (CM0131), TSAYE^{5,6}, or Oxoid Chromogenic *Listeria* Media.

Occasionally certain *Listeria* when grown on non-selective media form colonies that are difficult to remove with an inoculating loop. Ensure that there is sufficient colonial material on the inoculating loop before inoculating the reaction area. Failure to pick up sufficient material on the loop will result in a negative DALA test.

O.B.I.S. mono provides presumptive identification of *Listeria monocytogenes*, but does not replace full biochemical testing.

REFERENCES

- Kämpfer, P., Böttcher, S., Dott, W. and Rüden, H. (1991). Physiological characterization and identification of *Listeria* species. Zentralblatt für Bakteriologie 275, 423–435.
- Kämpfer, P. (1992). Differentiation of *Corynebacterium* spp., *Listeria* spp., and related organisms by fluorogenic substrates. Journal of Clinical Microbiology, 30, 1067–1071.
- Clark, A. G. and McLaughlin, J. (1997). Simple color tests based on an alanyl peptidase reaction which differentiate *Listeria monocytogenes* from other *Listeria* species. Journal of Clinical Microbiology 35, 2155–2156.
- Druggan, P., Roberts, P. B. and Swaine, D. (1999). A rapid chromogenic method for the differentiation of *Citrobacter* spp. and *Salmonella* spp. directly from enteric media. Abstract of the Annual Meeting of the American Society for Microbiology 1999, C444, p. 71. Washington: ASM Press.
- International Standards Organisation (1996). ISO 11290–1, Microbiology of food and animal feeding stuffs – horizontal method for the detection and enumeration of *Listeria monocytogenes* – part 1: detection method. London: British Standards Institute.
- FDA (1998). Appendix 3: Trypticase Soy Agar with 0.6% Yeast Extract (TSAYE). In Food and Drug Administration Bacteriological Analytical Manual, 8th edition, revision A ed. Jackson, G. J. pp. 54–55. Gaithersburg: AOAC International.

SYMBOL LEGEND

REF	Catalogue Number
IVD	In Vitro Diagnostic Medical Device
	Consult Instructions for Use (IFU)
	Temperature Limitations (Storage temp.)
	Contains sufficient for <N> tests
	Do not use if package is damaged
	Do not re-use
LOT	Batch Code (Lot Number)
	Use By (Expiration Date)
	Importer
EC REP	Authorized representative in the European Community
UK CA	UK Conformity Assessed
CE	European Conformity Assessment
	Manufacturer

IFU X6501E Revised May 2025

Printed in the UK

Oxoid Ltd
Wade Road, Basingstoke
Hants, RG24 8PW, UK

For technical assistance please contact your local distributor.

©2025 Thermo Fisher Scientific Inc. All rights reserved.

ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Biochemical Identification System (O.B.I.S.) mono

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LITTERATUR

- Kämpfer, P., Böttcher, S., Dott, W. and Rüden, H. (1991). Physiological characterization and identification of *Listeria* species. Zentralblatt für Bakteriologie 275, 423–435.
- Kämpfer, P. (1992). Differentiation of *Corynebacterium* spp., *Listeria* spp., and related organisms by fluorogenic substrates. Journal of Clinical Microbiology, 30, 1067–1071.
- Clark, A. G. and McLaughlin, J. (1997). Simple color tests based on an alanyl peptidase reaction which differentiate *Listeria* monocytogenes from other *Listeria* species. Journal of Clinical Microbiology 35, 2155–2156.
- Druggan, P., Roberts, P. B. and Swaine, D. (1999). A rapid chromogenic method for the differentiation of *Citrobacter* spp. and *Salmonella* spp. directly from enteric media. Abstract of the Annual Meeting of the American Society for Microbiology 1999, C444, p. 71. Washington: ASM Press.
- International Standards Organisation (1996). ISO 11290–1, Microbiology of food and animal feeding stuffs – horizontal method for the detection and enumeration of *Listeria monocytogenes* – part 1: detection method. London: British Standards Institute.
- FDA (1998). Appendix 3: Trypticase Soy Agar with 0.6% Yeast Extract (TSAYE). In Food and Drug Administration Bacteriological Analytical Manual, 8th edition, revision A ed. Jackson, G. J. pp. 54–55. Gaithersburg: AOAC International.

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 Oxoid Ltd
Wade Road, Basingstoke
Hampshire, RG24 8PW, UK

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©2025 Thermo Fisher Scientific Inc. Med enerett.
Alle andre varemerker tilhører Thermo Fisher Scientific Inc. og dets datterselskaper.