remel EN RapID™ ONE System

REF R8311006..... 20

1. INTENDED USE

Remel RapID™ ONE System is a qualitative micromethod using enzyme reactions to identify clinical isolates grown on agar of oxidase negative, Gram-negative Enterobacteriaceae. Used in a diagnostic workflow to aid clinicians in treatment options for patients suspected of having bacterial infections. The device is not automated, is for professional use only and is not a companion diagnostic.

A complete listing of the organisms addressed by the RapID ONE system is provided in the RapID ONE Differential Chart.

2. SUMMARY AND EXPLANATION

RapID ONE System is comprised of (1) RapID ONE Panels and (2) RapID ONE Reagent. RapID ONE Panels are disposable plastic trays with 18 reaction cavities, which contain dehydrated reactants. The panel allows for the simultaneous inoculation of each cavity with a predetermined amount of inoculum. A suspension of the test organism in RapID Inoculation Fluid is used as the inoculum which rehydrates and initiates test reactions. After incubation of the panel, each test cavity is examined for reactivity by noting the development of a color. In some cases, reagents must be added to the test cavities to provide a color change. The resulting pattern of positive and negative test scores is used as the basis for identification of the test isolate by comparison of test results with the probability values in the Differential Chart (Table 4), or by use of the RapID ERIC™ software.

3. PRINCIPLE

The tests used in the RapID ONE System are based on microbial degradation of specific substrates detected by various indicator systems. The reactions employed are a combination of conventional tests and single-substrate chromogenic tests, described in Table 1.

4. REAGENTS

RapID ONE Reagent (provided with kit)	(15 ml/Btl)
Reactive ingredient per liter:	
<i>p</i> -Dimethylaminocinnamaldehyde	0.06 g
RapID Inoculation Fluid (R8325106, supplied separately) (2 ml/tube)	
KCl	6.0 g
CaCl ₂	0.5 g
Deminerlized Water	1000.0 ml
RapID Spot Indole Reagent (R8309002, supplied separately) (15 ml/Btl)	
<i>p</i> -Dimethylaminocinnamaldehyde	10.0 g
Hydrochloric Acid	100.0 ml
Deminerlized Water	900.0 ml

5. PRECAUTIONS

This product is for *in vitro* diagnostic use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, media, and test panels after use. Directions should be read and followed carefully.

Do not use reagents beyond the printed expiration dates.

Do not use if there is evidence of contamination or other signs of deterioration.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

In the event of malfunction do not use device.

Caution!

- RapID ONE Reagent is toxic and may cause harm to the environment. Harmful by inhalation, contact with skin or eyes, or if swallowed. May impair fertility or cause harm to unborn child.
- RapID Spot Indole Reagent may cause irritation to skin, eyes, and respiratory system.
- Refer to Safety Data Sheet for detailed information on reagent chemicals.

Composition / information on ingredients

2-Methoxyethanol 109-86-4
Acetic acid 64-19-7
Hydrochloric acid 7647-01-0

WARNING! This product contains a chemical known in the State of California to cause birth defects or other reproductive harm.

Emergency Telephone Number

INFOTRAC - 24 Hour Number: 1-800-535-5053

Outside of the United States, call 24 Hour Number: 001-352-323-3500 (Call Collect)

DANGER		
	H315	Causes skin irritation
	H319	Causes serious eye irritation
US ONLY	H335	May cause respiratory irritation
	H336	May cause drowsiness or dizziness
	H360	May damage fertility. May damage the unborn child
	H373	May cause damage to organs through prolonged or repeated exposure
	P201	Obtain special instructions before use
	P202	Do not handle until all safety precautions have been read and understood
	P281	Use personal protective equipment as required
	P264	Wash face, hands and any exposed skin thoroughly after handling
	P280	Wear eye/face protection
	P260	Do not breathe dust/fume/gas/mist/vapors/spray
	P271	Use only outdoors or in a well-ventilated area
	P308+313	IF exposed or concerned: Get medical advice/attention
	P304+P340	IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.
	P302+P352	IF ON SKIN: Wash with plenty of soap and water
P332+P313	If skin irritation occurs: Get medical advice/attention	
P362	Take off contaminated clothing and wash before reuse	
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing	
P337+P313	If eye irritation persists: Get medical advice/attention	
P405	Store locked up	
P403+P233	Store in a well-ventilated place. Keep container tightly closed	
P501	Dispose of contents/container to an approved waste disposal plant	

6. STORAGE

RapID ONE System and RapID Spot Indole Reagent should be stored in their original containers at 2-8°C until use. Allow products to equilibrate to room temperature before use. DO NOT interchange reagents among different RapID systems. Remove only the number of panels necessary for testing. Reseal the plastic pouch and promptly return to 2-8°C. Panels must be used the same day they are removed from storage. RapID Inoculation Fluid should be stored in its original container at room temperature (20-25°C) until used.

7. PRODUCT DETERIORATION

This product should not be used if (1) the expiration date has passed, (2) the plastic tray is broken or the lid is compromised, or (3) there are other signs of deterioration.

8. SPECIMEN COLLECTION, STORAGE, AND TRANSPORT

Specimens should be collected and handled following recommended guidelines.^{8,9}

9. MATERIALS SUPPLIED

- 20 RapID ONE Panels
- 20 report forms
- RapID ONE Reagent (one plastic dropper-bottle containing reagent sufficient for 20 panels),
- 2 chipboard incubation trays
- 1 color guide
- Instructions for use (IFU).

10. CONTENTS SYMBOLS

ONE Panels	ONE Panels
Report Forms	RapID Report Forms
ONE Reagent	ONE Reagent
Incubation Trays	Incubation Trays

11. MATERIALS REQUIRED BUT NOT SUPPLIED

- Loop sterilization device
- Inoculating loop, swab, collection containers,
- Incubators, alternative environmental systems,
- Supplemental media
- Quality control organisms
- Gram stain reagents
- Microscope slides
- Oxidase reagent
- Cotton swabs
- RapID Inoculation Fluid - 2 ml (R8325106)
- McFarland #2 turbidity standard or equivalent (R20412)
- Pipettes
- RapID Spot Indole Reagent (R8309002)
- ERIC (R8323600) (optional).

12. PROCEDURE

Inoculum Preparation:

- Test organisms must be grown in pure culture and examined by Gram stain and oxidase prior to use in the system.

Notes:

- Only oxidase-negative, Gram-negative bacilli should be tested using the RapID ONE System. Oxidase-positive bacilli should be tested using the RapID NF Plus System (R8311005).
 - The oxidase test should be interpreted with caution when using bacterial growth from differential agars that contain dyes which may interfere with interpretation.
- Test organisms may be removed from a variety of selective and nonselective agar growth media. The following types of media are recommended:
 - Nonselective Media: Tryptic Soy Agar (TSA) with blood.
 - Differential or Selective Media: Hektoen Enteric (HE) Agar; MacConkey Agar; Eosin Methylene Blue (EMB) Agar; Desoxycholate Agar; Salmonella-Shigella (SS) Agar.

Notes:

- Plates used for inoculum preparation should preferably be 18-24 hours old. Slow-growing isolates may be tested using 48-hour plates.
- The use of media other than those recommended may compromise test performance.


- Using a cotton swab or inoculating loop, suspend sufficient growth from the agar plate culture in RapID Inoculation Fluid (2 ml) to achieve a visual turbidity equal to a #2 McFarland turbidity standard or equivalent.

Notes:

- Suspensions significantly less turbid than a #2 McFarland standard will result in aberrant reactions.
 - Bacterial suspensions that are slightly more turbid than a #2 McFarland standard will not affect test performance and are recommended for stock cultures and quality control strains. However, suspensions prepared with a turbidity far greater than a #2 McFarland standard will compromise test performance.
 - Suspensions should be mixed thoroughly and vortexed if required.
 - Suspensions should be used within 15 minutes of preparation.
- An agar plate may be inoculated for purity and any additional testing that may be required using a loopful of the test suspension from the inoculation fluid tube. Incubate the plate for at least 18-24 hours at 35-37°C.

Inoculation of RapID ONE Panels:


- Peel back the lid of the panel over the inoculation port by pulling the tab marked "Peel to Inoculate" up and to the left.
- Using a pipette, gently transfer the entire contents of the Inoculation Fluid tube into the upper right-hand corner of the panel. Reseal the inoculation port of the panel by pressing the peel-back tab back in place.
- After adding the test suspension, and while keeping the panel on a level surface, tilt the panel back away from the reaction cavities at approximately 45° (see below).

- While tilted back, gently rock the panel from side to side to evenly distribute the inoculum along the rear baffles as illustrated below.


Table 1. Principles and Components of the RapID ONE System

Cavity #	Test Code	Reactive Ingredient	Quantity	Principle	Bibliography #
1	URE	Urea	0.25%	Hydrolysis of urea produces basic products which raise the pH and change the indicator.	1, 3
2	ADH	Arginine	1.0%		
3	ODC	Ornithine	1.0%		
4	LDC	Lysine	1.0%		
5	TET	Aliphatic thiol	0.2%	Utilization of the thiol compound produces acidic products which lower the pH and change the indicator.	4
6	LIP	Fatty acid ester	1.0%	Hydrolysis of the fatty acid ester releases acidic products which lower the pH and change the indicator.	3, 4
7	KSF	Sugar aldehyde	1.0%	Utilization of the carbohydrate substrate produces acidic products which lower the pH and change the indicator.	3, 4
8	SBL	Sorbitol	1.0%		
9	GUR	<i>p</i> -Nitrophenyl-β, D-glucuronide	0.1%		
10	ONPG	<i>o</i> -Nitrophenyl-β, D-galactoside	0.1%		
11	βGLU	<i>p</i> -Nitrophenyl-β,D-glucoside	0.1%	Enzymatic hydrolysis of the colorless aryl-substituted glycoside or phosphoester releases yellow <i>o</i> - or <i>p</i> -nitrophenyl.	1, 3, 4-7
12	βXYL	<i>p</i> -Nitrophenyl-β,D-xyloside	0.1%		
13	NAG	<i>p</i> -Nitrophenyl-N-acetyl-β, D- glucosaminide	0.1%		
14	MAL	Malonate	0.5%	Utilization of malonate produces basic products which raise the pH and change the indicator	1,3
15	PRO	Proline-β-naphthylamide	0.1%	Enzymatic hydrolysis of the arylamide substrate releases free β-naphthylamine which is detected with the RapID ONE Reagent.	2, 7
16	GGT	γ-Glutamyl- β- naphthylamide	0.1%		
17	PYR	Pyrrrolidonyl- β- naphthylamide	0.1%		
18	ADON	Adonitol	1.0%	Utilization of the carbohydrate substrate produces acidic products which lower the pH and change the indicator.	1, 3
18	IND	Tryptophane	0.4%	Utilization of tryptophane results in the formation of indole which is detected with RapID Spot Indole Reagent.	1-3

- While maintaining a level, horizontal position (best achieved by using the bench top against the reaction cavity bottoms), slowly tilt the panel forward toward the reaction cavities until the inoculum flows along the baffles into the reaction cavities (see below). This should evacuate all of the inoculum from the rear portion of the panel.

Note: If the panel is tilted too quickly, air may be trapped at the test cavity junction, restricting fluid movement.

- Return the panel to a level position. If necessary, gently tap the panel on the bench top to remove any air trapped in the cavities.

Notes:

- Examine the test cavities, which should appear bubble-free and uniformly filled. Slight irregularities in test cavity fills are acceptable and will not affect test performance. If the panel is grossly misfilled, a new panel should be inoculated and the misfilled panel discarded.
- Complete the inoculation of each panel receiving inoculation fluid before inoculating additional panels.
- Do not allow the inoculum to rest in the back portion of the panel for prolonged periods without completing the procedure.

Incubation of RapID ONE Panels:

Incubate inoculated panels at 35-37°C in a non-CO₂ incubator for 4 hours. For ease of handling, panels may be incubated in the chipboard incubation trays provided with the kit.

RapID ONE Panel Test Location

Cavity #	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18
Test Code	URE	ADH	ODC	LDC	TET	LIP	KSF	SBL	GUR	ONPG	βGLU	βXYL	NAG	MAL	PRO	GGT	PYR	ADON
															RapID ONE Reagent			IND

Table 2. Interpretation of RapID ONE System Tests*

Cavity #	Test Code	Reagent	Reaction		Comments
			Positive	Negative	
1	URE	None	Red or violet	Yellow or orange	Only the development of a distinct red or violet color should be scored as positive.
2	ADH	None	Bright purple or blue	Yellow, gray, straw, or yellow-green	Only the development of a distinct, bright purple or blue color should be scored as positive. Shades of yellow, gray, straw, brown, or green should be scored as negative.
3	ODC				
4	LDC				
5	TET	None	Yellow	Red or orange	Only the development of a distinct yellow color throughout the well should be scored as positive. Layers of yellow color should not be interpreted as positive. Cavities may be mixed with an applicator stick to aid reading.
6	LIP				
7	KSF				
8	SBL	None	Yellow	Clear or tan	Only the development of a distinct yellow color should be scored as positive. Very faint yellow tints or equivocal colors should be scored as negative.
9	GUR				
10	ONPG				
11	βGLU				
12	βXYL				
13	NAG	None	Red	Yellow or orange	Only the development of a distinct red color should be scored as positive. Shades of orange should be scored as negative.
14	MAL				
15	PRO				
16	GGT	RapID ONE Reagent	Violet, purple, red, or dark pink	Clear, yellow, orange, or very pale pink	Any development of a violet, purple, or distinct red or dark pink color should be scored as positive. Pale orange or faint pinks should be scored as negative.
17	PYR				
18	ADON	None	Yellow or very light orange	Red or dark red-orange	Any development of a yellow or yellow-orange color throughout the cavity should be scored as positive.
19	IND	RapID Spot Indole Reagent	Brown, black, or purple	Orange or red	Any development of a brown, black, or purple color should be scored as positive.

*NOTE: Panels should be read by looking down through the reaction cavities against a white background.

Table 3. Quality Control Chart for RapID ONE Panels

Organism	URE	ADH	ODC	LDC	TET	LIP	KSF	SBL	GUR	ONPG	βGLU	βXYL	NAG	MAL	PRO	GGT	PYR	ADON	IND	OXI
<i>Proteus vulgaris</i> ATCC® 6380	+	V	-	-	+	-	-	-	-	-	V	-	-	-	-	+	-	-	+	-
<i>Escherichia coli</i> ® ATCC® 25922	-	-	+	(+)	-	-	-	(+)	+	+	-	-	-	-	-	-	-	-	+	-
<i>Pseudomonas aeruginosa</i> ATCC® 27853	V	+	V	V	V	+	-	-	-	-	-	-	-	V	+	+	V	-	-	+
<i>Klebsiella aerogenes</i> ® ATCC® 13048	-	-	+	+	-	-	+	+	-	+	+	+	+	+	V	V	+	+	-	-

+, positive; -, negative; V, variable; (+), usually positive

® Key indicator strains demonstrate acceptable performance of the most labile substrate in the system and reactivity in a significant number of wells, according to Clinical and Laboratory Standards Institute recommendations for streamlined quality control.¹⁴

14. QUALITY CONTROL

All lot numbers of the RapID ONE System have been tested using the following quality control organisms (Table 3) and have been found to be acceptable. Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported. Table 3. lists expected results for the selected battery of test organisms.

Notes:

- RapID reagent quality control is accomplished by obtaining the expected reactions for tests requiring the addition of the reagents (cavities 15-18).
- Organisms which have been repeatedly transferred on agar media for prolonged periods may provide aberrant results.
- Quality control strains should be stored frozen or lyophilized. Prior to use, quality control strains should be transferred 2-3 times from storage on an agar medium that is recommended for use with the RapID ONE System.

Formulations, additives, and ingredients of culture media vary from manufacturer to manufacturer and may vary from batch to batch. As a result, culture media may influence constitutive enzymatic activity of designated quality control strains. If quality control strain results differ from the patterns indicated, a subculture onto medium from a different batch or from another manufacturer will often resolve quality control discrepancies.

15. LIMITATIONS

- The use of RapID ONE System and the interpretation of results requires the knowledge of a competent microbiologist, familiar with laboratory procedures, who is trained in general microbiological methods and who judiciously makes use of training, experience, specimen information, and other pertinent procedures before reporting the identification obtained using this system.
- Specimen source, oxidase reaction, Gram stain characteristics, and growth on selective agars should be considered when using the RapID ONE System.

- The RapID ONE System must be used with pure cultures of test organisms. The use of mixed microbial populations or direct testing of clinical material without culture will result in aberrant results.
- The RapID ONE System is designed for use with the taxa listed in the RapID ONE Differential Chart. The use of organisms not specifically listed may lead to misidentifications.
- Expected values listed for RapID ONE System tests may differ from conventional test results or previously reported information.
- The accuracy of the RapID ONE System is based upon the statistical use of a multiplicity of specially designed tests and an exclusive, proprietary database. The use of any single test found in the RapID ONE System to establish the identification of a test isolate is subject to the error inherent in that test alone.

16. PERFORMANCE CHARACTERISTICS

The performance characteristics of RapID ONE System have been established by laboratory testing of reference and stock cultures at Remel, and by clinical evaluations using fresh clinical and stock isolates.¹⁰⁻¹³

17. BIBLIOGRAPHY

- Ewing, W.H. 1986. Edwards and Ewing's Identification of Enterobacteriaceae. 4th ed. Elsevier Science Publishing Company, Inc., New York, NY.
- Balows, A., W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy. 1991. Manual of Clinical Microbiology. 5th ed. ASM, Washington, D.C.
- MacFaddin, J.F. 2000. Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Lippincott Williams & Wilkins, Philadelphia, PA.
- Holt J.G., N.R. Krieg, P.H.A. Sneath, J.T. Staley, and S.T. Williams. 2000. Bergey's Manual of Determinative Bacteriology. 9th ed. Lippincott Williams and Wilkins. Philadelphia, PA.
- Eriqez, L.A., N.E. Hodinka, and K.A. Hornsby. 1993. Abstract C-294. Abstracts of the 93rd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.

- Kilian, M. and P. Bulow. 1976. Acta Pathol. Microbiol. Scand. B. 84:245-251.
- Humble, W.M., A. King, and I. Phillips. 1977. J. Clin. Pathol. 30:275-277.
- K.C. Carroll, M.A. Pfaller, M.L. Landry, A.J. McAdam, R. Patel, S.S. Richter, D.W. Warnock. 2019. Manual of Clinical Microbiology. 12th Edition. ASM Press.
- Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2007. Bailey and Scott's Diagnostic Microbiology. 12th ed. Mosby Elsevier, St. Louis, MO.
- Jerris, R., A. Kabant, S. Peroulas, T. Lee, K.A. Hornsby, and D. Lockhart. 1993. Abstract C-304. Abstracts of the 93rd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.

- Kitch, T., M.R. Jacobs, and P.C. Appelbaum. 1993. Abstract C-303. Abstracts of the 93rd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Schreckenberger, P., M. Montero, and N. Heldt. 1993. Abstract C-309. Abstracts of the 93rd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Eriqez, L.A., A.P. Jones, and N.E. Hodinka. 1992. Abstract C-5. Abstracts of the 92nd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Clinical and Laboratory Standards Institute (CLSI). 2008. Quality Control for Commercial Microbial Identification Systems; Approved Guideline. M50-A. CLSI, Wayne, PA.

18. PACKAGING

REF R8311006 RapID ONE System..... 20 Tests/Kit

19. SYMBOL LEGEND

	Catalogue Number
	In Vitro Diagnostic Medical Device
	Consult Instructions for Use (IFU)
	Temperature Limitations (Storage temp.)
	Contains sufficient for <N> tests
	Do not use if package is damaged
	Do not re-use
	Batch Code (Lot Number)
	Use By (Expiration Date)
	Importer
	Unique Device Identifier
	Authorized representative in the European Community
	UK Conformity Assessed
	European Conformity Assessment
	Manufacturer

RapID™ and ERIC™ are trademarks of Thermo Fisher Scientific and its subsidiaries.

ATCC® is a registered trademark of American Type Culture Collection.

Remel Inc., 12076 Santa Fe Trail Drive, Lenexa, KS 66215, USA

www.thermofisher.com/microbiology
Tel: (800) 255-6730 • International: (913) 888-0939

www.oxid.com/IFU
Europe +800 135 79 135 • US 1 855 2360 190
CA 1 855 805 8539 • ROW +31 20 794 7071

Version	Date of modifications introduced
IFU8311006	July 2025 Culture name changes. The strains themselves have not changed.

Table 4 - RapID ONE Differential Chart

Organism	URE	ADH	ODC	LDC	TET	LIP	KSF	SBL	GUR	ONPG	βGLU	βXYL	NAG	MAL	PRO	GGT	PYR	ADON	IND	OXI
<i>Acinetobacter calcoaceticus</i>	5	29	11	7	2	94	0	1	0	0	0	2	21	3	5	2	0	0	0	0
<i>Burkholderia cepacia</i> ^a	14	12	31	84	3	95	2	0	0	14	9	0	19	16	71	95	0	0	0	62
<i>Cedecea davisae</i> (EG-15)	2	27	51	0	0	71	88	0	0	51	98	96	99	92	1	82	0	0	0	0
<i>Cedecea lapagei</i>	0	38	0	7	0	90	95	0	0	90	92	0	99	98	1	98	0	0	0	0
<i>Cedecea neteri</i>	2	46	2	5	0	93	98	95	0	98	98	0	99	98	0	99	0	0	0	0
<i>Cedecea sp. 3</i>	0	90	0	0	0	95	90	0	0	91	97	0	98	0	0	97	0	0	0	0
<i>Cedecea sp. 5</i>	0	91	35	2	0	81	90	98	0	90	98	0	98	0	0	95	0	0	0	0
<i>Citrobacter amaloniticus</i>	11	11	96	3	98	0	95	96	0	95	83	1	9	2	2	92	98	2	98	0
<i>Citrobacter freundii</i>	11	18	24	0	88	0	94	92	0	92	5	1	0	5	3	88	96	1	1	0
<i>Citrobacter koseri</i> ^b	5	20	20	2	0	0	95	95	0	95	92	0	1	90	4	95	96	80	96	0
<i>Cronobacter sakazakii</i> ^c	0	89	92	3	0	0	98	1	0	96	84	96	98	19	0	95	38	0	20	0
<i>Edwardsiella hoshinae</i>	0	0	79	98	11	0	8	0	0	0	0	2	90	69	99	5	0	0	78	0
<i>Edwardsiella tarda</i>	0	0	83	89	90	0	2	0	0	3	0	2	96	0	98	4	0	0	92	0
<i>Enterobacter asburiae</i> (EG-17)	7	6	89	0	0	0	98	94	0	98	94	30	97	2	0	82	9	0	0	0
<i>Enterobacter cancerogenus</i> ^d	0	91	95	1	0	0	98	0	0	98	84	2	98	98	0	98	96	0	0	0
<i>Enterobacter cloacae</i>	9	92	92	1	0	2	95	96	0	98	82	92	94	76	1	94	21	17	0	0
<i>Escherichia coli</i>	1	15	75	84	2	0	1	90	92	95	1	0	0	1	3	44	1	2	98	0
<i>Escherichia ferqusonii</i>	0	5	88	96	0	0	96	0	0	86	2	0	0	14	2	90	98	82	94	0
<i>Escherichia hermannii</i> (EG-11)	0	2	87	2	98	0	98	3	0	85	67	0	1	0	4	96	92	0	96	0
<i>Ewingella americana</i> ^e	0	0	0	0	0	0	56	0	0	92	44	0	94	0	0	21	99	0	0	0
<i>Hafnia alvei</i>	0	13	92	96	5	1	87	0	0	31	7	0	77	20	97	96	11	0	2	0
<i>Klebsiella aerogenes</i> ^f	4	7	97	97	1	1	96	97	0	98	98	98	97	92	1	93	98	88	0	0
<i>Klebsiella oxytoca</i>	93	12	9	97	9	2	98	96	0	95	98	84	0	95	0	96	82	76	99	0
<i>Klebsiella pneumoniae subsp. ozaenae</i>	6	12	2	32	0	0	84	48	0	97	94	17	0	1	1	88	98	89	0	0
<i>Klebsiella pneumoniae</i>	96	13	9	97	4	1	96	96	0	97	98	94	0	91	0	98	92	88	0	0
<i>Klebsiella pneumoniae subsp. rhinoscleromatis</i>	0	0	0	0	0	0	5	27	0	0	14	0	0	29	0	2	99	4	0	0
<i>Kluyvera ascorbata</i>	0	0	98	90	0	0	90	40	0	98	96	96	0	93	0	91	8	0	89	0
<i>Kluyvera cryocrescens</i>	0	0	98	11	0	0	91	38	0	98	93	96	0	79	0	7	2	0	90	0
<i>Kluyvera intermedia</i> ^g	0	0	86	0	0	0	98	98	0	97	98	98	0	98	0	92	98	0	0	0
<i>Leclercia adecarboxylata</i> (EG-40)	20	0	0	0	0	0	42	0	0	95	92	94	96	86	0	5	93	91	94	0
<i>Lelliottia amniqena</i> ^h	0	2	78	0	9	0	90	21	0	98	97	98	96	93	0	28	96	0	0	0
<i>Leminorella grimontii</i> (EG-57)	0	0	7	6	88	0	5	0	0	0	0	0	99	0	0	98	0	0	0	0
<i>Leminorella richardii</i>	0	0	5	2	91	1	9	0	0	0	0	0	0	0	0	99	0	0	0	0
<i>Moellerella wisconsensis</i> (EG-46)	0	0	0	0	89	0	2	0	0	93	8	0	0	0	0	5	0	97	0	0
<i>Morganella morganii</i>	98	18	91	4	98	0	2	0	0	0	0	0	0	2	0	96	0	0	98	0
<i>Pantoea agglomerans</i> ⁱ	2	2	0	0	2	0	91	32	0	93	71	51	21	61	0	94	77	7	11	0
<i>Pluralibacter gergoviae</i> ^j	74	9	96	84	0	0	95	0	0	88	91	0	11	92	0	98	9	0	0	0
<i>Proteus mirabilis</i>	98	13	91	1	98	5	1	0	0	0	0	0	0	4	7	96	0	0	1	0
<i>Proteus penneri</i>	99	0	0	0	79	0	0	0	0	0	0	0	0	8	9	49	0	0	0	0
<i>Proteus vulgaris</i> Group 2	98	18	6	0	95	5	0	0	0	0	98	0	0	2	2	94	0	0	96	0
<i>Proteus vulgaris</i> Group 3	98	12	0	0	97	4	3	0	0	0	0	0	0	2	1	98	0	0	96	0
<i>Providencia alcalifaciens</i>	4	4	1	0	94	0	1	0	0	0	0	0	1	2	4	98	0	86	98	0
<i>Providencia rettgeri</i>	98	11	2	0	92	0	73	2	0	0	30	0	3	4	0	96	0	90	96	0
<i>Providencia rustigianii</i>	0	2	0	0	99	0	2	0	0	0	0	0	0	0	0	95	0	0	98	0
<i>Providencia stuartii</i>	22	3	2	1	96	0	2	0	0	0	1	0	97	0	1	98	2	2	95	0
<i>Pseudoescherichia vulneris</i> ^k	0	9	0	91	0	0	95	0	0	96	90	94	0	93	0	57	93	4	0	0
<i>Pseudomonas luteola</i> ^l	8	77	3	4	2	2	0	0	0	92	95	0	0	32	99	78	98	0	0	8
<i>Pseudomonas oryzihabitans</i> ^m	5	2	8	4	0	90	5	4	0	2	3	0	0	88	99	21	98	0	0	0
<i>Rahnella aquatilis</i> ⁿ	0	0	0	0	0	0	88	94	0	96	97	98	0	98	92	4	98	0	0	0
<i>Raoultella ornitholytica</i> ^k	98	8	99	98	0	0	98	98	1	98	98	98	98	98	0	98	98	93	98	0
<i>Raoultella planticola</i> (EG-47) ^o	98	11	2	98	0	0	98	90	0	99	98	15	90	98	0	98	98	90	11	0
Salmonella l - includes the following Salmonella choleraesuis serotypes:	1	54	92	94	91	0	2	91	20	1	1	0	0	1	8	97	0	0	1	0
Choleraesuis	0	60	98	94	56	0	0	90	11	0	0	0	0	0	7	95	0	0	0	0

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- Ewing, W.H. 1986. Edwards and Ewing's Identification of Enterobacteriaceae. 4th ed. Elsevier Science Publishing Company, Inc., New York, NY.
- Balows, A., W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy. 1991. Manual of Clinical Microbiology. 5th ed. ASM, Washington, D.C.
- MacFaddin, J.F. 2000. Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Lippincott Williams & Wilkins, Philadelphia, PA.
- Holt J.G., N.R. Krieg, P.H.A. Sneath, J.T. Staley, and S.T. Williams. 2000. Bergey's Manual of Determinative Bacteriology. 9th ed. Lippincott Williams and Wilkins. Philadelphia, PA.
- Eriquerz, L.A., N.E. Hodinka, and K.A. Hornsby. 1993. Abstract C-294. Abstracts of the 93rd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Kilian, M. and P. Bulow. 1976. Acta Pathol. Microbiol. Scand. B. 84:245-251.
- Humble, W.M., A. King, and I. Phillips. 1977. J. Clin. Pathol. 30:275-277.
- K.C. Carroll, M.A. Pfaller, M.L. Landry, A.J. McAdam, R. Patel, S.S. Richter, D.W. Warnock. 2019. Manual of Clinical Microbiology. 12th Edition. ASM Press.
- Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2007. Bailey and Scott's Diagnostic Microbiology. 12th ed. Mosby Elsevier, St. Louis, MO.
- Jerris, R, A. Kabant, S. Peroulas, T. Lee, K.A. Hornsby, and D. Lockhart. 1993. Abstract C-304. Abstracts of the 93rd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Kitch, T., M.R. Jacobs, and P.C. Appelbaum. 1993. Abstract C-303. Abstracts of the 93rd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Schreckenberger, P., M. Montero, and N. Heldt. 1993. Abstract C-309. Abstracts of the 93rd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Eriquerz, L.A., A.P. Jones, and N.E. Hodinka. 1992. Abstract C-5. Abstracts of the 92nd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Clinical and Laboratory Standards Institute (CLSI). 2008. Quality Control for Commercial Microbial Identification Systems; Approved Guideline. M50-A. CLSI, Wayne, PA.

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Remel Inc., 12076 Santa Fe Trail Drive, Lenexa, KS 66215, САЩ

www.thermofisher.com/microbiology

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- Jerris, R., A. Kabant, S. Peroulas, T. Lee, K.A. Hornsby, and D. Lockhart. 1993. Abstract C-304. Abstracts of the 93rd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Kitch, T., M.R. Jacobs, and P.C. Appelbaum. 1993. Abstract C-303. Abstracts of the 93rd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Schreckenberger, P., M. Montero, and N. Heldt. 1993. Abstract C-309. Abstracts of the 93rd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Eriquez, L.A., A.P. Jones, and N.E. Hodinka. 1992. Abstract C-5. Abstracts of the 92nd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Clinical and Laboratory Standards Institute (CLSI). 2008. Quality Control for Commercial Microbial Identification Systems; Approved Guideline. M50-A. CLSI, Wayne, PA.

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17. BIBLIOGRAFIA

- Ewing, W.H. 1986. Edwards and Ewing's Identification of Enterobacteriaceae. 4th ed. Elsevier Science Publishing Company, Inc., New York, NY.
- Balows, A., W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy, 1991. Manual of Clinical Microbiology. 5th ed. ASM, Washington, D.C.
- MacFaddin, J.F. 2000. Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Lippincott Williams & Wilkins, Philadelphia, PA.

- Holt J.G., N.R. Krieg, P.H.A. Sneath, J.T. Staley, and S.T. Williams. 2000. Bergey's Manual of Determinative Bacteriology. 9th ed. Lippincott Williams and Wilkins. Philadelphia, PA.
- Eriquez, L.A., N.E. Hodinka, and K.A. Hornsby. 1993. Abstract C-294. Abstracts of the 93rd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Kilian, M. and P. Bulow. 1976. Acta Pathol. Microbiol. Scand. B. 84:245-251.
- Humble, W.M., A. King, and I. Phillips. 1977. J. Clin. Pathol. 30:275-277.
- K.C. Carroll, M.A. Pfaller, M.L. Landry, A.J. McAdam, R. Patel, S.S. Richter, D.W. Warnock. 2019. Manual of Clinical Microbiology. 12th edition. ASM Press.
- Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2007. Bailey and Scott's Diagnostic Microbiology. 12th ed. Mosby Elsevier, St. Louis, MO.
- Jerris, R., A. Kabant, S. Peroulas, T. Lee, K.A. Hornsby, and D. Lockhart. 1993. Abstract C-304. Abstracts of the 93rd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Kitch, T., M.R. Jacobs, and P.C. Appelbaum. 1993. Abstract C-303. Abstracts of the 93rd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Schreckenberger, P., M. Montero, and N. Heldt. 1993. Abstract C-309. Abstracts of the 93rd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Eriquez, L.A., A.P. Jones, and N.E. Hodinka. 1992. Abstract C-5. Abstracts of the 92nd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Clinical and Laboratory Standards Institute (CLSI). 2008. Quality Control for Commercial Microbial Identification Systems; Approved Guideline. M50-A. CLSI, Wayne, PA.

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REF R8311006 Sistema RapID ONE..... 20 test/kit

19. LEGENDA DEI SIMBOLI

Table with 2 columns: Symbol and Description. Symbols include REF, IVD, book icon, thermometer icon, Sigma N, crossed-out person icon, crossed-out person icon, LOT, lightbulb icon, globe icon, UDI, EC REP, UK CA, CE, factory icon.

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