MICROLAB

LISTERIA

MONOCYTOGENES

1. INTRODUCTION AND SCOPE

Microlab *Listeria monocytogenes* is a qualitative test for the detection of *Listeria monocytogenes* in food and surfaces that integrates all the steps that would be carried out in a microbiological laboratory. ***L. innocua* or *L. marthii* could show positive results in foods at levels above 85 CFU/25 g.**

Microlab *Listeria monocytogenes* can be used in the food industry or in a laboratory. This test does not require qualified personnel or specific equipment.

Microlab *Listeria monocytogenes* allows to determine the absence of *Listeria monocytogenes* in 25 g of food, in accordance with the current legislation (Commission Regulation N° 2073/2005). Microlab *Listeria monocytogenes* integrates the steps of enrichment, detection and inactivation in a single device.

- » Enrichment step: *Listeria monocytogenes* grows to detectable levels in a specific culture medium (LOVETT media + selective supplements; ZE/MILIS/BROTH) .
- » Detection step: bacteria are detected by an immunochromatographic test with specific antibodies against *L. monocytogenes* (ZE/MILIS/STRIP).
- » Inactivation step: an inactivation solution eliminates any microorganism that may have grown in the culture medium (ZE/MILIS/INAC).

2. VALIDATION STATMENT

This method has been certified (AOAC-PTM 062502) for the use with cheddar cheese, ice cream, smoked salmon, deli ham and deli turkey using a 25 g test portion. Surface testing is not included in the AOAC-PTM certification claim.

3. KIT COMPONENTS

- » 6 Microlab Listeria monocytogenes devices.
- » Kit instructions.
- » Product certificate.

4. MATERIAL AND EQUIPMENT NOT INCLUDED

- » Incubator (37°C).
- » Weighing scale.

5. STABILITY AND STORAGE

The device should be stored at 4-10°C. kit shelf life is 12 months from production. See the kit expiry date on the package.

6. PRECAUTIONS

1. Storage conditions higher than 4-10°C may adversely affect performance of the test.
2. If a positive result is obtained, it should be confirmed. Sample extraction must be performed through the septum located at the bottom of the device, ALWAYS BEFORE THE INACTIVATION. Confirm positive results following ISO 11290-1 or other approved methods. Contact your distributor for more information.
3. Microlab Listeria monocytogenes is an in vitro diagnostic test. In analysis involving legal processes, the results should be re-evaluated with an official reference method. Zeulab does not assume any legal responsibility.

7. SAFETY INSTRUCTIONS

A safety data sheet is available from your local distributor or ZEULAB. Use this kit following good laboratory practices. Additional precautions and qualified personnel are not required.

8. TEST PROCEDURE

See point 9 (page 12) to identify the different parts of the device marked with letters. **Let the device reach 37°C prior to be used.**

1. Turn the green ring (F) anticlockwise holding the enrichment chamber (I) to open the device. Remove the protective band.

Sample preparation:

2.a. Food Analysis: Weigh 25 g of a representative and homogenised sample on the top of the enrichment chamber (I).

Surface analysis 2.b (not certified by AOAC-RI):

Cut the aluminum foil with a previously sterilized tool (scalpel, cutter...), hydrate the sponge with the culture medium, sample the surface and place the sponge in the enrichment chamber (I).

Surface sampling should be performed wearing gloves.

Surface analysis can be performed with premoistened sponges.

Enrichment:

3. Align the arrow (E) with the trigger (H) to assemble the device. While pressing the trigger (H), turn the ring (F) clockwise for two rounds.

4. Turn clockwise the upper cylinder while holding the enrichment chamber (I) until the aluminium foil tears.


5. Continue turning the ring (F) clockwise, without pressing the trigger (H). After 1,5 turns, the ring block (G) should be at the left of the trigger (H), causing the final click of the rounds.

6. Shake the device for 5 seconds to obtain a homogeneous mixture between the sample and the enrichment media (included).

7. Incubate the device at 37 °C ($\pm 1^\circ\text{C}$) for 48 hours (± 1 hour).

Detection:

8. Shake the device for 5 seconds to homogenize its contents.
9. Remove the rotator's security seal (K).
10. Turn the rotator (A) anticlockwise until the rotator's green marker (L) reaches position 1 (J).
11. Turn the device slowly upside down and wait for 10 seconds until the chamber marked as 1 is filled. In case it is not properly filled, shake it gently for 3 seconds in a vertical position.
12. Turn the device to the upright position. Continue turning the rotator (A) anticlockwise until the rotator's marker (L) reaches position 2 (C).
13. Wait for 15 minutes and read the results. Readings after 20 minutes will not be valid.

Positive: Two lines, red and blue, are observed in the white central zone.

Negative: Only a blue line is observed in the white central zone of the strip test.

Invalid: When no blue line is shown.

Inactivation:

14. Continue turning the rotator (A) until the rotator's marker (L) reaches the end (B).
15. Shake the device to homogenize and turn it upside down several times. The inactivation solution should reach all the chambers.
16. Leave the device in upright position for 30 minutes.
17. Turn the device upside down and leave 30 more minutes to complete the inactivation.
18. Discard the device as a solid urban waste.

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