MICROLAB **SALMONELLA**

PERFORMANCE TESTED

AOAC
RESEARCH INSTITUTE

LICENSE NUMBER 052401

Performance Tested MethodsSM number #052401.

1. OBJECTIVE

MicroLab Salmonella is a qualitative test for the detection of *Salmonella enterica* serogroups A-E in food and surfaces that integrates all the steps that would be carried out in a microbiological laboratory.

MicroLab Salmonella can be used in the food industry or in a laboratory. This test does not require qualified personnel or specific equipment.

MicroLab Salmonella allows to determine the absence of *Salmonella* in 25 g of food, in accordance with the current legislation (Commission Regulation (EC) No 2073/2005 and subsequent modifications).

NOTE: MicroLab Salmonella is an in vitro diagnostic test. In analysis involving legal processes, the results should be re-evaluated with an official reference method. GTZ MicroLab Detect S.L. does not assume any legal responsibility.

MicroLab Salmonella integrates the steps of enrichment, detection and inactivation in a single device.

- » Enrichment step: *Salmonella* grows to detectable levels in a ready-to-use culture medium.
- » Detection step: the bacteria are detected by an immunochromatographic test with specific antibodies against *Salmonella enterica* serogroups A-E.
- » Inactivation step: an inactivation solution eliminates any microorganism that may have grown in the culture medium.

2. VALIDATION STATMENT

This method has been certified (AOAC-PTM) for the use with raw ground beef, raw turkey (thermal processed marinated), ready-to-eat fresh cheese, pasteurized liquid eggs and ready-to-eat deli ham using a 25 g test portion.

3. SAFETY PRECAUTIONS

A safety data sheet is available from your local distributor or ZEULAB. Use this kit following good laboratory practices. Additional precautions and qualified personnel are not required.

4. KIT COMPONENTS

- » 6 MicroLab Salmonella devices.
- » Kit instructions.
- » Product certificate.

5. MATERIAL AND EQUIPMENT NOT INCLUDED

- » Incubator (37°C).
- » Weighing scale.

6. STABILITY AND STORAGE

The kit components should be stored at 4-25°C. The kit shelf life is 12 months from production. See the kit expiry date on the package.

7. TEST PROCEDURE

See point 8 (page 12) to identify the different parts of the device marked with letters. Let the device reach 37°C prior to be used.

1. Turn the green ring (F) anticlockwise holding the enrichment chamber (I) to open the device. Remove the protective band.


2. Sample preparation:

Surface analysis (2.a) is not certified by AOAC-PTM.


2.a. Surface analysis: Cut the aluminum foil with a previously sterilized tool (scalpel, cutter...), hydrate the sponge with the culture medium, sample the surface and place the sponge in the enrichment chamber (I).

Surface sampling should be performed wearing gloves.


Surface análisis can be performed with premoistened sponges.

2.b. Food Analysis: Weigh 25 g of a representative and homogenized sample on the top of the enrichment chamber (I).

Enrichment:


3. Align the arrow (E) with the trigger (H) to assemble the device. While pressing the trigger (H), turn the ring (F) clockwise for two rounds.

4. Turn clockwise the upper cylinder while holding the enrichment chamber (I) until the aluminium foil tears.

5. Continue turning the ring (F) clockwise, without pressing the trigger (H). After 1,5 turns, the ring block (G) should be at the left of the trigger (H), causing the final click of the rounds.

6. Shake the device for 5 seconds to obtain a homogeneous mixture between the sample and the enrichment media (included).

7. Incubate the device at 37 °C ($\pm 1^{\circ}\text{C}$) for 24 hours (± 1 hour).

Detection:

8. Shake the device for 5 seconds to homogenize its contents.


9. Remove the rotator's security seal (K).

10. Turn the rotator (A) anticlockwise until the rotator's green marker (L) reaches position 1 (J).

11. Turn the device slowly upside down and wait for 10 seconds until the chamber marked as 1 is filled. In case it is not properly filled, shake it gently for 3 seconds in a vertical position.

12. Turn the device to the upright position. Continue turning the rotator (A) anticlockwise until the rotator's marker (L) reaches position 2 (C).

13. Wait for 15 minutes and read the results. Readings after 20 minutes will not be valid.


Positive: Two lines, red and blue, are observed in the white central zone.

Negative: Only a blue line is observed in the white central zone of the strip test.

Invalid: When no blue line is shown.

Inactivation:

NOTE: In case the confirmation of a result is necessary, the extraction of the sample should be done ALWAYS BEFORE THE INACTIVATION. Perform the confirmation following the ISO 6579-1 or other approved methods. Contact your distributor for more information.

14. Continue turning the rotator (A) until the rotator's marker (L) reaches the end (B).

15. Shake the device to homogenize and turn it upside down several times. The inactivation solution should reach all the chambers.

16. Leave the device in upright position for 30 minutes.

17. Turn the device upside down and leave 30 more minutes to complete the inactivation.

18. Discard the device as a solid urban waste.

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