



## Heathrow Scientific Microtube Storage Boxes



- Made from polypropylene
- Durable three point hinges and a clasp that locks for secure lid closure
- Featuring moulded grid lines on the lid
- Imprinted coordinates on the base, and imprinted and raised coordinates on the box
- User friendly and stackable
- Black, light tight version available for light sensitive samples
- Autoclavable
- Suitable for storage down to  $-80^{\circ}\text{C}$

Supplied in packs of 5

Code	Alt Ref	Capacity	Colour	Dims, w x d x h, mm
<b>RAC0470</b>	HS120032	50 x 1.5 - 2.0mL tubes	Natural	141 x 92 x 56
<b>RAC0472</b>	HS120033	50 x 1.5 - 2.0mL tubes	Assorted (blue, green, purple, yellow, orange)	141 x 92 x 56
<b>RAC0484</b>	HS120043	100 x 1.5 - 2.0mL tubes	Black	141 x 151 x 57
<b>RAC0486</b>	HS120201	100 x 1.5 - 2.0mL tubes	Blue	141 x 151 x 57
<b>RAC0487</b>	HS120202	100 x 1.5 - 2.0mL tubes	Green	141 x 151 x 57
<b>RAC0488</b>	HS120203	100 x 1.5 - 2.0mL tubes	Yellow	141 x 151 x 57
<b>RAC0489</b>	HS120204	100 x 1.5 - 2.0mL tubes	Orange	141 x 151 x 57
<b>RAC0480</b>	HS120042	100 x 1.5 - 2.0mL tubes	Natural	141 x 151 x 57
<b>RAC0482</b>	HS120044	100 x 1.5 - 2.0mL tubes	Assorted (blue, green, purple, yellow, orange)	141 x 151 x 57



## Bravo Storage/Cryo Boxes

- Water repellent fibreboard boxes in a choice of colours
- Two heights: 75 or 50mm
- Available with plastic coating for additional waterproofing
- Does not include dividers
- Dimensions, w x d: 133 x 133mm



Code	Colour	Height, mm	Description	Pack
<b>FRE2130</b>	White	50	Standard	12
<b>FRE2132</b>	White	75	Standard	12
<b>FRE2246</b>	Red	50	Standard	Each
<b>FRE2247</b>	Red	50	Plastic coated	Each
<b>FRE2216</b>	Red	75	Standard	Each
<b>FRE2217</b>	Red	75	Plastic coated	Each
<b>FRE2240</b>	Yellow	50	Standard	Each
<b>FRE2241</b>	Yellow	50	Plastic coated	Each
<b>FRE2210</b>	Yellow	75	Standard	Each
<b>FRE2211</b>	Yellow	75	Plastic coated	Each
<b>FRE2242</b>	Blue	50	Standard	Each
<b>FRE2243</b>	Blue	50	Plastic coated	Each
<b>FRE2212</b>	Blue	75	Standard	Each
<b>FRE2213</b>	Blue	75	Plastic coated	Each
<b>FRE2244</b>	Green	50	Standard	Each
<b>FRE2245</b>	Green	50	Plastic coated	Each
<b>FRE2214</b>	Green	75	Standard	Each
<b>FRE2215</b>	Green	75	Plastic coated	Each