

eShield® Sterile Cover Smartphone Sizing Chart

NOTE: If you have a protective case on your smartphone, you may need the larger size cover to protect your device.

SKU: EC2100



Overall Dimensions: 9" x 14" (23 cm x 36 cm)
Functional Dimensions: 9" x 4.5" (23 cm x 11 cm)

Manufacturer	Smartphone Model	Height (in.)	Width (in.)	Depth (in.)
Apple	iPhone 16	5.81	2.82	0.31
Apple	iPhone 16 Plus	6.34	3.06	0.31
Apple	iPhone 16 Pro	5.89	2.81	0.32
Apple	iPhone 16 Pro Max	6.42	3.06	0.32
Apple	iPhone 15	5.81	2.82	0.31
Apple	iPhone 15 Plus	6.33	3.06	0.31
Apple	iPhone 15 Pro	5.77	2.78	0.33
Apple	iPhone 15 Pro Max	6.30	3.02	0.33
Apple	iPhone 14	5.78	2.82	0.31
Apple	iPhone 14 Pro	5.81	2.81	0.31
Apple	iPhone 13	5.78	2.82	0.30
Apple	iPhone 13 Pro	5.78	2.82	0.30
Apple	iPhone 12 Mini	5.18	2.53	0.29
Apple	iPhone 12	5.78	2.82	0.29
Apple	iPhone 12 Pro	5.78	2.82	0.29
Apple	iPhone 12 Pro Max	6.33	3.07	0.29
Apple	iPhone 11	5.94	2.98	0.33
Apple	iPhone 11 Pro	5.66	2.81	0.32
Apple	iPhone 11 Pro Max	6.22	3.06	0.32
Apple	iPhone 5	4.87	2.3	0.30
Apple	iPhone 5c	4.89	2.33	0.35
Apple	iPhone 5s	4.87	2.3	0.30
Apple	iPhone 6	5.43	2.63	0.27
Apple	iPhone 6s	5.44	2.64	0.28
Apple	iPhone 6s Plus	6.22	3.06	0.29
Apple	iPhone 7	5.44	2.64	0.28
Apple	iPhone 7 Plus	6.22	3.06	0.29
Apple	iPhone 8	5.44	2.64	0.29
Apple	iPhone 8 Plus	6.23	3.07	0.30
Apple	iPhone SE	4.87	2.3	0.30
Apple	iPhone X	5.65	2.79	0.30
Apple	iPhone XR	5.94	2.98	0.33
Apple	iPhone XS	5.65	2.79	0.30
Apple	iPhone XS Max	6.20	3.04	0.30
Google	Pixel 9 Pro	6.40	3.00	0.33
Google	Pixel 8 Pro	6.40	3.01	0.35
Google	Pixel 5a	6.10	2.90	0.30
Google	Pixel 4a	5.67	2.73	0.32
Google	Pixel (2016)	5.66	2.74	0.33

Manufacturer	Smartphone Model	Height (in.)	Width (in.)	Depth (in.)
Motorola	Razr V3	3.86	2.10	0.54
Motorola	Moto Z4	6.22	2.95	0.29
Motorola	Moto E4	5.70	2.85	0.37
Motorola	Edge 40	6.24	2.83	0.30
Motorola	Edge 40 Pro	6.35	2.95	0.34
Samsung	Galaxy Z Flip 4 (open/closed)	6.50/3.35	2.83	0.27/0.67
Samsung	Galaxy Z Flip 5 (open/closed)	6.50/3.35	2.83	0.27/0.60
Samsung	Galaxy Z Flip 7 (open/closed)	6.53/3.36	2.96	0.25/0.54
Samsung	Galaxy Z Fold 5	6.10	5.11/2.64	0.24/0.52
Samsung	Galaxy Z Fold 7 (open/closed)	6.23	5.64/2.86	0.16/0.35
Samsung	Galaxy S25	5.78	2.78	0.28
Samsung	Galaxy S25 Ultra	6.41	3.06	0.32
Samsung	Galaxy S24	5.79	2.78	0.30
Samsung	Galaxy S24 Ultra	6.39	3.11	0.34
Samsung	Galaxy S24 Plus	6.24	2.99	0.30
Samsung	Galaxy S23	5.76	2.79	0.30
Samsung	Galaxy S23 Ultra	6.43	3.07	0.35
Samsung	Galaxy S23 Plus	6.21	3.00	0.30
Samsung	Galaxy S22	5.75	2.78	0.30
Samsung	Galaxy S22 Ultra	6.43	3.07	0.35
Samsung	Galaxy S22 Plus	6.19	2.98	0.30
Samsung	Galaxy S21	5.97	2.80	0.31
Samsung	Galaxy S21 Ultra	6.50	2.98	0.35
Samsung	Galaxy S21 Plus	6.35	2.98	0.31
Samsung	Galaxy S21 FE	6.13	2.93	0.31
Samsung	Galaxy S20	5.97	2.72	0.31
Samsung	Galaxy S20 Ultra	6.57	3.00	0.35
Samsung	Galaxy S20 Plus	6.37	2.90	0.31
Samsung	Galaxy S20 FE	6.30	2.93	0.33
Samsung	Galaxy S10	5.90	2.77	0.31
Samsung	Galaxy S10+	6.20	2.92	0.31
Samsung	Galaxy S10e	5.59	2.75	0.31
Samsung	Galaxy S9	5.81	2.7	0.33
Samsung	Galaxy S9+	6.22	2.9	0.33
Samsung	Galaxy S8	5.86	2.68	0.31
Samsung	Galaxy S8+	6.28	2.89	0.32
Samsung	Galaxy S7	5.60	2.74	0.31
Samsung	Galaxy S7 Edge	5.94	2.85	0.30
Samsung	Galaxy S6	5.64	2.77	0.27
Samsung	Galaxy S5	5.59	2.85	0.32
Samsung	Galaxy S5 Mini	5.16	2.55	0.36
Samsung	Galaxy Note 4	6.04	3.09	0.33
Samsung	Galaxy S III	5.37	2.77	0.34
Samsung	Galaxy S III Mini	4.78	2.48	0.39
Samsung	Galaxy S IV	5.37	2.74	0.31

eShield® Sterile Cover Tablet Sizing Chart

NOTE: If you have a protective case on your tablet, you may need the larger size cover to protect your device.

SKU: EC2300



Overall Dimensions: 14" x 18" (36 cm x 46 cm)
Functional Dimensions: 14" x 8.5" (36 cm x 22 cm)

Manufacturer	Tablet Model	Height (in.)	Width (in.)	Depth (in.)
Amazon	Fire 7	7.60	4.5	0.40
Amazon	Fire HD 8 (2022)	7.94	5.4	0.37
Amazon	Fire HD 8	8.40	5	0.40
Amazon	Fire HD 10 (2021, 9th Gen)	9.73	6.5	0.36
Amazon	Fire HD 10	10.30	6.3	0.40
Amazon	Fire Max 11	10.20	6.5	0.40
Apple	iPad (1st Gen)	9.56	7.47	0.50
Apple	iPad (2nd Gen)	9.50	7.31	0.35
Apple	iPad (3rd Gen)	9.50	7.31	0.37
Apple	iPad (4th Gen)	9.50	7.31	0.37
Apple	iPad (5th Gen)	9.40	6.67	0.30
Apple	iPad (6th Gen)	9.40	6.6	0.29
Apple	iPad (7th Gen)	9.87	6.85	0.30
Apple	iPad 4 (Retina Display)	9.46	7.31	0.37
Apple	iPad Air (1st Gen)	9.40	6.67	0.30
Apple	iPad Air (2nd Gen)	9.40	6.67	0.24
Apple	iPad Air (3rd Gen)	9.87	6.85	0.24
Apple	iPad Air 13-inch (M2)	11.04	8.46	0.24
Apple	iPad Mini (1st Gen)	7.87	5.3	0.28
Apple	iPad Mini (2nd Gen)	7.90	5.3	0.30
Apple	iPad Mini (3rd Gen)	7.90	5.3	0.30
Apple	iPad Mini (4th Gen)	8.00	5.3	0.24
Apple	iPad Mini (5th Gen)	8.00	5.3	0.24
Apple	iPad Mini (6th Gen)	7.69	5.3	0.25
Apple	iPad Pro (3rd Gen)	9.74	7.02	0.23
Apple	iPad Pro 9.7" (1st Gen)	9.40	6.67	0.24
Apple	iPad Pro 10.5" (2nd Gen)	9.80	6.8	0.24
Apple	iPad Pro 11" (A16, 2025)	9.79	7.07	0.24
Apple	iPad Pro 11" (M2)	9.74	7.02	0.24
Apple	iPad Pro 12.9" (3rd Gen)	11.04	8.46	0.23
Apple	iPad Pro 12.9" (4th Gen)	11.04	8.46	0.23
Apple	iPad Pro 12.9" (5th Gen)	11.05	8.46	0.25
Google	Pixel	10.20	6.7	0.30

Manufacturer	Tablet Model	Height (in.)	Width (in.)	Depth (in.)
HP	Chromebook x2	11.50	8.32	0.33
HP	Elite X2 1012 G1	11.80	8.4	0.30
Samsung	Galaxy Tab A	8.35	4.89	0.35
Samsung	Galaxy Tab S3	9.30	6.7	0.20
Samsung	Galaxy Tab S6	9.63	6.28	0.22
Samsung	Galaxy Tab S7	9.99	6.5	0.25
Samsung	Galaxy Tab S7+	11.22	7.28	0.22
Samsung	Galaxy Tab S7FE	11.21	7.28	0.25
Samsung	Galaxy Tab S9	10.01	6.53	0.23
Samsung	Galaxy Tab S10+	11.24	7.3	0.22
Samsung	Galaxy Tab S10 Ultra	12.85	8.21	0.21
Samsung	Galaxy Tab S10FE	10.01	6.53	0.24
Samsung	Galaxy Tab S10FE+	11.83	7.67	0.24