


C-Chip_{4ch}

*Count your cells,
for the most affordable price!*



C-Chip_{4ch}



Scan me



- ✓ New 4channel type disposable hemocytometer
- ✓ Created by the best-known microfluidic slide manufacturer
- ✓ High Quality, Low Price

C-Chip 4channel is a Disposable Plastic Hemocytometer

- No need for coverslips
- Disposable, eliminates the need for washing
- Reduces exposure to infectious samples and risk of contamination
- Guaranteed reproducibility and reliability
- Four tests
- Made of sturdy quartz grade optical plastic

Models

No.	Grid	Depth
DHC-N04	Neubauer Improved	100 μ m
DHC-B04	Burker Turk	100 μ m
DHC-F04	Fuchs Rosenthal	200 μ m




Applications

Blood analysis (Hematology) : Blood cell counting
Cell culture : Cell concentration measurement / Cell viability test
Microbiology : Fungal spore counting

Specifications

Dimension : 25 (W) x 75 (L) x Thickness 1.6 mm
Chamber depth : 100 μ m, 200 μ m
Chamber volume : 10 μ L (DHC-N04, B04)
20 μ L (DHC-F04)



website www.nanoentek.com
e-mail sales@nanoentek.com

NESCT-N04-001E (V.1.3)

NanoEntek, Inc.

Head Office

12F, 5, Digital-ro 26-gil, Guro-gu, Seoul, 08389, Korea
Tel +82-2-6220-7940 / Fax +82-2-6220-7999

NanoEntek America, Inc.

220 Bear Hill Road, Suite 102, Waltham, MA 02451, USA
Tel +1-781-472-2558 / Fax +1-781-790-5649

NanoEntek Europe | mts med-tech supplies GmbH

Lochhamerstr. 4a, 82152 Martinsried, Germany
Tel +49-89-21-55-38-43 / Fax +49-89-99-95-46-60