



Per- and Polyfluoroalkylated Substances (PFAS) Reference Materials

LGC Quality

ISO 17034 | ISO/IEC 17025 | ISO 9001

To learn more about our full range of
PFAS reference materials, get in touch
or visit: lgcstandards.com



PFAS Reference
Materials

PFAS

What you need to know

Per- and Polyfluoroalkylated Substances (PFAS) have garnered much attention in recent years as legislators around the world have begun to ban or limit their use.

The impact in Europe

Several PFAS compounds have come under the scrutiny of the Stockholm Convention on Persistent Organic Pollutants (POPs) including PFOA, PFOS, PFOSF, PFHxS, associated salts and other related compounds.

Towards the end of 2020, food regulators in the European Union set a limit for combined exposure to PFOA, PFOS and two other PFAS compounds in food – Perfluorononanoic Acid (PFNA) and Perfluorohexane Sulphonic Acid (PFHxS).

More regulation followed to regulate the presence of PFAS in drinking water and groundwater within the European Water Framework Directive to safeguard consumer and environmental health.

And most recently, authorities in Denmark, Germany, the Netherlands, Norway and Sweden are calling for further restriction or banning of PFAS compounds under REACH.

The impact in America

The U.S. Senate's PFAS Release Disclosure and Protection Act of 2019 requires the regulation of PFAS in relation to safe drinking water and management of toxic chemicals. The bill also required 172 specified PFAS to be immediately included in the Toxics Release Inventory (TRI) of the Environmental Protection Agency (EPA).

In 2020, the EPA outlined its strategic roadmap, including timelines and key actions, to safeguard public health and the environment, and hold polluters accountable.

Introducing our new ¹³C-labelled single solutions & mixtures with specific analytes for EPA Methods 533, 537.1, 1633 and the EU Drinking Water Directive.

Order now [!\[\]\(9c2e8d1b5bd77cb5c9f83b7a9cff79fd_img.jpg\)](#)





History & environmental impact

PFAS have been in use in manufacturing since the 1950s, for cookware coatings, food packaging, firefighting foam, textiles, clothing and other household goods. Unfortunately, the same properties that make PFAS so useable – hardness, oil and water resistance – also make these compounds incredibly environmentally persistent, lending them the nickname ‘forever chemicals’. This tenacity is largely down to the carbon and fluorine bonds that make up PFAS molecules, meaning that soil, air and water contamination may be widespread.

One report by the Centers for Disease Control (CDC) detected PFOS, PFOA, PFNA and PFHxS in the blood of more than 95% of the US population over 12 years old.

A growing global legal issue

The legal ramifications of the emerging PFAS problem are significant. To date, over 5000 PFAS-contaminated sites have been identified in the United States, and legal actions have been successfully brought against multiple manufacturers. Meanwhile, despite the Australian government's claim that there is “limited to no evidence” linking PFAS to human disease, one law firm still reached an in-principle agreement with the Australian Department of Defence. They secured a reported \$212.5 million AUD for three communities affected by PFAS contamination linked to PFAS compounds on army bases, fire stations and industrial sites from the 1970s to the 2000s.

Countless compounds, endless research

There are thousands of different types of PFAS compounds, all varying in the threat they pose to human, animal and environmental health. So the research and monitoring that's needed to set standards on a case-by-case basis is lengthy and highly resource-intensive. To date, a core group of more widely used compounds has dominated research and become regulated by bodies such as the EPA and European Union. For those that the industry is still gaining intelligence on, there are simply guidance measures in place.

References:

- 1 - <https://echa.europa.eu/hot-topics/perfluoroalkyl-chemicals-pfas>
- 2 - <https://chm.pops.int/TheConvention/ThePOPs/TheNewPOPs/tabid/2511/Default.aspx>
- 3 - <https://echa.europa.eu/-/next-steps-for-pfas-restriction-proposal>
- 4 - <https://www.epa.gov/newsreleases/epa-requires-toxics-releaseinventory-reporting-seven-additional-pfas>
- 5 - <https://www.epa.gov/pfas/press-releases-related-pfas>
- 6 - <https://www.epa.gov/pfas/pfas-strategic-roadmap-epas-commitments-action-2021-2024>
- 7 - https://www.ewg.org/interactive-maps/pfas_contamination/

Our PFAS portfolio

Over the following pages, we'll introduce you to a selection of our PFAS-related products, including both neat and single solutions. Our EPA method tailored solutions improve the efficiency of your laboratory, while our Stable Isotopically Labelled products (SIL) give you a range of options for your testing needs. You'll also see our latest PFAS mixtures which support the EU Drinking Water Directive, UCMR 5 and EPA Method 1633.

To see our complete range of reference materials, visit lgcstandards.com

Speak to our team 

PFAS reference materials

Product Listings

LGC brings you a comprehensive range of reference materials for PFAS testing through our Dr. Ehrenstorfer portfolio, and matrix materials through our preferred partners The National Institute of Standards and Technology (NIST) and Joint Research Centre (JRC).

Ehrenstorfer Quality

At Dr. Ehrenstorfer we place an emphasis on quality as part of our commitment to providing you with products you can trust. Here is what 'Ehrenstorfer Quality' means to us:

Producing to the highest standard

Dr. Ehrenstorfer reference materials are produced to the highest quality, with all analytical measurements performed under our ISO/IEC 17025 scope of accreditation and a leading portfolio of products produced according to our ISO 17034 accreditation. We use the most advanced analytical techniques to characterise our reference materials so that you can rely on the scientific integrity of the data contained in your Certificate of Analysis (CoA).

Ensuring confidence, from characterisation to implementation

We use real-time stability testing and expiry date management to give you confidence in your Dr. Ehrenstorfer reference materials and ensure you receive your products as certified, ready for your analysis. Our careful packaging choices protect your product during delivery and storage, and are made with your convenience and safety in mind.

Offering a unique and extensive portfolio

We produce an unrivalled portfolio of reference materials for food and environmental analysis, including unique substances, stable isotope labelled compounds and metabolites. Dr. Ehrenstorfer continues to be a global leader in pesticides and our range also features pharmaceutical and veterinary compounds, food related compounds, dyes, food packaging contaminants and more. We offer multiple formats including neat, single and multicomponent solutions.

Understanding your analytical needs

Through direct interactions with our customers and our expertise in the latest scientific and regulatory developments, we are able to quickly adapt our portfolio of reference materials to address your needs. We are committed to providing you with trusted solutions, today and tomorrow.

Providing expert support

At Dr. Ehrenstorfer we combine experience with continuous training to ensure that the latest knowledge and skills are being applied to producing your reference materials. As part of the LGC family, we are proud to connect with our customers across a global network, with dedicated local teams able to support your reference materials decisions and the implementation of our products in your analytical testing.

Our heritage, our vision, your guarantee

Dr. Ehrenstorfer is built upon 50 years of history in planning, developing, producing, analysing, packaging and delivering high quality reference materials to our customers around the world with speed and reliability. We are passionate about our work, which supports you in your science for a safer world.

Take a look at our evolution [↗](#)

PFAS reference materials

Product Listings



Since 1975, Dr. Ehrenstorfer™ has led the way in producing pesticide reference standards. Today, our portfolio has expanded to adapt to changing regulations and technology.

- Over 280 PFAS reference materials
- 50% of PFAS products produced to our ISO 17034 accreditation
- Popular tailor-made mixtures for EPA and other regulatory methods
- Centres of excellence in Germany and the United States
- As part of the LGC family, we're proud to connect with our customers across a global network, with dedicated local teams ready to support your needs



Order at
lgcstandards.com/drehrenstorfer

Connect with us now [🔗](#)

Introducing our latest PFAS mixtures to support the EU Drinking Water & Groundwater Directives, UK DWI, UCMR 5 and EPA Method 1633.

Improve efficiency

Specifically designed to meet testing requirements, our UCMR 5 mixture has 29 analytes, our EU Drinking Water Directive mixture has 20 analytes, EU Groundwater mixture 23 analytes, UK DWI mixture 48 analytes and our EPA Method 1633 mixture has 40 analytes.

- No need to conduct stability studies or develop your own mixtures
- Reduce your calibration runs and save preparation time

Free up resources

A higher concentration to give you clear signal intensity, even with the most difficult sample matrices.

- Increase your efficiency with higher sample throughput
- Save time compared to preparing your solutions in-house

Have confidence in your results

A robust ISO 17034 mixture designed by our team of expert chemists.

- Decrease your likelihood of downtime
- Reduce issues
- Avoid error-prone preparations

Order now DRE-A50000753MW, DRE-A50000754MW, DRE-A30000051MW, DRE-A30000079MW, DRE-K30000087 or contact us for more information. [📄](#)

New 13C-labelled single solutions and mixtures to support the EU Drinking Water Directive and EPA Method 1633

Improve efficiency

13C-labelled mixture of 20 analytes specifically designed to meet EU Drinking Water Directive testing requirements.

Product Number	Product Description	Pack size	Analytes
DRE-A30000095MW	EU Drinking Water Directive Labelled PFAS Mixture 2 µg/mL in Methanol:Water	1ml	20

13C-labelled mixtures of 24 and 7 analytes specifically for your EPA Method 1633 testing.

Product Number	Product Description	Pack size	Analytes
DRE-A50000792MW	EPA Method 1633 EIS Labelled PFAS Mixture 0.25-5 µg/mL in Methanol:Water	1ml	24
DRE-A50000794MW	EPA Method 1633 NIS Labelled PFAS Mixture 0.25-1 µg/mL in Methanol:Water	1ml	7

ES (Extraction Standard) – Add before extraction to track sample preparation recovery.

IS (Injection Standard) – Add before injection to correct for instrument response.

Have confidence in results

Robust ISO 17034 13C-labelled mixtures designed for maximum stability by our team of expert chemists.

Order now DRE-A30000095MW, DRE-A50000792MW and DRE-A50000794MW



To learn more, visit
lgcstandards.com/drehrenstorfer

Speak to our team [📄](#)

Contents

Mixtures	8
Stable Isotopes	31
Per- and Polyfluorocarboxylic Acids	33
Per- and Polyfluorosulfonic Acids	36
Fluorotelomer Alcohols	39
Perfluoroalkyl Sulfonamides	40
Perfluoroalkyl Acrylates	41
Polyfluorinated Phosphate Esters	41
Per- and Polyfluoroalkyl Ether Carboxylic Acids	42
Other PFAS	43
Matrix Materials	44
PFAS Proficiency Testing	45



Mixtures

Part Number	Part Description	Isotope label	No. of analytes
Water & Soil Analysis - EPA and US methods			
DRE-A50000753MW	EPA Method 1633 mixture 10 µg/mL in Methanol:Water	Unlabelled	40
DRE-A50000766MW	EPA Method 1633 PFAS Mixture 0.25-6.25 µg/mL in Methanol:Water	Unlabelled	40
DRE-A50000792MW	EPA Method 1633 EIS Labelled PFAS Mixture 0.25-5 µg/mL in Methanol:Water	13C, 2H	24
DRE-A50000794MW	EPA Method 1633 NIS Labelled PFAS Mixture 0.25-1 µg/mL in Methanol:Water	13C, 18O	7
DRE-A50000751MW	EPA Method 533 PFAS Mixture 10 µg/mL in Methanol:Water	Unlabelled	25
DRE-A50000765MW	EPA Method 533 Native Analyte Primary Dilution Standard 0.5 µg/mL in Methanol:Water	Unlabelled	25
DRE-A50000795MW	EPA Method 533 Isotope Performance Standard 1-3 µg/mL in Methanol:Water	13C	3
DRE-A50000790MW	EPA Method 533 Isotope Dilution Standard 0.5-2 µg/mL in Methanol:Water	13C	16
DRE-A50000767MW	EPA Method 537 Analyte Primary Dilution Standard 2 µg/mL in Methanol:Water	Unlabelled	18
DRE-A50000768MW	EPA Method 537 Analyte Primary Dilution Standard (branched/linear mix) 2 µg/mL in Methanol:Water	Unlabelled	18
DRE-A50000796MW	EPA Method 537 Internal Standard Primary Dilution Standard 1-4 µg/mL in Methanol:Water	13C, 2H	3
DRE-A50000152MW	EPA Method 537.1 PFAS Mixture 152 100 µg/mL in Methanol:Water	Unlabelled	18
DRE-A50000791MW	EPA Method 537.1 Surrogate Primary Dilution Standard 1-4 µg/mL in Methanol:Water	13C, 2H	4
DRE-A50000752MW	EPA Method 8327 PFAS Mixture 10 µg/mL in Methanol:Water	Unlabelled	24
DRE-A50000754MW	UCMR5 PFAS Mixture 10 µg/mL in Methanol:Water	Unlabelled	29
DRE-A50000738MW	Massachusetts PFAS Mixture 2 µg/mL in Methanol:Water	Unlabelled	6
Water Analysis - EU/UK			
DRE-A30000051MW	EU Drinking Water Directive PFAS Mixture 10 µg/mL in Methanol:Water	Unlabelled	20
DRE-A30000095MW	EU Drinking Water Directive Labelled PFAS Mixture 2 µg/mL in Methanol:Water	13C	20
DRE-K30000087	UK Drinking Water Inspectorate (DWI) PFAS-48 Kit 10 µg/mL in Methanol:Water	Unlabelled	48
DRE-A30000080MW	DIN EN 17892:2024-08 PFAS Mixture 10 µg/mL in Methanol:Water	Unlabelled	29
DRE-A30000079MW	EU Groundwater Testing PFAS Mixture (w/o C6O4) 10-100 µg/mL in Methanol:Water	Unlabelled	23
Food Analysis			
DRE-A30000094MW	AOAC 2023.003 PFAS Mixture 1-100 µg/mL in Methanol:Water	Unlabelled	30
DRE-A30000074MW	PFAS Mixture Food Testing EU 2022/1431 10 µg/mL in Methanol:Water	Unlabelled	28
DRE-A30000081MW	EU 2023/915 Food Testing PFAS-4 Mixture (PFOA, PFOS, PFNA, PFHxS) 1 µg/mL in Methanol:Water	Unlabelled	4
Cosmetics Analysis			
DRE-A30000083MW	Korea cosmetics testing PFAS Mixture 1 µg/mL in Methanol:Water	Unlabelled	5
DRE-A30000084MW	Korea cosmetics testing Labelled PFAS Mixture 1 µg/mL in Methanol:Water	13C	4
Other Mixtures			
DRE-A50000764MW	Native PFAS Precision and Recovery Standard Solution 1 µg/mL in Methanol:Water	Unlabelled	30
DRE-A50000647MW	PFASiMix 100 µg/mL in Methanol:Water (96:4)	Unlabelled	27
DRE-A30000070MW	PFAS EF-28 Native Mixture 1 µg/mL in Methanol:Water	Unlabelled	28
DRE-A50000760MW	PFAS Mixture (PFCAs, PFSA, FTs, FOSAs) 1-4 µg/mL in Methanol:Water	Unlabelled	25
DRE-A30000071MW	PFAS Superfund Mixture 2 µg/mL in Methanol:Water	Unlabelled	24
DRE-A50000151MW	PFAS Mixture 151 100 µg/mL in Methanol:Water	Unlabelled	24
DRE-GS09000218MW	PFAA Mixture 218 100 µg/mL in Methanol:Water	Unlabelled	11
DRE-A30000072MW	PFOA (C4-C14) Native Mixture 2 µg/mL in Methanol:Water	Unlabelled	11
DRE-A30000064MW	Ultrashort-chain-PFAS-7 (C1-C4) Mixture 50 µg/mL in Methanol:Water	Unlabelled	7
DRE-A30000082MW	Ultrashort-chain-PFAS-10 (C1-C4) Mixture 50 µg/mL in Methanol:Water	Unlabelled	10
DRE-A50000761MW	PFAS Mixture (DONA, 9CI-PF3ONS, 11CI-PF3OUdS, HFPO-DA) 2 µg/mL in Methanol:Water	Unlabelled	4
DRE-A50000762MW	PFAS Mixture (PF4OPeA, PF5OHxA, 3,6-OPFHpA, PFEEA) 2 µg/mL in Methanol:Water	Unlabelled	4

PFAS reference materials

Product Listings

Mixtures. Continued from previous page

Part Number	Part Description	Isotope label	No. of analytes
DRE-A50000763MW	PFAS Mixture (N-MeFOSA, N-EtFOSA, N-MeFOSE, N-EtFOSE) 1-10 µg/mL in Methanol:Water	Unlabelled	4
DRE-A50000759MW	PFAS Native X:3 FTCAs Mixture 4-20 µg/mL in Methanol:Water	Unlabelled	3
DRE-A30000088MW	Mass-Labelled PFC Injection Standards Solution 2 µg/mL in Methanol:Water	13C	4
DRE-A30000089MW	PFAS-C-ES Mixture 2 µg/mL in Methanol:Water	13C	13
DRE-A50000793MW	Mass-Labelled PFAS-24ES Extraction Mixture 1 µg/mL in Methanol:Water	13C, 2H	19

Product Number	Product Description			Pack size	Analytes
DRE-A50000753MW	EPA Method 1633 mixture 10 µg/mL in Methanol:Water			1ml	40 analytes
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluoro(2-ethoxyethane) sulfonic acid		PFEESA	113507-82-7	4	10 µg/mL
Perfluorobutanoic acid		PFBA	375-22-4	4	10 µg/mL
Perfluorobutanesulfonic acid		PFBS	375-73-5	4	10 µg/mL
Perfluoro-3-methoxypropanoic acid (PFMOPrA)		PFMPA	377-73-1	4	10 µg/mL
Perfluoro-3,6-dioxahexanoic acid		NFDHA	151772-58-6	5	10 µg/mL
Perfluoropentanoic acid		PFPeA	2706-90-3	5	10 µg/mL
Perfluoropentanesulfonic acid		PFPeS	2706-91-4	5	10 µg/mL
Perfluoro-4-methoxybutanoic acid (PFMOBA)		PFMBA	863090-89-5	5	10 µg/mL
Perfluoro-2-propoxypropanoic acid (PFPrOPrA)		PFPrOPrA	13252-13-6	6	10 µg/mL
Perfluorohexanoic acid		PFHxA	307-24-4	6	10 µg/mL
Perfluorohexanesulfonic acid		PFHxS	355-46-4	6	10 µg/mL
2H,2H,3H,3H-Perfluorohexanoic acid		3:3 FTCA	356-02-5	6	10 µg/mL
1H,1H,2H,2H-Perfluorohexanesulfonic acid		4:2 FTS	757124-72-4	6	10 µg/mL
Perfluoroheptanoic acid		PFHpA	375-85-9	7	10 µg/mL
Perfluoroheptanesulfonic acid		PFHpS	375-92-8	7	10 µg/mL
3H-Perfluoro-4,8-dioxanonanoic acid		DONA	919005-14-4	7	10 µg/mL
N-Ethyl-N-(2-hydroxyethyl)perfluorooctanesulfonamide		N-EtFOSE	1691-99-2	8	10 µg/mL
Perfluorooctane sulfonic acid		PFOS	1763-23-1	8	10 µg/mL
2-(N-Methylperfluorooctanesulfonamido)acetic acid		N-MeFOSAA	2355-31-9	8	10 µg/mL
N-(2-Hydroxyethyl)-N-methylperfluorooctanesulfonamide		N-MeFOSE	24448-09-7	8	10 µg/mL
1H,1H,2H,2H-Perfluorooctane sulfonic acid		6:2 FTS	27619-97-2	8	10 µg/mL
2-(N-Ethylperfluorooctanesulfonamido)acetic acid		N-EtFOSAA	2991-50-6	8	10 µg/mL
N-Methylperfluorooctanesulfonamide		MeFOSA	31506-32-8	8	10 µg/mL
Perfluorooctanoic acid		PFOA	335-67-1	8	10 µg/mL
N-Ethylperfluorooctanesulfonamide		EtFOSA	4151-50-2	8	10 µg/mL
<i>9-Chlorohexadecafluoro-3-oxanonane-1-sulfonic acid potassium</i>		9Cl-PF3ONS	73606-19-6	8	10 µg/mL
Perfluorooctane sulfonamide		PFOSA	754-91-6	8	10 µg/mL
2H,2H,3H,3H-Perfluorooctanoic acid		5:3 FTCA	914637-49-3	8	10 µg/mL
Perfluorononanoic acid		PFNA	375-95-1	9	10 µg/mL
<i>Perfluorononanesulfonic acid sodium</i>		PFNS	98789-57-2	9	10 µg/mL
Perfluorodecanoic acid		PFDA	335-76-2	10	10 µg/mL
Perfluorodecanesulfonic acid		PFDS	335-77-3	10	10 µg/mL
1H,1H,2H,2H-Perfluorodecanesulfonic acid		8:2 FTS	39108-34-4	10	10 µg/mL
1H,1H,2H,2H-Perfluorodecanoic acid		7:3 FTCA	812-70-4	10	10 µg/mL
<i>11-Chloroeicosafluoro-3-oxaundecane-1-sulfonic acid potassium</i>		11Cl-PF3OUdS	83329-89-9	10	10 µg/mL
Perfluoroundecanoic acid		PFUnDA	2058-94-8	11	10 µg/mL
<i>Perfluorododecanesulfonic acid sodium</i>		PFDoDS	1260224-54-1	12	10 µg/mL
Perfluorododecanoic acid		PFDoDA	307-55-1	12	10 µg/mL

PFAS reference materials

Product Listings

Mixtures. Continued from previous page

Analytes	Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluorotridecanoic acid	PFTDA	72629-94-8	13	10 µg/mL
Perfluorotetradecanoic acid	PFTeDA	376-06-7	14	10 µg/mL
Perfluorotetradecanoic acid	PFTeDA	376-06-7	14	10 µg/mL

Product Number	Product Description	Pack size	Analytes
DRE-A50000766MW	EPA Method 1633 PFAS Mixture 0.25-6.25 µg/mL in Methanol:Water	1ml	40 analytes

Analytes	Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluoro-4-oxapentanoic acid	PF4OPeA	377-73-1	4	0.5 µg/mL
Potassium perfluoro-1-butanefulfonate	PFBS	375-73-5	4	0.25 µg/mL
Potassium perfluoro(2-ethoxyethane)sulfonate	PFEESA	113507-82-7	4	0.5 µg/mL
Perfluoro-n-butanefulfonate	PFBA	375-22-4	5	1 µg/mL
Perfluoro-n-pentanoic acid	PFPeA	2706-90-3	5	0.5 µg/mL
Sodium perfluoro-1-pentanesulfonate	PFPeS	630402-22-1	5	0.25 µg/mL
Perfluoro-5-oxahexanoic acid	PF5OHxA	863090-89-5	6	0.5 µg/mL
Perfluoro-n-hexanoic acid	PFHxA	307-24-4	6	0.25 µg/mL
2,3,3,3-Tetrafluoro-2-(1,1,2,2,3,3,3-heptafluoropropoxy) propanoic acid	HFPO-DA	13252-13-6	6	1 µg/mL
3-Perfluoropropyl propanoic acid	FPrPA	356-02-5	6	1.25 µg/mL
Potassium perfluorohexanesulfonate (linear and branched isomers)	br-PFHxSK	355-46-4	6	0.25 µg/mL
Sodium 1H,1H,2H,2H-perfluorohexanesulfonate	4:2FTS	757124-72-4	6	1 µg/mL
Perfluoro-3,6-dioxahexanoic acid	3,6-OPFHpA	151772-58-6	7	0.5 µg/mL
Perfluoro-n-heptanoic acid	PFHpA	375-85-9	7	0.25 µg/mL
Sodium dodecafluoro-3H-4,8-dioxanonanoate	NaDONA	919005-14-4	7	1 µg/mL
Sodium perfluoro-1-heptanesulfonate	PFHpS	21934-50-9	7	0.25 µg/mL
Perfluorooctanoic acid (linear and branched isomers)	br-PFOA	335-67-1	8	0.25 µg/mL
Perfluorooctanesulfonamide (linear and branched isomers)	br-FOSA	754-91-6	8	0.25 µg/mL
N-Methylperfluorooctanesulfonamide (linear and branched isomers)	br-NMeFOSA	31506-32-8	8	0.25 µg/mL
N-Ethylperfluorooctanesulfonamide (linear and branched isomers)	br-NEtFOSA	4151-50-2	8	0.25 µg/mL
N-Methylperfluorooctanesulfonamidoacetic acid (linear and branched isomers)	br-NMeFOSAA	2355-31-9	8	0.25 µg/mL
N-Ethylperfluorooctanesulfonamidoacetic acid (linear and branched isomers)	br-NEtFOSAA	2991-50-6	8	0.25 µg/mL
2-(N-Methylperfluorooctanesulfonamido)ethanol (linear and branched isomers)	br-NMeFOSE	24448-09-7	8	2.5 µg/mL
2-(N-Ethylperfluorooctanesulfonamido)ethanol (linear and branched isomers)	br-NEtFOSE	1691-99-2	8	2.5 µg/mL
3-Perfluoropentyl propanoic acid	FPePA	914637-49-3	8	6.25 µg/mL
Potassium perfluorooctanesulfonate (linear and branched isomers)	br-PFOSK	1763-23-1	8	0.25 µg/mL
Sodium 1H,1H,2H,2H-perfluorooctanesulfonate	6:2FTS	27619-97-2	8	1 µg/mL
Potassium 9-chlorohexadecafluoro-3-oxanonane-1-sulfonate	9Cl-PF3ONS	73606-19-6	8	1 µg/mL
Perfluorononanoic acid (linear and branched isomers)	br-PFNA	375-95-1	9	0.25 µg/mL
Sodium perfluoro-1-nonanesulfonate	PFNS	98789-57-2	9	0.25 µg/mL
Perfluoro-n-decanoic acid	PFDA	335-76-2	10	0.25 µg/mL
3-Perfluoroheptyl propanoic acid	FHpPA	812-70-4	10	6.25 µg/mL
Sodium perfluoro-1-decanesulfonate	PFDS	2806-15-7	10	0.25 µg/mL
Sodium 1H,1H,2H,2H-perfluorodecanesulfonate	8:2FTS	39108-34-4	10	1 µg/mL
Potassium 11-chloroeicosafluoro-3-oxaundecane-1-sulfonate	11Cl-PF3OUdS	83329-89-9	10	1 µg/mL
Perfluoro-n-undecanoic acid	PFUnDA	2058-94-8	11	0.25 µg/mL
Perfluoro-n-dodecanoic acid	PFDoDA	307-55-1	12	0.25 µg/mL
Sodium perfluoro-1-dodecanesulfonate	PFDoS	1260224-54-1	12	0.25 µg/mL
Perfluoro-n-tridecanoic acid	PFTDA	72629-94-8	13	0.25 µg/mL
Perfluoro-n-tetradecanoic acid	PFTeDA	376-06-7	14	0.25 µg/mL

*Concentrations are always given for the free acids, which is also stated on our CoA

lgcstandards.com

Continued on next page

PFAS reference materials

Product Listings

Mixtures. Continued from previous page

Product Number	Product Description			Pack size	Analytes
DRE-A50000792MW	EPA Method 1633 EIS Labelled PFAS Mixture 0.25-5 µg/mL in Methanol:Water			1ml	24 analytes
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluoro-n-(13C4)butanoic acid		13C4-PFBA	1017281-29-6	4	2 µg/mL
Sodium perfluoro-1-(2,3,4-13C3)butanesulfonate		13C3-PFBS	2708218-84-0	4	0.5 µg/mL
Perfluoro-n-(13C5)pentanoic acid		13C5-PFPeA	2283397-79-3	5	1 µg/mL
Perfluoro-n-(1,2,3,4,6-13C5)hexanoic acid		13C5-PFHxA	2328024-54-8	6	0.5 µg/mL
Perfluoro-2-propoxypropanoic acid 13C3		13C3-PFPrOPrA	3030247-97-0	6	2 µg/mL
Sodium perfluoro-1-(1,2,3-13C3)hexanesulfonate		13C3-PFHxS	2708218-86-2	6	0.5 µg/mL
Sodium 1H,1H,2H,2H-perfluoro(1,2-13C2)hexanesulfonate		13C2-4:2 FTS	2708218-88-4	6	1 µg/mL
Perfluoro-n-(1,2,3,4-13C4)heptanoic acid		13C4-PFHpA	2328024-55-9	7	0.5 µg/mL
Perfluoro-n-(13C8)octanoic acid		13C8-PFOA	1350614-84-4	8	0.5 µg/mL
Perfluoro-1-(13C8)octanesulfonamide		13C8-FOSA	1365803-60-6	8	0.5 µg/mL
N-methyl-d3-perfluoro-1-octanesulfonamide		D3-N-MeFOSA	936109-37-4	8	0.5 µg/mL
N-ethyl-d5-perfluoro-1-octanesulfonamide		D5-N-EtFOSA	936109-40-9	8	0.5 µg/mL
Sodium perfluoro-1-(13C8)octanesulfonate		13C8-PFOS	2522762-16-7	8	0.5 µg/mL
Sodium 1H,1H,2H,2H-perfluoro(1,2-13C2)octanesulfonate		13C2-6:2 FTS	2708218-89-5	8	1 µg/mL
N-methyl-d3-perfluoro-1-octanesulfonamidoacetic acid		D3-N-MeFOSAA	1400690-70-1	8	1 µg/mL
N-ethyl-d5-perfluoro-1-octanesulfonamidoacetic acid		D5-N-EtFOSAA	1265205-97-7	8	1 µg/mL
2-(N-methyl-d3-perfluoro-1-octanesulfonamido)ethan-d4-ol		D7-N-MeFOSE	1265205-95-5	8	5 µg/mL
2-(N-ethyl-d5-perfluoro-1-octanesulfonamido)ethan-d4-ol		D9-N-EtFOSE	1265205-96-6	8	5 µg/mL
Perfluoro-n-(13C9)nonanoic acid		13C9-PFNA	2283397-80-6	9	0.25 µg/mL
Perfluoro-n-(1,2,3,4,5,6-13C6)decanoic acid		13C6-PFDA	2328024-56-0	10	0.25 µg/mL
Sodium 1H,1H,2H,2H-perfluoro(1,2-13C2)decanesulfonate		13C2-8:2 FTS	2708218-90-8	10	1 µg/mL
Perfluoro-n-(1,2,3,4,5,6,7-13C7)undecanoic acid		-	-	11	0.25 µg/mL
Perfluoro-n-(1,2-13C2)dodecanoic acid		13C2-PFDoA	960315-52-0	12	0.25 µg/mL
Perfluoro-n-(1,2-13C2)tetradecanoic acid		13C2-PFTeDA	2708218-82-8	14	0.25 µg/mL

Product Number	Product Description			Pack size	Analytes
DRE-A50000794MW	EPA Method 1633 NIS Labelled PFAS Mixture 0.25-1 µg/mL in Methanol:Water			1ml	7 analytes
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluoro-n-(2,3,4-13C3)butanoic acid		13C3-PFBA	2483735-33-5	4	1 µg/mL
Perfluoro-n-(1,2-13C2)hexanoic acid		13C2-PFHxA	960315-47-3	6	0.5 µg/mL
Sodium perfluoro-1-hexane(18O2)sulfonate		-	1585941-14-5	6	0.5 µg/mL
Sodium perfluoro-1-(1,2,3,4-13C4)octanesulfonate		13C4-PFOS	960315-53-1	8	0.5 µg/mL
Perfluoro-n-(1,2,3,4-13C4)octanoic acid		13C4-PFOA	960315-48-4	8	0.5 µg/mL
Perfluoro-n-(1,2,3,4,5-13C5)nonanoic acid		13C5-PFNA	960315-49-5	9	0.25 µg/mL
Perfluoro-n-(1,2-13C2)decanoic acid		13C2-PFDA	960315-50-8	10	0.25 µg/mL

PFAS reference materials

Product Listings

Mixtures. Continued from previous page

Product Number	Product Description		Pack size	Analytes	
DRE-A50000751MW	EPA Method 533 PFAS Mixture 10 µg/mL in Methanol:Water		1ml	25 analytes	
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluorobutanesulfonic Acid		PFBS	375-73-5	4	10 µg/mL
Perfluoro-3-methoxypropanoic Acid		PFMPA	377-73-1	4	10 µg/mL
Perfluoro(2-ethoxyethane) Sulfonic Acid		PFEESA	113507-82-7	4	10 µg/mL
Perfluorobutanoic Acid		PFBA	375-22-4	4	10 µg/mL
Perfluoropentanoic Acid		PFPeA	2706-90-3	5	10 µg/mL
Perfluoropentanesulfonic Acid		PFPeS	2706-91-4	5	10 µg/mL
Perfluoro-4-methoxybutanoic Acid		PFMBA	863090-89-5	5	10 µg/mL
Perfluoro-3,6-dioxaheptanoic Acid		NFDHA	151772-58-6	5	10 µg/mL
Perfluorohexanesulfonic Acid		PFHxS	355-46-4	6	10 µg/mL
1H,1H,2H,2H-Perfluorohexanesulfonic Acid		4:2 FTS	757124-72-4	6	10 µg/mL
Perfluoro-2-propoxypropanoic acid (PFPrOPrA)		PFPrOPrA	13252-13-6	6	10 µg/mL
Perfluorohexanoic Acid		PFHxA	307-24-4	6	10 µg/mL
Perfluoroheptanesulfonic Acid		PFHpS	375-92-8	7	10 µg/mL
Perfluoroheptanoic Acid		PFHpA	375-85-9	7	10 µg/mL
3H-Perfluoro-4,8-dioxanonanoic Acid		DONA	919005-14-4	7	10 µg/mL
9Cl-Perfluoro-3-oxanonanesulfonic Acid		PFHxA	756426-58-1	8	10 µg/mL
1H,1H,2H,2H-Perfluorooctanesulfonic Acid		6:2 FTS	27619-97-2	8	10 µg/mL
Perfluorooctanesulfonic Acid		PFOS	1763-23-1	8	10 µg/mL
Perfluorooctanoic Acid		PFOA	335-67-1	8	10 µg/mL
Perfluorononanoic Acid		PFNA	375-95-1	9	10 µg/mL
11Cl-Perfluoro-3-oxaundecanesulfonic Acid		11Cl-PF3OUdS	763051-92-9	10	10 µg/mL
1H,1H,2H,2H-Perfluorodecane sulfonic Acid		8:2 FTS	39108-34-4	10	10 µg/mL
Perfluorodecanoic Acid		PFDA	335-76-2	10	10 µg/mL
Perfluoroundecanoic Acid		PFUnDA	2058-94-8	11	10 µg/mL
Perfluorododecanoic Acid		PFDoDA	307-55-1	12	10 µg/mL

Product Number	Product Description		Pack size	Analytes	
DRE-A50000765MW	EPA Method 533 Native Analyte Primary Dilution Standard 0.5 µg/mL in Methanol:Water		1ml	25 analytes	
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluoro-n-butanoic acid		PFBA	375-22-4	4	0.5 µg/mL
Perfluoro-4-oxapentanoic acid		PF4OPeA	377-73-1	4	0.5 µg/mL
Potassium perfluoro-1-butanefulfonate		L-PFBS	29420-49-3	4	0.5 µg/mL
Potassium perfluoro(2-ethoxyethane)sulfonate		PFEESA	117205-07-9	4	0.5 µg/mL
Perfluoro-n-pentanoic acid		PFPeA	2706-90-3	5	0.5 µg/mL
Perfluoro-5-oxahexanoic acid		PF5OHxA	863090-89-5	5	0.5 µg/mL
Perfluoro-3,6-dioxahexanoic acid		3,6-OPFHpA	151772-58-6	5	0.5 µg/mL
Sodium perfluoro-1-pentanesulfonate		L-PFPeS	-	5	0.5 µg/mL
Perfluoro-n-hexanoic acid		PFHxA	307-24-4	6	0.5 µg/mL
Perfluoro-2-propoxypropanoic acid (PFPrOPrA)		PFPPrOPrA	13252-13-6	6	0.5 µg/mL
Potassium perfluorohexanesulfonate		PFHxSK	-	6	0.5 µg/mL
Sodium 1H,1H,2H,2H-perfluorohexanesulfonate		4:2FTS	27619-93-8	6	0.5 µg/mL
Perfluoro-n-heptanoic acid		PFHpA	375-85-9	7	0.5 µg/mL
Sodium perfluoro-1-heptanesulfonate		L-PFHpS	21934-50-9	7	0.5 µg/mL
Perfluorooctanoic acid		PFOA	335-67-1	8	0.5 µg/mL
Potassium perfluorooctanesulfonate		PFOSK	-	8	0.5 µg/mL
Sodium 1H,1H,2H,2H-perfluorooctanesulfonate		6:2FTS	27619-94-9	8	0.5 µg/mL

*Concentrations are always given for the free acids, which is also stated on our CoA

lgcstandards.com

Continued on next page

PFAS reference materials

Product Listings

Mixtures. Continued from previous page

Analytes	Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Potassium 9-chlorohexadecafluoro-3-oxanonane-1-sulfonate	9Cl-PF3ONS	73606-19-6	8	0.5 µg/mL
Perfluorononanoic acid	PFNA	375-95-1	9	0.5 µg/mL
Sodium dodecafluoro-3H-4,8-dioxanonanoate	NaDONA	2250081-67-3	9	0.5 µg/mL
Perfluoro-n-decanoic acid	PFDA	335-76-2	10	0.5 µg/mL
Sodium 1H,1H,2H,2H-perfluorodecanesulfonate	8:2FTS	27619-96-1	10	0.5 µg/mL
Potassium 11-chloroeicosafuoro-3-oxaundecane-1-sulfonate	11Cl-PF3OUdS	83329-89-9	10	0.5 µg/mL
Perfluoro-n-undecanoic acid	PFUnDA	2058-94-8	11	0.5 µg/mL
Perfluoro-n-dodecanoic acid	PFDoDA	307-55-1	12	0.5 µg/mL

Product Number	Product Description			Pack size	Analytes
DRE-A50000795MW	EPA Method 533 Isotope Performance Standard 1-3 µg/mL in Methanol:Water			1ml	3 analytes
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluoro-n-(2,3,4-13C3) butanoic acid		13C3-PFBA	2483735-33-5	4	1 µg/mL
Perfluoro-n-(1,2-13C2) octanoic acid		13C2-PFOA	864071-08-9	8	1 µg/mL
Sodium perfluoro-1-(1,2,3,4-13C4) octanesulfonate		13C4-PFOS	960315-53-1	8	3 µg/mL

Product Number	Product Description		Pack size	Analytes	
DRE-A50000790MW	EPA Method 533 Isotope Dilution Standard 0.5-2 µg/mL in Methanol:Water		1ml	16 analytes	
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluoro-n-(13C4) butanoic acid		13C4-PFBA	1017281-29-6	4	0.5 µg/mL
Sodium perfluoro-1-(2,3,4-13C3) butanesulfonate		13C3-PFBS	2708218-84-0	4	0.5 µg/mL
Perfluoro-n-(13C5) pentanoic acid		13C5-PFPeA	2283397-79-3	5	0.5 µg/mL
Perfluoro-n-(1,2,3,4,6-13C5) hexanoic acid		13C5-PFHxA	2328024-54-8	6	0.5 µg/mL
Perfluoro-2-propoxypropanoic acid 13C3		13C3-PFPrOPrA	3030247-97-0	6	0.5 µg/mL
Sodium perfluoro-1-(1,2,3-13C3) hexanesulfonate		13C3-PFHxS	2708218-86-2	6	0.5 µg/mL
Sodium 1H,1H,2H,2H-perfluoro(1,2-13C2) hexanesulfonate		13C2-4:2 FTS	2708218-88-4	6	2 µg/mL
Perfluoro-n-(1,2,3,4-13C4) heptanoic acid		13C4-PFHpA	2328024-55-9	7	0.5 µg/mL
Perfluoro-n-(13C8) octanoic acid		13C8-PFOA	1350614-84-4	8	0.5 µg/mL
Sodium perfluoro-1-(13C8) octanesulfonate		13C8-PFOSNa	2522762-16-7	8	0.5 µg/mL
Sodium 1H,1H,2H,2H-perfluoro(1,2-13C2) octanesulfonate		13C2-6:2 FTS	2708218-89-5	8	2 µg/mL
Perfluoro-n-(13C9) nonanoic acid		13C9-PFNA	2283397-80-6	9	0.5 µg/mL
Perfluoro-n-(1,2,3,4,5,6-13C6) decanoic acid		13C6-PFDA	2328024-56-0	10	0.5 µg/mL
Sodium 1H,1H,2H,2H-perfluoro(1,2-13C2) decanesulfonate		13C2-8:2 FTS	2708218-90-8	10	2 µg/mL
Perfluoro-n-(1,2,3,4,5,6,7-13C7) undecanoic acid		13C7-PFUnDA	-	11	0.5 µg/mL
Perfluoro-n-(1,2-13C2) dodecanoic acid		13C3-PFDoDA	960315-52-0	12	0.5 µg/mL

Product Number	Product Description		Pack size	Analytes	
DRE-A50000767MW	EPA Method 537 Analyte Primary Dilution Standard 2 µg/mL in Methanol:Water		1ml	18 analytes	
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluorobutanesulfonic acid		PFBS	375-73-5	4	2 µg/mL
Perfluoro-2-propoxypropanoic acid (PFPrOPrA)		PFPrOPrA	13252-13-6	6	2 µg/mL
Perfluorohexanesulfonic acid		PFHxS	355-46-4	6	2 µg/mL
Perfluorohexanoic acid		PFHxA	307-24-4	6	2 µg/mL
3H-Perfluoro-4,8-dioxanonanoic acid		DONA	919005-14-4	7	2 µg/mL

*Concentrations are always given for the free acids, which is also stated on our CoA

lgcstandards.com

Continued on next page

PFAS reference materials

Product Listings

Mixtures. Continued from previous page

Analytes	Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluoroheptanoic acid	PFHpA	375-85-9	7	2 µg/mL
2-(N-Ethylperfluorooctanesulfonamido)acetic acid	N-EtFOSAA	2991-50-6	8	2 µg/mL
2-(N-Methylperfluorooctanesulfonamido)acetic acid	N-MeFOSAA	2355-31-9	8	2 µg/mL
<i>9-Chloro-perfluoro-3-oxanonanesulfonic acid potassium</i>	9Cl-PF3ONS	73606-19-6	8	2 µg/mL
Perfluorooctane sulfonic acid	PFOS	1763-23-1	8	2 µg/mL
Perfluorooctanoic acid	PFOA	335-67-1	8	2 µg/mL
Perfluorononanoic acid	PFNA	375-95-1	9	2 µg/mL
<i>11-Chloroeicosafluoro-3-oxaundecane-1-sulfonic acid potassium</i>	11Cl-PF3OUdS	83329-89-9	10	2 µg/mL
Perfluorodecanoic acid	PFDA	335-76-2	10	2 µg/mL
Perfluoroundecanoic acid	PFUnDA	2058-94-8	11	2 µg/mL
Perfluorododecanoic acid	PFDoDA	307-55-1	12	2 µg/mL
Perfluorotridecanoic acid	PFTrDA	72629-94-8	13	2 µg/mL
Perfluorotetradecanoic acid	PFTeDA	376-06-7	14	2 µg/mL

Product Number	Product Description	Pack size	Analytes
DRE-A50000768MW	EPA Method 537 Analyte Primary Dilution Standard (branched/75linear mix) 2 µg/mL in Methanol:Water	1ml	18 analytes

Analytes	Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluorobutanesulfonic acid	PFBS	375-73-5	4	2 µg/mL
Perfluoro-2-propoxypropanoic acid (PFPrOPrA)	PFPrOPrA	13252-13-6	6	2 µg/mL
Perfluorohexanesulfonic acid	PFHxS	355-46-4	6	2 µg/mL
Perfluorohexanoic acid	PFHxA	307-24-4	6	2 µg/mL
3H-Perfluoro-4,8-dioxanonanoic acid	DONA	919005-14-4	7	2 µg/mL
Perfluoroheptanoic acid	PFHpA	375-85-9	7	2 µg/mL
2-(N-Ethylperfluorooctanesulfonamido)acetic acid	N-EtFOSAA	2991-50-6	8	2 µg/mL
2-(N-Methylperfluorooctanesulfonamido)acetic acid	N-MeFOSAA	2355-31-9	8	2 µg/mL
<i>9-Chloro-perfluoro-3-oxanonanesulfonic acid potassium</i>	9Cl-PF3ONS	73606-19-6	8	2 µg/mL
Perfluorooctane sulfonic acid	PFOS	1763-23-1	8	2 µg/mL
Perfluorooctanoic acid	PFOA	335-67-1	8	2 µg/mL
Perfluorononanoic acid	PFNA	375-95-1	9	2 µg/mL
<i>11-Chloroeicosafluoro-3-oxaundecane-1-sulfonic acid potassium</i>	11Cl-PF3OUdS	83329-89-9	10	2 µg/mL
Perfluorodecanoic acid	PFDA	335-76-2	10	2 µg/mL
Perfluoroundecanoic acid	PFUnDA	2058-94-8	11	2 µg/mL
Perfluorododecanoic acid	PFDoDA	307-55-1	12	2 µg/mL
Perfluorotridecanoic acid	PFTrDA	72629-94-8	13	2 µg/mL
Perfluorotetradecanoic acid	PFTeDA	376-06-7	14	2 µg/mL

Product Number	Product Description	Pack size	Analytes
DRE-A50000796MW	EPA Method 537 Internal Standard Primary Dilution Standard 1-4 µg/mL in Methanol:Water	1ml	3 analytes

Analytes	Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluoro-n-(1,2-13C2)octanoic acid	13C2-PFOA	864071-08-9	8	1 µg/mL
Sodium perfluoro-1-(1,2,3,4-13C4)octanesulfonate	13C4-PFOS	960315-53-1	8	3 µg/mL
N-methyl-d3 -perfluoro-1-octanesulfonamidoacetic acid	D3-N-MeFOSAA	1400690-70-1	10	4 µg/mL

*Concentrations are always given for the free acids, which is also stated on our CoA

PFAS reference materials

Product Listings

Mixtures. Continued from previous page

Product Number	Product Description		Pack size	Analytes	
DRE-A50000152MW	EPA Method 537.1 PFAS Mixture 152 100 µg/mL in Methanol:Water		1ml	18 analytes	
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluorobutanesulfonic acid		PFBS	375-73-5	4	100 µg/mL
Perfluoro-2-propoxypropanoic acid (PFPrOPrA)		PFPrOPrA	13252-13-6	6	100 µg/mL
Perfluorohexanesulfonic acid		PFHxS	355-46-4	6	100 µg/mL
Perfluorohexanoic acid		PFHxA	307-24-4	6	100 µg/mL
3H-Perfluoro-4,8-dioxanonanoic acid		DONA	919005-14-4	7	100 µg/mL
Perfluoroheptanoic acid		PFHpA	375-85-9	7	100 µg/mL
2-(N-Ethylperfluorooctanesulfonamido)acetic acid		N-EtFOSAA	2991-50-6	8	100 µg/mL
2-(N-Methylperfluorooctanesulfonamido)acetic acid		N-MeFOSAA	2355-31-9	8	100 µg/mL
9-Chloro-perfluoro-3-oxanonanesulfonic acid potassium		9Cl-PF3ONS	73606-19-6	8	100 µg/mL
Perfluorooctane sulfonic acid		PFOS	1763-23-1	8	100 µg/mL
Perfluorooctanoic acid		PFOA	335-67-1	8	100 µg/mL
Perfluorononanoic acid		PFNA	375-95-1	9	100 µg/mL
11-Chloroeicosafuoro-3-oxaundecane-1-sulfonic acid potassium		11Cl-PF3OUdS	83329-89-9	10	100 µg/mL
Perfluorodecanoic acid		PFDA	335-76-2	10	100 µg/mL
Perfluoroundecanoic acid		PFUnDA	2058-94-8	11	100 µg/mL
Perfluorododecanoic acid		PFDoDA	307-55-1	12	100 µg/mL
Perfluorotridecanoic acid		PFTrDA	72629-94-8	13	100 µg/mL
Perfluorotetradecanoic acid		PFTeDA	376-06-7	14	100 µg/mL

Product Number	Product Description			Pack size	Analytes
DRE-A50000791MW	EPA Method 537.1 Surrogate Primary Dilution Standard 1-4 µg/mL in Methanol:Water			1ml	4 analytes
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluoro-n-(1,2-13C2)hexanoic acid		13C2-PFHxA	960315-47-3	6	1 µg/mL
Perfluoro-2-propoxypropanoic acid 13C3		13C3-PFPrOPrA	3030247-97-0	6	1 µg/mL
N-ethyl-d5 -perfluoro-1-octanesulfonamidoacetic acid		D5-N-EtFOSAA	1265205-97-7	8	4 µg/mL
Perfluoro-n-(1,2-13C2)decanoic acid		13C2-PFDA	960315-50-8	10	1 µg/mL

Product Number	Product Description		Pack size	Analytes	
DRE-A50000752MW	EPA Method 8327 PFAS Mixture 10 µg/mL in Methanol:Water		1ml	24 analytes	
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluorobutanoic Acid		PFBA	375-22-4	4	10 µg/mL
Perfluorobutanesulfonic Acid		PFBS	375-73-5	4	10 µg/mL
Perfluoropentanesulfonic Acid		PFPeS	2706-91-4	5	10 µg/mL
Perfluoropentanoic Acid		PFPeA	2706-90-3	5	10 µg/mL
Perfluorohexanoic Acid		PFHxA	307-24-4	6	10 µg/mL
Perfluorohexanesulfonic Acid		PFHxS	355-46-4	6	10 µg/mL
1H,1H,2H,2H-Perfluorohexanesulfonic Acid		4:2 FTS	757124-72-4	6	10 µg/mL
Perfluoroheptanoic Acid		PFHpA	375-85-9	7	10 µg/mL
Perfluoroheptanesulfonic Acid		PFHpS	375-92-8	7	10 µg/mL
2-(N-Ethylperfluorooctanesulfonamido)acetic Acid		N-EtFOSAA	2991-50-6	8	10 µg/mL
Perfluorooctanesulfonic Acid		PFOS	1763-23-1	8	10 µg/mL

PFAS reference materials

Product Listings

Mixtures. Continued from previous page

Analytes	Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluorooctanoic Acid	PFOA	335-67-1	8	10 µg/mL
2-(N-Methylperfluorooctanesulfonamido)acetic Acid	N-MeFOSAA	2355-31-9	8	10 µg/mL
Perfluorooctane Sulfonamide	PFOSA	754-91-6	8	10 µg/mL
1H,1H,2H,2H-Perfluorooctanesulfonic Acid	6:2 FTS	27619-97-2	8	10 µg/mL
Perfluorononanoic Acid	PFNA	375-95-1	9	10 µg/mL
Perfluorononanesulfonic Acid	PFNS	68259-12-1	9	10 µg/mL
1H,1H,2H,2H-Perfluorodecanesulfonic Acid	8:2 FTS	39108-34-4	10	10 µg/mL
Perfluorodecanoic Acid	PFDA	335-76-2	10	10 µg/mL
Perfluorodecanesulfonic Acid	PFDS	335-77-3	10	10 µg/mL
Perfluoroundecanoic Acid	PFUnDA	2058-94-8	11	10 µg/mL
Perfluorododecanoic Acid	PFDODA	307-55-1	12	10 µg/mL
Perfluorotridecanoic Acid	PFTDA	72629-94-8	13	10 µg/mL
Perfluorotetradecanoic Acid	PFTeDA	376-06-7	14	10 µg/mL

Product Number	Product Description		Pack size	Analytes	
DRE-A50000754MW	UCMR5 PFAS Mixture 10 µg/mL in Methanol:Water		1ml	29 analytes	
Analytes	Abbreviation	CAS No.	Carbon Chain Length	*Concentration	
Perfluorobutanoic Acid	PFBA	375-22-4	4	10 µg/mL	
Perfluorobutanesulfonic Acid	PFBS	375-73-5	4	10 µg/mL	
Perfluoro-3-methoxypropanoic Acid	PFMPA	377-73-1	4	10 µg/mL	
Perfluoro(2-ethoxyethane) Sulfonic Acid	PFEESA	113507-82-7	4	10 µg/mL	
Perfluoropentanoic Acid	PFPeA	2706-90-3	5	10 µg/mL	
Perfluoropentanesulfonic Acid	PFPeS	2706-91-4	5	10 µg/mL	
Perfluoro-4-methoxybutanoic Acid	PFMBA	863090-89-5	5	10 µg/mL	
Perfluoro-3,6-dioxahexanoic Acid	NFDHA	151772-58-6	5	10 µg/mL	
Perfluoro-2-propoxypropanoic acid (PFPrOPrA)	PFPrOPrA	13252-13-6	6	10 µg/mL	
Perfluorohexanoic Acid	PFHxA	307-24-4	6	10 µg/mL	
1H,1H,2H,2H-Perfluorohexanesulfonic Acid	4:2 FTS	757124-72-4	6	10 µg/mL	
Perfluorohexanesulfonic Acid	PFHxS	355-46-4	6	10 µg/mL	
Perfluoroheptanoic Acid	PFHpA	375-85-9	7	10 µg/mL	
3H-Perfluoro-4,8-dioxanonanoic Acid	DONA	919005-14-4	7	10 µg/mL	
Perfluoroheptanesulfonic Acid	PFHpS	375-92-8	7	10 µg/mL	
Perfluorooctanoic Acid	PFOA	335-67-1	8	10 µg/mL	
1H,1H,2H,2H-Perfluorooctanesulfonic Acid	6:2 FTS	27619-97-2	8	10 µg/mL	
2-(N-Methylperfluorooctanesulfonamido)acetic Acid	N-MeFOSAA	2355-31-9	8	10 µg/mL	
2-(N-Ethylperfluorooctanesulfonamido)acetic Acid	N-EtFOSAA	2991-50-6	8	10 µg/mL	
Perfluorooctanesulfonic Acid	PFOS	1763-23-1	8	10 µg/mL	
9CI-Perfluoro-3-oxanonanesulfonic Acid	9CI-PF3ONS	756426-58-1	8	10 µg/mL	
Perfluorononanoic Acid	PFNA	375-95-1	9	10 µg/mL	
11CI-Perfluoro-3-oxaundecanesulfonic Acid	11CI-PF3OUdS	763051-92-9	10	10 µg/mL	
Perfluorodecanoic Acid	PFDA	335-76-2	10	10 µg/mL	
1H,1H,2H,2H-Perfluorodecanesulfonic Acid	8:2 FTS	39108-34-4	10	10 µg/mL	
Perfluoroundecanoic Acid	PFUnDA	2058-94-8	11	10 µg/mL	
Perfluorododecanoic Acid	PFDODA	307-55-1	12	10 µg/mL	
Perfluorotridecanoic Acid	PFTDA	72629-94-8	13	10 µg/mL	
Perfluorotetradecanoic Acid	PFTeDA	376-06-7	14	10 µg/mL	

*Concentrations are always given for the free acids, which is also stated on our CoA

lgcstandards.com

Continued on next page

PFAS reference materials

Product Listings

Mixtures. Continued from previous page

Product Number	Product Description		Pack size	Analytes	
DRE-A50000738MW	Massachusetts PFAS Mixture 2 µg/mL in Methanol:Water		1ml	6 analytes	
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluorohexanesulfonic Acid		PFHxS	355-46-4	6	2 µg/mL
Perfluoroheptanoic Acid		PFHpA	375-85-9	7	2 µg/mL
Perfluorooctanesulfonic Acid		PFOS	1763-23-1	8	2 µg/mL
Perfluorooctanoic Acid		PFOA	335-67-1	8	2 µg/mL
Perfluorononanoic Acid		PFNA	375-95-1	9	2 µg/mL
Perfluorodecanoic Acid		PFDA	335-76-2	10	2 µg/mL

Product Number	Product Description	Pack size		Analytes	
DRE-A30000051MW	EU Drinking Water Directive PFAS Mixture 10 µg/mL in Methanol:Water	1ml		20 analytes	
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluorobutanoic Acid		PFBA	375-22-4	4	10 µg/mL
Perfluorobutanesulfonic Acid		PFBS	375-73-5	4	10 µg/mL
Perfluoropentanoic Acid		PFPeA	2706-90-3	5	10 µg/mL
Perfluoropentanesulfonic Acid		PFPeS	2706-91-4	5	10 µg/mL
Perfluorohexanesulfonic Acid		PFHxS	355-46-4	6	10 µg/mL
Perfluorohexanoic Acid		PFHxA	307-24-4	6	10 µg/mL
Perfluoroheptanoic Acid		PFHpA	375-85-9	7	10 µg/mL
Perfluoroheptanesulfonic Acid		PFHpS	375-92-8	7	10 µg/mL
Perfluorooctanesulfonic Acid		PFOS	1763-23-1	8	10 µg/mL
Perfluorooctanoic Acid		PFOA	335-67-1	8	10 µg/mL
Perfluorononanoic Acid		PFNA	375-95-1	9	10 µg/mL
Perfluorononanesulfonic Acid		PFNS	68259-12-1	9	10 µg/mL
Perfluorodecanoic Acid		PFDA	335-76-2	10	10 µg/mL
Perfluorodecanesulfonic Acid		PFDS	335-77-3	10	10 µg/mL
Perfluoroundecanoic Acid		PFUnDA	2058-94-8	11	10 µg/mL
Perfluoroundecanesulfonic Acid		PFUnDS	749786-16-1	11	10 µg/mL
Perfluorododecanoic Acid		PFDoDA	307-55-1	12	10 µg/mL
Perfluorododecanesulfonic Acid		PFDoDS	79780-39-5	12	10 µg/mL
Perfluorotridecanoic Acid		PFTTrDA	72629-94-8	13	10 µg/mL
Perfluorotridecanesulfonic Acid		PFTTrDS	791563-89-8	13	10 µg/mL

Product Number	Product Description			Pack size	Analytes
DRE-A30000095MW	EU Drinking Water Directive Labelled PFAS Mixture 2 µg/mL in Methanol:Water			1ml	20 analytes
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluorobutanesulfonic acid 13C3 (2,3,4-13C3)		13C3-PFBS	2708218-83-9	4	2 µg/mL
Perfluorobutanoic acid 13C3 (2,3,4-13C3)		13C3-PFBA	2483735-33-5	4	2 µg/mL
Perfluoropentanesulfonic acid 13C3 (1,2,3-13C3)		13C3-PFPeS	-	5	2 µg/mL
Perfluoropentanoic acid 13C5		13C5-PFPeA	2283397-79-3	5	2 µg/mL
Perfluorohexanesulfonic acid 13C3 (1,2,3-13C3)		13C3-PFHxS	2708218-85-1	6	2 µg/mL
Perfluorohexanoic acid 13C5 (1,2,3,4,6-13C5)		13C5-PFHxA	2328024-54-8	6	2 µg/mL

*Concentrations are always given for the free acids, which is also stated on our CoA

lgcstandards.com

Continued on next page

PFAS reference materials

Product Listings

Mixtures. Continued from previous page

Analytes	Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluoroheptanesulfonic acid 13C3 (1,2,3-13C3)	13C3-PFHpS	-	7	2 µg/mL
Perfluoroheptanoic acid 13C4 (1,2,3,4-13C4)	13C4-PFHpA	2328024-55-9	7	2 µg/mL
Perfluorooctanesulfonic acid 13C4 (1,2,3,4-13C4)	13C4-PFOS	960354-19-2	8	2 µg/mL
Perfluorooctanoic acid 13C8	13C8-PFOA	1350614-84-4	8	2 µg/mL
Perfluorononanesulfonic acid 13C3 (1,2,3-13C3)	13C3-PFNS	-	9	2 µg/mL
Perfluorononanoic acid 13C9	13C9-PFNA	2283397-80-6	9	2 µg/mL
Perfluorodecanesulfonic acid 13C3 (1,2,3-13C3)	13C3-PFDS	-	10	2 µg/mL
Perfluorodecanoic acid 13C6 (1,2,3,4,5,6-13C6)	13C6-PFDA	2328024-56-0	10	2 µg/mL
Perfluoroundecanesulfonic acid 13C3 (1,2,3-13C3)	13C3-PFUnDS	-	11	2 µg/mL
Perfluoroundecanoic acid 13C7 (1,2,3,4,5,6,7-13C7)	13C7-PFUnDA	-	11	2 µg/mL
Perfluorododecanesulfonic acid 13C3 (1,2,3-13C3)	13C3-PFDoDS	-	12	2 µg/mL
Perfluorododecanoic acid 13C2 (1,2-13C2)	13C2-PFDoDA	960315-52-0	12	2 µg/mL
Perfluorotridecanesulfonic acid 13C3 (1,2,3-13C3)	13C3-PFTrDS	-	13	2 µg/mL
Perfluorotridecanoic acid 13C2 (1,2-13C2)	13C2-PFTrDA	-	13	2 µg/mL

Product Number	Product Description			Pack size	Analytes
DRE-K30000087	UK Drinking Water Inspectorate (DWI) PFAS-48 Kit 10 µg/mL in Methanol:Water			3 x 1ml	48
DRE-A30000085MW	UK Drinking Water Inspectorate (DWI) PFAS-48 Mixture A 10 µg/mL in Methanol:Water			1ml	25
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluorobutanoic acid		PFBA	375-22-4	4	10 µg/mL
Perfluorobutanesulfonic acid		PFBS	375-73-5	4	10 µg/mL
Perfluoropentanoic acid		PFPeA	2706-90-3	5	10 µg/mL
Perfluoropentanesulfonic acid		PFPeS	2706-91-4	5	10 µg/mL
Perfluorohexanoic acid		PFHxA	307-24-4	6	10 µg/mL
Perfluorohexanesulfonic acid		PFHxS	355-46-4	6	10 µg/mL
1H,1H,2H,2H-Perfluorohexanesulfonic acid		4:2 FTS	757124-72-4	6	10 µg/mL
Perfluoroheptanoic acid		PFHpA	375-85-9	7	10 µg/mL
Perfluoroheptanesulfonic acid		PFHpS	375-92-8	7	10 µg/mL
Perfluorooctanoic acid		PFOA	335-67-1	8	10 µg/mL
Perfluorooctane sulfonic acid		PFOS	1763-23-1	8	10 µg/mL
1H,1H,2H,2H-Perfluorooctane sulfonic acid		6:2 FTS	27619-97-2	8	10 µg/mL
Perfluorononanoic acid		PFNA	375-95-1	9	10 µg/mL
Perfluorononanesulfonic acid		PFNS	98789-57-2	9	10 µg/mL
Perfluorodecanoic acid		PFDA	335-76-2	10	10 µg/mL
Perfluorodecanesulfonic acid		PFDS	335-77-3	10	10 µg/mL
1H,1H,2H,2H-Perfluorodecanesulfonic acid		8:2 FTS	39108-34-4	10	10 µg/mL
Perfluoroundecanoic acid		PFUnDA	2058-94-8	11	10 µg/mL
Perfluoroundecanesulfonic acid		PFUnDS	3016286-99-7	11	10 µg/mL
Perfluorododecanoic acid		PFDoDA	307-55-1	12	10 µg/mL
Perfluorododecanesulfonic acid		PFDoDS	1260224-54-1	12	10 µg/mL
Perfluorotridecanoic acid		PFTTrDA	72629-94-8	13	10 µg/mL
Perfluorotetradecanoic acid		PFTeDA	376-06-7	14	10 µg/mL
Perfluorohexadecanoic acid		PFHxDA	67905-19-5	16	10 µg/mL
Perfluorooctadecanoic acid		PFODA	16517-11-6	18	10 µg/mL

*Concentrations are always given for the free acids, which is also stated on our CoA

lgcstandards.com

Continued on next page

PFAS reference materials

Product Listings

Mixtures. Continued from previous page

Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluorotetradecanoic acid		PFTeDA	376-06-7	14	10 µg/mL
Perfluorohexadecanoic acid		PFHxDA	67905-19-5	16	10 µg/mL
Perfluorooctadecanoic acid		PFODA	16517-11-6	18	10 µg/mL
Product Number	Product Description			Pack size	Analytes
DRE-A30000086MW	UK Drinking Water Inspectorate (DWI) PFAS-48 Mixture B 10 µg/mL in Methanol:Water			1ml	20
Perfluoro-3-methoxypropanoic acid (PFMOPrA)		PFMPA	377-73-1	4	10 µg/mL
Perfluoro(2-ethoxyethane) sulfonic acid		PFEESA	113507-82-7	4	10 µg/mL
Perfluorobutane sulfonamide		PFBSA	30334-69-1	4	10 µg/mL
Perfluoro-3,6-dioxahexanoic acid		NFDHA	151772-58-6	5	10 µg/mL
Perfluoro-4-methoxybutanoic acid (PFMOBA)		PFMOBA	863090-89-5	5	10 µg/mL
Perfluoro-2-propoxypropanoic acid (PFPrOPrA)		PFPrOPrA	13252-13-6	6	10 µg/mL
Perfluorohexane sulfonamide		PFHxSA	41997-13-1	6	10 µg/mL
3H-Perfluoro-4,8-dioxanonanoic acid		DONA	919005-14-4	7	10 µg/mL
Perfluoro-4-ethylcyclohexylsulfonic acid potassium		PFECHS	646-83-3	8	10 µg/mL
9-Chlorohexadecafluoro-3-oxanonane-1-sulfonic acid		9CI-PF3ONS	756426-58-1	8	10 µg/mL
Perfluorooctane sulfonamide		PFOSA	754-91-6	8	10 µg/mL
N-Methylperfluorooctanesulfonamide		MeFOSA	31506-32-8	8	10 µg/mL
Sulfluramid		EtFOSA	4151-50-2	8	10 µg/mL
N-(2-Hydroxyethyl)-N-methylperfluorooctanesulfonamide		N-MeFOSE	24448-09-7	8	10 µg/mL
N-Ethyl-N-(2-hydroxyethyl)perfluorooctanesulfonamide		N-EtFOSE	1691-99-2	8	10 µg/mL
2-(N-Methylperfluorooctanesulfonamido)acetic acid		N-MeFOSAA	2355-31-9	8	10 µg/mL
2-(N-Ethylperfluorooctanesulfonamido)acetic acid		N-EtFOSAA	2991-50-6	8	10 µg/mL
6:2 fluorotelomer sulfonamide alkylbetaine		-	34455-29-3	8	10 µg/mL
Perfluoro-2,5-dimethyl-3,6-dioxanonanoic acid (HFPO-TA)		HFPO-TA	13252-14-7	9	10 µg/mL
11-Chloroeicosafluoro-3-oxaundecane-1-sulfonic acid		11CI-PF3OUdS	763051-92-9	10	10 µg/mL
Product Number	Product Description			Pack size	Analytes
DRE-A30000098MW	UK Drinking Water Inspectorate (DWI) PFAS-48 Mixture C 10 µg/mL in Methanol:Water			1ml	3
2H,2H,3H,3H-Perfluorohexanoic acid		3:3 FTCA	356-02-5	6	10 µg/mL
2H,2H,3H,3H-Perfluorooctanoic acid		5:3 FTCA	914637-49-3	8	10 µg/mL
1H,1H,2H,2H-Perfluorodecanoic acid		7:3 FTCA	812-70-4	10	10 µg/mL

Product Number	Product Description			Pack size	Analytes
DRE-A30000080MW	DIN EN 17892:2024-08 PFAS Mixture 10 µg/mL in Methanol:Water			1ml	29 analytes
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluorobutanoic acid		PFBA	375-22-4	4	10 µg/mL
Perfluorobutanesulfonic acid		PFBS	375-73-5	4	10 µg/mL
Perfluoro-3-methoxypropanoic acid		PFMPA	377-73-1	4	10 µg/mL
Perfluoropentanoic acid		PFPeA	2706-90-3	5	10 µg/mL
Perfluoropentanesulfonic acid		PFPeS	2706-91-4	5	10 µg/mL
Perfluorohexanoic acid		PFHxA	307-24-4	6	10 µg/mL
Perfluoro-2-propoxypropanoic acid (PFPrOPrA)		PFPrOPrA	13252-13-6	6	10 µg/mL
Perfluorohexanesulfonic acid		PFHxS	355-46-4	6	10 µg/mL
1H,1H,2H,2H-Perfluorohexanesulfonic acid		4:2 FTS	757124-72-4	6	10 µg/mL

*Concentrations are always given for the free acids, which is also stated on our CoA

lgcstandards.com

Continued on next page

PFAS reference materials

Product Listings

Mixtures. Continued from previous page

Analytes	Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluoroheptanoic acid	PFHpA	375-85-9	7	10 µg/mL
Perfluoroheptanesulfonic acid	PFHpS	375-92-8	7	10 µg/mL
3H-Perfluoro-4,8-dioxanonanoic acid	DONA	919005-14-4	7	10 µg/mL
Perfluorooctanesulfonic acid	PFOS	1763-23-1	8	10 µg/mL
1H,1H,2H,2H-Perfluorooctanesulfonic acid	6:2 FTS	27619-97-2	8	10 µg/mL
2-(N-Ethylperfluorooctanesulfonamido)acetic acid	N-EtFOSAA	2991-50-6	8	10 µg/mL
Perfluorooctanoic acid	PFOA	335-67-1	8	10 µg/mL
Perfluorooctanesulfonamide	PFOSA	754-91-6	8	10 µg/mL
9CI-Perfluoro-3-oxanonanesulfonic acid	9CI-PF3ONS	756426-58-1	8	10 µg/mL
Perfluorononanoic acid	PFNA	375-95-1	9	10 µg/mL
Perfluorononanesulfonic acid	PFNS	68259-12-1	9	10 µg/mL
Perfluorodecanoic acid	PFDA	335-76-2	10	10 µg/mL
Perfluorodecanesulfonic acid	PFDS	335-77-3	10	10 µg/mL
1H,1H,2H,2H-Perfluorodecanesulfonic acid	8:2 FTS	39108-34-4	10	10 µg/mL
Perfluoroundecanoic acid	PFUnDA	2058-94-8	11	10 µg/mL
Perfluoroundecanesulfonic acid	PFUnDS	749786-16-1	11	10 µg/mL
Perfluorododecanoic acid	PFDoDA	307-55-1	12	10 µg/mL
Perfluorododecanesulfonic acid	PFDoDS	79780-39-5	12	10 µg/mL
Perfluorotridecanoic acid	PFTTrDA	72629-94-8	13	10 µg/mL
Perfluorotridecanesulfonic acid	PFTTrDS	791563-89-8	13	10 µg/mL

Product Number	Product Description		Pack size	Analytes	
DRE-A30000079MW	EU Groundwater Testing PFAS Mixture (w/o C6O4) 10-100 µg/mL in Methanol:Water		1ml	23	
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluorobutane sulfonic acid		PFBS	375-73-5	4	10 µg/mL
Perfluorobutanoic acid		PFBA	375-22-4	4	10 µg/mL
Perfluoropentanoic acid		PFPeA	2706-90-3	5	10 µg/mL
Perfluoropentane sulfonic acid		PFPeS	2706-91-4	5	10 µg/mL
Perfluorohexane sulfonic acid		PFHxS	355-46-4	6	10 µg/mL
Perfluorohexanoic acid		PFHxA	307-24-4	6	10 µg/mL
Perfluoro-2-propoxypropanoic acid (PFPrOPrA)		PFPrOPrA	13252-13-6	6	10 µg/mL
Perfluoroheptanoic acid		PFHpA	375-85-9	7	10 µg/mL
Perfluoroheptane sulfonic acid		PFHpS	375-92-8	7	10 µg/mL
3H-Perfluoro-4,8-dioxanonanoic acid		DONA	919005-14-4	7	10 µg/mL
Perfluorooctanoic acid		PFOA	335-67-1	8	10 µg/mL
Perfluorooctane sulfonic acid		PFOS	1763-23-1	8	10 µg/mL
2-(Perfluorohexyl)ethyl alcohol (6:2 FTOH)		6:2 FTOH	647-42-7	8	100 µg/mL
Perfluorononanoic acid		PFNA	375-95-1	9	10 µg/mL
Perfluorodecanoic acid		PFDA	335-76-2	10	10 µg/mL
Perfluorodecane sulfonic acid		PFDS	335-77-3	10	10 µg/mL
2-(Perfluorooctyl)ethanol (8:2 FTOH)		8:2 FTOH	678-39-7	10	100 µg/mL
Perfluoroundecanoic acid		PFUnDA	2058-94-8	11	10 µg/mL
Perfluorododecanoic acid		PFDoDA	307-55-1	12	10 µg/mL
Perfluorotridecanoic acid		PFTTrDA	72629-94-8	13	10 µg/mL

*Concentrations are always given for the free acids, which is also stated on our CoA

lgcstandards.com

Continued on next page

PFAS reference materials

Product Listings

Mixtures. Continued from previous page

Analytes	Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluorotetradecanoic acid	PFTeDA	376-06-7	14	10 µg/mL
Perfluorohexadecanoic acid	PFHxDA	67905-19-5	16	10 µg/mL
Perfluorooctadecanoic acid	PFODA	16517-11-6	18	10 µg/mL

Product Number	Product Description		Pack size	Analytes	
DRE-A30000094MW	AOAC 2023.003 PFAS Mixture 1-100 µg/mL in Methanol:Water		1ml	30	
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluorobutanoic acid		PFBA	375-22-4	4	100 µg/mL
Perfluorobutanesulfonic acid		PFBS	375-73-5	4	10 µg/mL
Perfluoropentanoic acid		PFPeA	2706-90-3	5	100 µg/mL
Perfluoropentanesulfonic acid		PFPeS	2706-91-4	5	10 µg/mL
Perfluorohexanoic acid		PFHxA	307-24-4	6	10 µg/mL
Perfluoro-2-propoxypropanoic acid (PFPrOPrA)		PFPrOPrA	13252-13-6	6	10 µg/mL
Perfluorohexanesulfonic acid		PFHxS	355-46-4	6	1 µg/mL
1H,1H,2H,2H-Perfluorohexanesulfonic acid		4:2 FTS	757124-72-4	6	10 µg/mL
Perfluoroheptanoic acid		PFHpA	375-85-9	7	10 µg/mL
3H-Perfluoro-4,8-dioxanonoic acid		DONA	919005-14-4	7	10 µg/mL
Perfluoroheptanesulfonic acid		PFHpS	375-92-8	7	10 µg/mL
Perfluorooctanoic acid		PFOA	335-67-1	8	1 µg/mL
Perfluorooctane sulfonamide		PFOSA	754-91-6	8	10 µg/mL
Perfluorooctanesulfonic acid		PFOS	1763-23-1	8	1 µg/mL
1H,1H,2H,2H-Perfluorooctanesulfonic acid		6:2 FTS	27619-97-2	8	10 µg/mL
9-Chlorohexadecafluoro-3-oxanonane-1-sulfonic acid potassium		9Cl-PF3ONS	756426-58-1	8	10 µg/mL
Perfluorononanoic acid		PFNA	375-95-1	9	1 µg/mL
Perfluorononanesulfonic acid		PFNS	791563-89-8	9	10 µg/mL
Perfluorodecanoic acid		PFDA	335-76-2	10	10 µg/mL
Perfluorodecanesulfonic acid		PFDS	335-77-3	10	10 µg/mL
1H,1H,2H,2H-Perfluorodecanesulfonic acid		8:2 FTS	39108-34-4	10	10 µg/mL
11-Chloroeicosafluoro-3-oxaundecane-1-sulfonic acid		11Cl-PF3OUdS	763051-92-9	10	10 µg/mL
Perfluoroundecanoic acid		PFUnDA	2058-94-8	11	10 µg/mL
Perfluoroundecanesulfonic acid		PFUnDS	749786-16-1	11	10 µg/mL
1H,1H,2H,2H-Perfluorododecane sulfonic acid		10:2 FTS	120226-60-0	12	10 µg/mL
Perfluorododecanoic acid		PFDoDA	307-55-1	12	10 µg/mL
Perfluorododecanesulfonic acid		PFDoDS	335-77-3	12	10 µg/mL
Perfluorotridecanoic acid		PFTTrDA	72629-94-8	13	10 µg/mL
Perfluorotridecanesulfonic acid		PFTTrDS	791563-89-8	13	10 µg/mL
Perfluorotetradecanoic acid		PFTeDA	376-06-7	14	10 µg/mL

Product Number	Product Description		Pack size	Analytes	
DRE-A30000074MW	PFAS Mixture Food Testing EU 2022/1431 10 µg/mL in Methanol:Water		1ml	28 analytes	
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluorobutanesulfonic acid		PFBS	73606-19-6	4	10 µg/mL
Perfluorobutanoic acid		PFBA	375-22-4	4	10 µg/mL
Perfluoropentanesulfonic acid		PFPeS	2706-91-4	5	10 µg/mL
Perfluoropentanoic acid		PFPeA	2706-90-3	5	10 µg/mL
Perfluorohexanesulfonic acid		PFHxS	355-46-4	6	10 µg/mL
Perfluorohexanoic acid		PFHxA	307-24-4	6	10 µg/mL
Perfluoro-2-propoxypropanoic acid (PFPrOPrA)		HFPO-DA	13252-13-6	6	10 µg/mL
3H-Perfluoro-4,8-dioxanonoic acid		DONA	919005-14-4	7	10 µg/mL

*Concentrations are always given for the free acids, which is also stated on our CoA

lgcstandards.com

Continued on next page

PFAS reference materials

Product Listings

Mixtures. Continued from previous page

Analytes	Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluoroheptanesulfonic acid	PFHpS	375-92-8	7	10 µg/mL
Perfluoroheptanoic acid	PFHpA	375-85-9	7	10 µg/mL
Capstone Product A	-	80475-32-7	8	10 µg/mL
Capstone Product B	-	34455-29-3	8	10 µg/mL
<i>9Cl-Perfluoro-3-oxanonanesulfonic acid sodium</i>	9Cl-PF3ONS	756426-58-1	8	10 µg/mL
Perfluorooctanesulfonamide	PFOSA	754-91-6	8	10 µg/mL
Perfluorooctanesulfonic acid	PFOS	1763-23-1	8	10 µg/mL
Perfluorooctanoic acid	PFOA	335-67-1	8	10 µg/mL
Perfluorononanoic acid	PFNA	375-95-1	9	10 µg/mL
Perfluorononanesulfonic acid	PFNS	98789-57-2	9	10 µg/mL
<i>11Cl-Perfluoro-3-oxaundecanesulfonic acid sodium</i>	11Cl-PF3OUdS	83329-89-9	10	10 µg/mL
Perfluorodecanesulfonic acid	PFDS	335-77-3	10	10 µg/mL
Perfluorodecanoic acid	PFDA	335-76-2	10	10 µg/mL
<i>Perfluoroundecanesulfonic acid sodium</i>	PFUnDS	3016286-99-7	11	10 µg/mL
Perfluoroundecanoic acid	PFUnDA	2058-94-8	11	10 µg/mL
Perfluorododecanoic acid	PFDoDA	307-55-1	12	10 µg/mL
<i>Perfluorododecanesulfonic acid sodium</i>	PFDoDS	1260224-54-1	12	10 µg/mL
<i>Perfluorotridecanesulfonic acid sodium</i>	PFTTrDS	174675-49-1	13	10 µg/mL
Perfluorotridecanoic acid	PFTTrDA	72629-94-8	13	10 µg/mL
Perfluorotetradecanoic acid	PFTeDA	376-06-7	14	10 µg/mL

Product Number	Product Description			Pack size	Analytes
DRE-A30000081MW	EU 2023/915 Food Testing PFAS-4 Mixture (PFOA, PFOS, PFNA, PFHxS) 1 µg/mL in Methanol:Water			1ml	4 analytes
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluorohexanesulfonic acid		PFHxS	355-46-4	6	1 µg/mL
Perfluorooctanesulfonic acid		PFOS	1763-23-1	8	1 µg/mL
Perfluorooctanoic acid		PFOA	335-67-1	8	1 µg/mL
Perfluorononanoic acid		PFNA	375-95-1	9	1 µg/mL

Product Number	Product Description	Pack size	Analytes		
DRE-A30000083MW	Korea cosmetics testing PFAS Mixture 1 µg/mL in Methanol:Water	1ml	5 analytes		
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Tetraconazole		-	112281-77-3	2	1 µg/mL
Perfluorooctanoic acid		PFOA	335-67-1	8	1 µg/mL
Perfluorooctanesulfonic acid		PFOS	1763-23-1	8	1 µg/mL
Perfluorononanoic acid		PFNA	375-95-1	9	1 µg/mL
Perfluorodecanoic acid		PFDA	335-76-2	10	1 µg/mL

PFAS reference materials

Product Listings

Mixtures. Continued from previous page

Product Number	Product Description			Pack size	Analytes
DRE-A30000084MW	Korea cosmetics testing Labelled PFAS Mixture 1 µg/mL in Methanol:Water			1ml	4 analytes
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluorooctanoic acid 13C4 (1,2,3,4-13C4)		13C4-PFOA	960315-48-4	8	1 µg/mL
Perfluorooctanesulfonic acid 13C4 (1,2,3,4-13C4) sodium		13C4-PFOS	960315-53-1	8	1 µg/mL
Perfluorononanoic acid 13C5 (1,2,3,4,5-13C5)		13C5-PFNA	960315-49-5	9	1 µg/mL
Perfluorododecanoic acid 13C2 (1,2-13C2)		13C2-PFDoDA	960315-52-0	12	1 µg/mL

Product Number	Product Description			Pack size	Analytes
DRE-A50000764MW	Native PFAS Precision and Recovery Standard Solution 1 µg/mL in Methanol:Water			1ml	30
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluoro(2-ethoxyethane) sulfonic acid		PFEESA	113507-82-7	4	100 µg/mL
Perfluoro-3-methoxypropanoic acid		PFMOPrA	377-73-1	4	100 µg/mL
Perfluorobutanesulfonic acid		PFBS	375-73-5	4	100 µg/mL
Perfluorobutanoic acid		PFBA	375-22-4	4	100 µg/mL
Perfluoro-3,6-dioxahexanoic acid		PFDDHA	151772-58-6	5	100 µg/mL
Perfluoro-4-methoxybutanoic acid		PFMOBA	863090-89-5	5	100 µg/mL
Perfluoropentanesulfonic Acid		PFPeS	2706-91-4	5	100 µg/mL
Perfluoropentanoic acid		PFPeA	2706-90-3	5	100 µg/mL
1H,1H,2H,2H-Perfluorohexanesulfonic acid sodium		4:2 FTS	27619-93-8	6	100 µg/mL
Perfluoro-2-propoxypropanoic acid		PFPPrOPrA	13252-13-6	6	100 µg/mL
Perfluorohexanesulfonic acid		PFHxS	355-46-4	6	100 µg/mL
Perfluorohexanoic acid		PFHxA	307-24-4	6	100 µg/mL
3H-Perfluoro-4,8-dioxanonanoic acid		DONA	919005-14-4	7	100 µg/mL
Perfluoroheptanesulfonic acid		PFHpS	375-92-8	7	100 µg/mL
Perfluoroheptanoic acid		PFHpA	375-85-9	7	100 µg/mL
1H,1H,2H,2H-Perfluorooctane sulfonic acid sodium		6:2 FTS	27619-94-9	8	100 µg/mL
2-(N-Ethylperfluorooctanesulfonamido)acetic acid		N-EtFOSAA	2991-50-6	8	100 µg/mL
2-(N-Methylperfluorooctanesulfonamido)acetic acid		N-MeFOSAA	2355-31-9	8	100 µg/mL
Perfluorooctane sulfonic acid		PFOS	1763-23-1	8	100 µg/mL
Perfluorooctanoic acid		PFOA	335-67-1	8	100 µg/mL
Perfluorononanoic acid		PFNA	375-95-1	9	100 µg/mL
1,1,2,2H-Perfluorodecanesulfonic acid sodium salt		8:2 FTS	27619-96-1	10	100 µg/mL
Perfluorodecanoic acid		PFDA	335-76-2	10	100 µg/mL
Perfluoroundecanoic acid		PFUnDA	2058-94-8	11	100 µg/mL
Perfluorododecanoic acid		PFDoDA	307-55-1	12	100 µg/mL
Perfluorotridecanoic acid		PFTTrDA	72629-94-8	13	100 µg/mL
Perfluorotetradecanoic acid		PFTeDA	376-06-7	14	100 µg/mL

*Concentrations are always given for the free acids, which is also stated on our CoA

PFAS reference materials

Product Listings

Mixtures. Continued from previous page

Product Number	Product Description		Pack size	Analytes	
DRE-A50000647MW	PFASiMix 100 µg/mL in Methanol: Water		1ml	27 analytes	
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluoro(2-ethoxyethane) sulfonic acid		PFEESA	113507-82-7	4	100 µg/mL
Perfluoro-3-methoxypropanoic acid		PFMOPrA	377-73-1	4	100 µg/mL
Perfluorobutanesulfonic acid		PFBS	375-73-5	4	100 µg/mL
Perfluorobutanoic acid		PFBA	375-22-4	4	100 µg/mL
Perfluoro-3,6-dioxahexanoic acid		PFDA	151772-58-6	5	100 µg/mL
Perfluoro-4-methoxybutanoic acid		PFMOBA	863090-89-5	5	100 µg/mL
Perfluoropentanesulfonic Acid		PFPeS	2706-91-4	5	100 µg/mL
Perfluoropentanoic acid		PFPeA	2706-90-3	5	100 µg/mL
1H,1H,2H,2H-Perfluorohexanesulfonic acid sodium		4:2 FTS	27619-93-8	6	100 µg/mL
Perfluoro-2-propoxypropanoic acid		PFPPrOPrA	13252-13-6	6	100 µg/mL
Perfluorohexanesulfonic acid		PFHxS	355-46-4	6	100 µg/mL
Perfluorohexanoic acid		PFHxA	307-24-4	6	100 µg/mL
3H-Perfluoro-4,8-dioxanonanoic acid		DONA	919005-14-4	7	100 µg/mL
Perfluoroheptanesulfonic acid		PFHpS	375-92-8	7	100 µg/mL
Perfluoroheptanoic acid		PFHpA	375-85-9	7	100 µg/mL
1H,1H,2H,2H-Perfluorooctane sulfonic acid sodium		6:2 FTS	27619-94-9	8	100 µg/mL
2-(N-Ethylperfluorooctanesulfonamido)acetic acid		N-EtFOSAA	2991-50-6	8	100 µg/mL
2-(N-Methylperfluorooctanesulfonamido)acetic acid		N-MeFOSAA	2355-31-9	8	100 µg/mL
Perfluorooctane sulfonic acid		PFOS	1763-23-1	8	100 µg/mL
Perfluorooctanoic acid		PFOA	335-67-1	8	100 µg/mL
Perfluorononanoic acid		PFNA	375-95-1	9	100 µg/mL
1,1,2,2H-Perfluorodecanesulfonic acid sodium salt		8:2 FTS	27619-96-1	10	100 µg/mL
Perfluorodecanoic acid		PFDA	335-76-2	10	100 µg/mL
Perfluoroundecanoic acid		PFUnDA	2058-94-8	11	100 µg/mL
Perfluorododecanoic acid		PFDoDA	307-55-1	12	100 µg/mL
Perfluorotridecanoic acid		PFTTrDA	72629-94-8	13	100 µg/mL
Perfluorotetradecanoic acid		PFTeDA	376-06-7	14	100 µg/mL

Product Number	Product Description		Pack size	Analytes	
DRE-A30000070MW	PFAS EF-28 Native Mixture 1 µg/mL in Methanol:Water		1ml	28 analytes	
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluorobutanesulfonic acid		PFBS	375-73-5	4	1 µg/mL
Perfluorobutanoic acid		PFBA	375-22-4	4	1 µg/mL
Perfluoropentanesulfonic acid		PFPeS	2706-91-4	5	1 µg/mL
Perfluoropentanoic acid		PFPeA	2706-90-3	5	1 µg/mL
1H,1H,2H,2H-Perfluorohexanesulfonic acid		4:2 FTS	757124-72-4	6	1 µg/mL
Perfluorohexanesulfonic acid		PFHxS	355-46-4	6	1 µg/mL
Perfluorohexanoic acid		PFHxA	307-24-4	6	1 µg/mL
Perfluoroheptanesulfonic acid		PFHpS	375-92-8	7	1 µg/mL
Perfluoroheptanoic acid		PFHpA	375-85-9	7	1 µg/mL
1H,1H,2H,2H-Perfluorooctane sulfonic acid		6:2 FTS	27619-97-2	8	1 µg/mL

*Concentrations are always given for the free acids, which is also stated on our CoA

lgcstandards.com

Continued on next page

PFAS reference materials

Product Listings

Mixtures. Continued from previous page

Analytes	Abbreviation	CAS No.	Carbon Chain Length	*Concentration
2-(N-Ethylperfluorooctanesulfonamido)acetic acid	N-EtFOSAA	2991-50-6	8	1 µg/mL
2-(N-Methylperfluorooctanesulfonamido)acetic acid	N-MeFOSAA	2355-31-9	8	1 µg/mL
N-Methylperfluorooctanesulfonamide	MeFOSA	31506-32-8	8	1 µg/mL
Perfluorooctane sulfonamide	PFOSA	754-91-6	8	1 µg/mL
Perfluorooctanoic acid	PFOA	335-67-1	8	1 µg/mL
Perfluorooctane sulfonic acid	PFOS	1763-23-1	8	1 µg/mL
Perfluorononanoic acid	PFNA	375-95-1	9	1 µg/mL
Perfluorodecanesulfonic acid	PFDS	335-77-3	10	1 µg/mL
1H,1H,2H,2H-Perfluorodecanesulfonic acid	8:2 FTS	39108-34-4	10	1 µg/mL
Perfluorodecanoic acid	PFDA	335-76-2	10	1 µg/mL
Perfluoroundecanoic acid	PFUnDA	2058-94-8	11	1 µg/mL
1H,1H,2H,2H-Perfluorododecane sulfonic acid	10:2 FTS	120226-60-0	12	1 µg/mL
Perfluorododecanoic acid	PFDoDA	307-55-1	12	1 µg/mL
Perfluorotridecanoic acid	PFTriDA	72629-94-8	13	1 µg/mL
Perfluorotetradecanoic acid	PFTeDA	376-06-7	14	1 µg/mL
Perfluorohexadecanoic acid	PFHxDA	67905-19-5	16	1 µg/mL
Perfluorooctadecanoic acid	PFODA	16517-11-6	18	1 µg/mL
Bis[2-(perfluorooctyl)ethyl] phosphate	8:2 Di-PAPS	678-41-1	10/10	1 µg/mL

Product Number	Product Description	Pack size	Analytes		
DRE-A50000760MW	PFAS Mixture (PFCAs, PFSA, FTSs, FOSAAs) 1-4 µg/mL in Methanol:Water	1ml	25 analytes		
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluorobutanoic acid		PFBA	375-22-4	4	4 µg/mL
Perfluorobutanesulfonic acid		PFBS	375-73-5	4	1 µg/mL
Perfluoropentanoic acid		PFPeA	2706-90-3	5	2 µg/mL
Perfluoropentanesulfonic acid		PFPeS	2706-91-4	5	1 µg/mL
Perfluorohexanoic acid		PFHxA	307-24-4	6	1 µg/mL
Perfluorohexanesulfonic acid		PFHxS	355-46-4	6	1 µg/mL
1H,1H,2H,2H-Perfluorohexane sulfonic acid		4:2FTS	757124-72-4	6	4 µg/mL
Perfluoroheptanoic acid		PFHpA	375-85-9	7	1 µg/mL
Perfluoroheptanesulfonic acid		PFHpS	375-92-8	7	1 µg/mL
Perfluorooctanoic acid		PFOA	335-67-1	8	1 µg/mL
Perfluorooctane sulfonamide		FOSA	754-91-6	8	1 µg/mL
2-(N-Methylperfluorooctanesulfonamido)acetic acid		N-MeFOSAA	2355-31-9	8	1 µg/mL
2-(N-Ethylperfluorooctanesulfonamido)acetic acid		N-EtFOSAA	2991-50-6	8	1 µg/mL
Perfluorooctane sulfonic acid		PFOS	1763-23-1	8	1 µg/mL
1H,1H,2H,2H-Perfluorooctane sulfonic acid		6:2FTS	27619-97-2	8	4 µg/mL
Perfluorononanoic acid		PFNA	375-95-1	9	1 µg/mL
Perfluorononanesulfonic acid		PFNS	68259-12-1	9	1 µg/mL
Perfluorodecanoic acid		PFDA	335-76-2	10	1 µg/mL
Perfluorodecanesulfonic acid		PFDS	335-77-3	10	1 µg/mL
Perfluoroundecanoic acid		PFUnDA	2058-94-8	11	1 µg/mL
Perfluorododecanoic acid		PFDoDA	307-55-1	12	1 µg/mL
Perfluorododecanesulfonic acid		PFDoDS	79780-39-5	12	1 µg/mL

PFAS reference materials

Product Listings

Mixtures. Continued from previous page

Analytes	Abbreviation	CAS No.	Carbon Chain Length	*Concentration
1H,1H,2H,2H-Perfluorododecane sulfonic acid	10:2FTS	120226-60-0	12	4 µg/mL
Perfluorotridecanoic acid	PFTTrDA	72629-94-8	13	1 µg/mL
Perfluorotetradecanoic acid	PFTeDA	376-06-7	14	1 µg/mL

Product Number	Product Description		Pack size	Analytes	
DRE-A30000071MW	PFAS Superfund Mixture 2 µg/mL in Methanol:Water		1ml	24 analytes	
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluorobutanesulfonic acid		PFBS	375-73-5	4	2 µg/mL
Perfluorobutanoic acid		PFBA	375-22-4	4	2 µg/mL
Perfluoropentanoic acid		PFPeA	2706-90-3	5	2 µg/mL
Perfluorohexanesulfonic acid		PFHxS	355-46-4	6	2 µg/mL
1H,1H,2H,2H-Perfluorohexanesulfonic acid		4:2 FTS	757124-72-4	6	2 µg/mL
Perfluorohexanoic acid		PFHxA	307-24-4	6	2 µg/mL
Perfluoro-2-propoxypropanoic acid		PFPPrOPrA	13252-13-6	6	2 µg/mL
Perfluoroheptanoic acid		PFHpA	375-85-9	7	2 µg/mL
N-Methylperfluorooctanesulfonamide		MeFOSA	31506-32-8	8	2 µg/mL
2-(N-Methylperfluorooctanesulfonamido)acetic acid		N-MeFOSAA	2355-31-9	8	2 µg/mL
2-(N-Ethylperfluorooctanesulfonamido)acetic acid		N-EtFOSAA	2991-50-6	8	2 µg/mL
Perfluorooctane sulfonamide		PFOSA	754-91-6	8	2 µg/mL
Perfluorooctane sulfonic acid		PFOS	1763-23-1	8	2 µg/mL
1H,1H,2H,2H-Perfluorooctane sulfonic acid		6:2 FTS	27619-97-2	8	2 µg/mL
Perfluorooctanoic acid		PFOA	335-67-1	8	2 µg/mL
Perfluorononanoic acid		PFNA	375-95-1	9	2 µg/mL
1H,1H,2H,2H-Perfluorodecanesulfonic acid		8:2 FTS	39108-34-4	10	2 µg/mL
Perfluoro-n-decanoic acid		PFDA	335-76-2	10	2 µg/mL
Perfluoroundecanoic acid		PFUnDA	2058-94-8	11	2 µg/mL
Perfluoro-n-dodecanoic acid		PFDoDA	307-55-1	12	2 µg/mL
1H,1H,2H,2H-Perfluorododecane sulfonic acid		10:2 FTS	120226-60-0	12	2 µg/mL
Perfluorotetradecanoic acid		PFTeDA	376-06-7	14	2 µg/mL
Perfluorohexadecanoic acid		PFHxDA	67905-19-5	16	2 µg/mL
Perfluorooctadecanoic acid		PFODA	16517-11-6	18	2 µg/mL

Product Number	Product Description		Pack size	Analytes	
DRE-A50000151MW	PFAS Mixture 151 100 µg/mL in Methanol:Water		1ml	24 analytes	
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluorobutanesulfonic acid potassium salt		PFBSK	29420-49-3	4	100 µg/mL
Perfluorobutanoic acid		PFBA	375-22-4	4	100 µg/mL
Perfluoropentanesulfonic acid sodium salt		PFPeS	630402-22-1	5	100 µg/mL
Perfluoropentanoic acid		PFPeA	2706-90-3	5	100 µg/mL
1,1,2,2H-Perfluorohexanesulfonic acid		4:2 FTS	757124-72-4	6	100 µg/mL
Perfluorohexanesulfonic acid potassium salt		PFHxS	3871-99-6	6	100 µg/mL
Perfluorohexanoic acid		PFHxA	307-24-4	6	100 µg/mL
Perfluoroheptanesulfonic acid		PFHpS	375-92-8	7	100 µg/mL

*Concentrations are always given for the free acids, which is also stated on our CoA

lgcstandards.com

Continued on next page

PFAS reference materials

Product Listings

Mixtures. Continued from previous page

Analytes	Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluoroheptanoic acid	PFHpA	375-85-9	7	100 µg/mL
1,1,2,2H-Perfluorooctanesulfonic acid	6:2 FTS	27619-97-2	8	100 µg/mL
2-(N-Ethylperfluorooctanesulfonamido)acetic acid	N-EtFOSAA	2991-50-6	8	100 µg/mL
2-(N-Methylperfluorooctanesulfonamido)acetic acid	N-MeFOSAA	2355-31-9	8	100 µg/mL
Perfluorooctane sulfonamide PFOSA	PFOSA	754-91-6	8	100 µg/mL
<i>Perfluorooctane sulfonic acid potassium salt (PFOS)</i>	PFOS	2795-39-3	8	100 µg/mL
Perfluorooctanoic acid	PFOA	335-67-1	8	100 µg/mL
<i>Perfluorononanesulfonic acid sodium salt</i>	PFNS	98789-57-2	9	100 µg/mL
Perfluorononanoic acid	PFNA	375-95-1	9	100 µg/mL
1,1,2,2H-Perfluorodecanesulfonic acid	8:2 FTS	39108-34-4	10	100 µg/mL
<i>Perfluorodecanesulfonic acid sodium salt</i>	PFDS	2806-15-7	10	100 µg/mL
Perfluorodecanoic acid	PFDA	335-76-2	10	100 µg/mL
Perfluoroundecanoic acid	PFUnDA	2058-94-8	11	100 µg/mL
Perfluorododecanoic acid	PFDODA	307-55-1	12	100 µg/mL
Perfluorotridecanoic acid	PFTTrDA	72629-94-8	13	100 µg/mL
Perfluorotetradecanoic acid	PFTeDA	376-06-7	14	100 µg/mL

Product Number	Product Description		Pack size	Analytes	
DRE-GS09000218MW	PFAA Mixture 218 100 µg/mL in Methanol:Water		1ml	11 analytes	
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluoro-1-butanedisulfonic Acid		PFBS	375-73-5	4	100 µg/mL
Perfluorohexanoic Acid		PFHxA	307-24-4	6	100 µg/mL
Perfluoroheptanoic Acid		PFHpA	375-85-9	7	100 µg/mL
Pentadecafluorooctanoic Acid Hydrate		PFOA	335-67-1	8	100 µg/mL
Perfluorooctanesulfonic Acid		PFOS	1763-23-1	8	100 µg/mL
Perfluorononanoic Acid		PFNA	375-95-1	9	100 µg/mL
Perfluorodecanoic Acid		PFDA	335-76-2	10	100 µg/mL
Perfluoroundecanoic Acid		PFUnDA	2058-94-8	11	100 µg/mL
Perfluorododecanoic Acid		PFDODA	307-55-1	12	100 µg/mL
Perfluorotridecanoic Acid		PFTTrDA	72629-94-8	13	100 µg/mL
Perfluorotetradecanoic Acid		PFTeDA	376-06-7	14	100 µg/mL

Product Number	Product Description		Pack size	Analytes	
DRE-A30000072MW	PFCA (C4-C14) Native Mixture 2 µg/mL in Methanol:Water		1ml	11 analytes	
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluorobutyric acid		PFBA	375-22-4	4	2 µg/mL
Perfluoropentanoic acid		PFPeA	2706-90-3	5	2 µg/mL
Perfluorohexanoic acid		PFHxA	307-24-4	6	2 µg/mL
Perfluoroheptanoic acid		PFHpA	375-85-9	7	2 µg/mL
Perfluoro-n-octanoic acid		PFOA	335-67-1	8	2 µg/mL
Perfluorononanoic acid		PFNA	375-95-1	9	2 µg/mL
Perfluorodecanoic acid		PFDA	335-76-2	10	2 µg/mL
Perfluoroundecanoic acid		PFUA	2058-94-8	11	2 µg/mL

PFAS reference materials

Product Listings

Mixtures. Continued from previous page

Analytes	Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluorododecanoic acid	PFDODA	307-55-1	12	2 µg/mL
Perfluorotridecanoic acid	PFTTrDA	72629-94-8	13	2 µg/mL
Perfluorotetradecanoic acid	PFTeDA	376-06-7	14	2 µg/mL

Product Number	Product Description		Pack size	Analytes	
DRE-A30000064MW	Ultrashort-chain-PFAS-7 (C1-C4) Mixture 50 µg/mL in Methanol:Water		1ml	7 analytes	
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluoromethanesulfonic acid		TFMS	1493-13-6	1	50 µg/mL
Perfluoroethanesulfonic acid		PFEtS	354-88-1	2	50 µg/mL
Trifluoroacetic acid		-	76-05-1	2	50 µg/mL
Perfluoropropanesulfonic acid		PFPrS	423-41-6	3	50 µg/mL
Perfluoropropionic acid		PFPrA	422-64-0	3	50 µg/mL
Perfluorobutanoic acid		PFBA	375-22-4	4	50 µg/mL
Perfluorobutanesulfonic acid		PFBS	375-73-5	4	50 µg/mL

Product Number	Product Description		Pack size	Analytes	
DRE-A30000082MW	Ultrashort-chain-PFAS-10 (C1-C4) Mixture 50 µg/mL in Methanol:Water		1ml	10 analytes	
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluoromethanesulfonic acid		TFMS	1493-13-6	1	50 µg/mL
Perfluoroethanesulfonic acid		PFEtS	354-88-1	2	50 µg/mL
Trifluoroacetic acid		-	76-05-1	2	50 µg/mL
Perfluoropropanesulfonic acid sodium		PFPrS	359868-82-9	3	50 µg/mL
2H-Perfluoropropionic Acid		-	359-49-9	3	50 µg/mL
Perfluoropropionic acid		PFPrA	422-64-0	3	50 µg/mL
Perfluoro-2-methoxyacetic acid		PFMOAA	674-13-5	3	50 µg/mL
Flupropanate		-	756-09-2	3	50 µg/mL
Perfluorobutanoic acid		PFBA	375-22-4	4	50 µg/mL
Perfluorobutanesulfonic acid		PFBS	375-73-5	4	50 µg/mL

Product Number	Product Description		Pack size	Analytes	
DRE-A50000761MW	PFAS Mixture (DONA, 9CI-PF3ONS, 11CI-PF3OUdS, HFPO-DA) 2 µg/mL in Methanol:Water		1ml	4 analytes	
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluoro-2-propoxypropanoic acid		PFPrOPrA	13252-13-6	6	2 µg/mL
3H-Perfluoro-4,8-dioxanonanoic acid		DONA	919005-14-4	7	2 µg/mL
9-Chlorohexadecafluoro-3-oxanonane-1-sulfonic acid		9CI-PF3ONS	756426-58-1	8	2 µg/mL
11-Chloroeicosafluoro-3-oxaundecane-1-sulfonic acid		11CI-PF3OUdS	763051-92-9	10	2 µg/mL

PFAS reference materials

Product Listings

Mixtures. Continued from previous page

Product Number	Product Description		Pack size	Analytes	
DRE-A50000762MW	PFAS Mixture (PF4OPeA, PF5OHxA, 3,6-OPFHpA, PFEESA) 2 µg/mL in Methanol:Water		1ml	4 analytes	
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluoro(2-ethoxyethane) sulfonic acid		PFEESA	113507-82-7	4	2 µg/mL
Perfluoro-3-methoxypropanoic acid		PFMOPrA	377-73-1	4	2 µg/mL
Perfluoro-4-methoxybutanoic acid		PFMOBA	863090-89-5	5	2 µg/mL
Perfluoro-3,6-dioxaheptanoic acid		NFDHA	151772-58-6	5	2 µg/mL

Product Number	Product Description		Pack size	Analytes	
DRE-A50000763MW	PFAS Mixture (N-MeFOSA, N-EtFOSA, N-MeFOSE, N-EtFOSE) 1-10 µg/mL in Methanol:Water		1ml	4 analytes	
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
N-Methylperfluorooctanesulfonamide		MeFOSA	31506-32-8	8	1 µg/mL
Sulfluramid		EtFOSA	4151-50-2	8	1 µg/mL
N-(2-Hydroxyethyl)-N-methylperfluorooctanesulfonamide		N-MeFOSE	24448-09-7	8	10 µg/mL
N-Ethyl-N-(2-hydroxyethyl)perfluorooctanesulfonamide		N-EtFOSE	1691-99-2	8	10 µg/mL

Product Number	Product Description			Pack size	Analytes
DRE-A50000759 MW	PFAS Native X:3 FTCAs Mixture 4-20 µg/mL in Methanol:Water			1ml	3 analytes
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
2H,2H,3H,3H-Perfluorohexanoic acid		3:3 FTCA	356-02-5	6	4 µg/mL
2H,2H,3H,3H-Perfluorooctanoic acid		5:3 FTCA	914637-49-3	8	20 µg/mL
2H,2H,3H,3H-Perfluorodecanoic acid		7:3 FTCA	812-70-4	10	20 µg/mL

Product Number	Product Description		Pack size	Analytes	
DRE-A30000088MW	Mass-Labelled PFC Injection Standards Solution 2 µg/mL in Methanol:Water		1ml	4	
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluoro-n-(2,3,4-13C3)butanoic acid		13C3-PFBA	2483735-33-5	4	2 µg/mL
Perfluoro-n-(1,2-13C2)octanoic acid		13C2-PFOA	864071-08-9	8	2 µg/mL
Sodium perfluoro-1-(1,2,3,4-13C4)octanesulfonate		13C4-PFOS	960315-53-1	8	2 µg/mL
Perfluoro-n-(1,2-13C2)decanoic acid		13C2-PFDA	960315-50-8	10	2 µg/mL

Product Number	Product Description			Pack size	Analytes
DRE-A30000089MW	PFAS-C-ES Mixture 2 µg/mL in Methanol:Water			1ml	13
Analytes		Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluoro-n-(13C4)butanoic acid		13C4-PFBA	1017281-29-6	4	2 µg/mL
Sodium perfluoro-1-(2,3,4-13C3)butanesulfonate		13C3-PFBS	2708218-84-0	4	2 µg/mL
Perfluoro-n-(13C5)pentanoic acid		13C5-PFPeA	2283397-79-3	5	2 µg/mL
Perfluoro-n-(1,2,3,4,6-13C5)hexanoic acid		13C5-PFHxA	2328024-54-8	6	2 µg/mL
Sodium perfluoro-1-(1,2,3-13C3)hexanesulfonate		13C3-PFHxS	2708218-86-2	6	2 µg/mL
Perfluoro-n-(1,2,3,4-13C4)heptanoic acid		13C4-PFHpA	2328024-55-9	7	2 µg/mL

PFAS reference materials

Product Listings

Mixtures. Continued from previous page

Analytes	Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluoro-n-(13C8)octanoic acid	13C8-PFOA	1350614-84-4	8	2 µg/mL
Sodium perfluoro-1-(13C8)octanesulfonate	13C8-PFOS	2522762-16-7	8	2 µg/mL
Perfluoro-n-(13C9)nonanoic acid	13C9-PFNA	2283397-80-6	9	2 µg/mL
Perfluoro-n-(1,2,3,4,5,6-13C6)decanoic acid	13C6-PFDA	2328024-56-0	10	2 µg/mL
Perfluoro-n-(1,2,3,4,5,6,7-13C7)undecanoic acid	13C7-PFUnDA	-	11	2 µg/mL
Perfluoro-n-(1,2-13C2)dodecanoic acid	13C2-PFDoDA	960315-52-0	12	2 µg/mL
Perfluoro-n-(1,2-13C2)tetradecanoic acid	13C2-PFTeDA	2708218-82-8	14	2 µg/mL

Product Number	Product Description	Pack size	Analytes
DRE-A50000793MW	Mass-Labelled PFAS-24ES Extraction Mixture 1 µg/mL in Methanol:Water	1ml	19

Analytes	Abbreviation	CAS No.	Carbon Chain Length	*Concentration
Perfluoro-n-(13C4)butanoic acid	13C4-PFBA	1017281-29-6	4	1 µg/mL
Sodium perfluoro-1-(2,3,4-13C3)butanesulfonate	13C3-PFBS	2708218-84-0	4	1 µg/mL
Perfluoro-n-(13C5)pentanoic acid	13C5-PFPeA	2283397-79-3	5	1 µg/mL
Perfluoro-n-(1,2,3,4,6-13C5)hexanoic acid	13C5-PFHxA	2328024-54-8	6	1 µg/mL
Sodium perfluoro-1-(1,2,3-13C3)hexanesulfonate	13C3-PFHxS	2708218-86-2	6	1 µg/mL
Sodium 1H,1H,2H,2H-perfluoro(1,2-13C2)hexanesulfonate	13C2-4:2 FTS	2708218-88-4	6	1 µg/mL
Perfluoro-n-(1,2,3,4-13C4)heptanoic acid	13C4-PFHpA	2328024-55-9	7	1 µg/mL
Perfluoro-n-(13C8)octanoic acid	13C8-PFOA	1350614-84-4	8	1 µg/mL
Perfluoro-1-(13C8)octanesulfonamide	13C8-FOSA	1365803-60-6	8	1 µg/mL
Sodium perfluoro-1-(13C8)octanesulfonate	13C8-PFOS	2522762-16-7	8	1 µg/mL
Sodium 1H,1H,2H,2H-perfluoro(1,2-13C2)octanesulfonate	13C2-6:2 FTS	2708218-89-5	8	1 µg/mL
N-methyl-d3 -perfluoro-1-octanesulfonamidoacetic acid	D3-N-MeFOSAA	1400690-70-1	8	1 µg/mL
N-ethyl-d5 -perfluoro-1-octanesulfonamidoacetic acid	D5-N-EtFOSAA	1265205-97-7	8	1 µg/mL
Perfluoro-n-(13C9)nonanoic acid	13C9-PFNA	2283397-80-6	9	1 µg/mL
Perfluoro-n-(1,2,3,4,5,6-13C6)decanoic acid	13C6-PFDA	2328024-56-0	10	1 µg/mL
Sodium 1H,1H,2H,2H-perfluoro(1,2-13C2)decanesulfonate	13C2-8:2 FTS	2708218-90-8	10	1 µg/mL
Perfluoro-n-(1,2,3,4,5,6,7-13C7)undecanoic acid	13C7-PFUnDA	-	11	1 µg/mL
Perfluoro-n-(1,2-13C2)dodecanoic acid	13C2-PFDoA	960315-52-0	12	1 µg/mL
Perfluoro-n-(1,2-13C2)tetradecanoic acid	13C2-PFTeDA	2708218-82-8	14	1 µg/mL

Stable Isotopes

Our range of PFAS reference materials include ¹³C stable isotope labelled PFAS compounds to support regulatory testing requirements and improved analytical accuracy, including when working with complex matrices. Available in single solution format and accredited to ISO 17034 to provide you with confidence in your results.

Product Number	Product Description	Pack Size	Abbreviation	Labelled CAS Number	Unlabelled CAS Number	Carbon Chain Length
Labelled Carboxylic Acids						
<i>DRE-C17844405</i>	Trifluoroacetic acid ¹³ C2 sodium	10 mg	13C2-TFA	1794767-05-7	76-05-1	2
<i>DRE-A15987241MW-50</i>	Perfluoropropionic acid ¹³ C3 50 µg/mL in Methanol:Water	1ml	13C3-PFPrA	-	422-64-0	3
<i>DRE-A15987253MW-50</i>	Perfluoro-2-propoxypropanoic acid ¹³ C3 (propanoic acid ¹³ C3) 50 µg/mL in Methanol:Water	1ml	13C3-PFPrOPrA	3030247-97-0	13252-13-6	3/3
<i>DRE-A15986523MW-50</i>	Perfluorobutanoic acid ¹³ C3 (2,3,4- ¹³ C3) 50 µg/mL in Methanol:Water	1ml	13C3-PFBA	2483735-33-5	375-22-4	4
<i>DRE-A15986521MW-50</i>	Perfluorobutanoic acid ¹³ C4 50 µg/mL in Methanol:Water	1ml	13C4-PFBA	1017281-29-6	375-22-4	4
<i>DRE-A15987201MW-50</i>	Perfluoropentanoic acid ¹³ C5 50 µg/mL in Methanol:Water	1ml	13C5-PFPeA	2283397-79-3	2706-90-3	5
<i>DRE-A15986951MW-50</i>	Perfluoro-4-methoxybutanoic acid ¹³ C2 (1,2- ¹³ C2) 50 µg/mL in Methanol:Water	1ml	13C2-PFMBA	-	863090-89-5	5
<i>DRE-A15986914MW-50</i>	Perfluorohexanoic acid ¹³ C2 (1,2- ¹³ C2) 50 µg/mL in Methanol:Water	1ml	13C2-PFHxA	960315-47-3	307-24-4	6
<i>DRE-A15986911MW-50</i>	Perfluorohexanoic acid ¹³ C5 (1,2,3,4,6- ¹³ C5) 50 µg/mL in Methanol:Water	1ml	13C5-PFHxA	2328024-54-8	307-24-4	6
<i>DRE-A15986891MW-50</i>	Perfluoroheptanoic acid ¹³ C4 (1,2,3,4- ¹³ C4) 50 µg/mL in Methanol:Water	1ml	13C4-PFHpA	2328024-55-9	375-85-9	7
<i>DRE-A15986875MW-50</i>	Perfluoroheptanoic acid ¹³ C5 (1,2,3,4,5- ¹³ C5) 50 µg/mL in Methanol:Water	1ml	13C5-PFHpA	2328024-55-9	375-85-9	7
<i>DRE-A11041320MW-50</i>	Capstone product B ¹³ C3 (dimethyl ¹³ C2, glycine-beta- ¹³ C1) 50 µg/mL in Methanol:Water	1ml	13C4-PFOA	-	34455-29-3	8
<i>DRE-A15987154MW-50</i>	Perfluorooctanoic acid ¹³ C4 (1,2,3,4- ¹³ C4) 50 µg/mL in Methanol:Water	1ml	13C4-PFOA	960315-48-4	335-67-1	8
<i>DRE-A15987149MW-50</i>	Perfluorooctanoic acid ¹³ C2 (1,2- ¹³ C2) sodium 50 µg/mL in Methanol:Water	1ml	13C2-PFOA	-	335-95-5	8
<i>DRE-A15987151MW-50</i>	Perfluorooctanoic acid ¹³ C8 50 µg/mL in Methanol:Water	1ml	13C8-PFOA	1350614-84-4	335-67-1	8
<i>DRE-A15987001MW-50</i>	Perfluorononanoic acid ¹³ C9 50 µg/mL in Methanol:Water	1ml	13C9-PFNA	2283397-80-6	375-95-1	9
<i>DRE-A15987002MW-50</i>	Perfluorononanoic acid ¹³ C5 (1,2,3,4,5- ¹³ C5) 50 µg/mL in Methanol:Water	1ml	13C5-PFNA	960315-49-5	375-95-1	9
<i>DRE-A15986599MW-50</i>	Perfluorodecanoic acid ¹³ C6 (1,2,3,4,5,6- ¹³ C6) 50 µg/mL in Methanol:Water	1ml	13C6-PFDA	2328024-56-0	335-76-2	10
<i>DRE-A15986596MW-50</i>	Perfluorodecanoic acid ¹³ C2 (1,2- ¹³ C2) 50 µg/mL in Methanol:Water	1ml	13C2-PFDA	960315-50-8	335-76-2	10
<i>DRE-A15989002MW-50</i>	Perfluoroundecanoic acid ¹³ C2 (1,2- ¹³ C2) 50 µg/mL in Methanol:Water	1ml	13C2-PFUnDA	960315-51-9	2058-94-8	11
<i>DRE-A15989001MW-50</i>	Perfluoroundecanoic acid ¹³ C7 (1,2,3,4,5,6,7- ¹³ C7) 50 µg/mL in Methanol:Water	1ml	13C7-PFUnDA	-	2058-94-8	11
<i>DRE-A15986628MW-50</i>	Perfluorododecanoic acid ¹³ C2 (1,2- ¹³ C2) 50 µg/mL in Methanol:Water	1ml	13C2-PFDoDA	960315-52-0	307-55-1	12
<i>DRE-A15988003MW-50</i>	Perfluorotridecanoic acid ¹³ C2 (1,2- ¹³ C2) 50 µg/mL in Methanol:Water	1ml	13C2-PFTrDA	-	72629-94-8	13
<i>DRE-A15987410MW-50</i>	Perfluorotetradecanoic acid ¹³ C2 (1,2- ¹³ C2) 50 µg/mL in Methanol:Water	1ml	13C2-PFTeDA	2708218-82-8	376-06-7	14
<i>DRE-A15986899MW-50</i>	Perfluorohexadecanoic acid ¹³ C2 (1,2- ¹³ C2) 50 µg/mL in Methanol:Water	1ml	-	-	67905-19-5	16
Labelled Sulfonic Acids - Sulfonates & Salts						
<i>DRE-A15986510MW-50</i>	Perfluorobutanesulfonic acid ¹³ C4 sodium 50 µg/mL in Methanol:Water	1ml	13C4-PFBS	-	29420-49-3	4
<i>DRE-A15986518MW-50</i>	Perfluorobutanesulfonic acid ¹³ C3 (2,3,4- ¹³ C3) sodium 50 µg/mL in Methanol:Water	1ml	13C3-PFBS	2708218-84-0	375-73-5	4
<i>DRE-A15987207MW-50</i>	Perfluoropentanesulfonic acid ¹³ C3 (1,2,3- ¹³ C3) sodium 50 µg/mL in Methanol:Water	1ml	13C3-PFPeS	-	2732349-38-9	5
<i>DRE-A15986906MW-50</i>	Perfluorohexanesulfonic acid ¹³ C3 (1,2,3- ¹³ C3) sodium 50 µg/mL in Methanol:Water	1ml	13C3-PFHxS	2708218-86-2	355-46-4	6
<i>DRE-A15986902MW-50</i>	Perfluorohexanesulfonic acid 18O2 sodium 50 µg/mL in Methanol:Water	1ml	18O2-PFHxS	-	82382-12-5	6
<i>DRE-A15986629MW-50</i>	1H,1H,2H,2H-Perfluorohexanesulfonic acid ¹³ C2 (1,2- ¹³ C2) sodium 50 µg/mL in Methanol:Water	1ml	13C2-4:2 FTS	2708218-88-4	27619-93-8	6
<i>DRE-A15986883MW-50</i>	Perfluoroheptanesulfonic acid ¹³ C3 (1,2,3- ¹³ C3) sodium 50 µg/mL in Methanol:Water	1ml	13C3-PFHpS	-	21934-50-9	7

PFAS reference materials

Product Listings

Stable Isotopes. Continued from previous page

Product Number	Product Description	Pack Size	Abbreviation	Labelled CAS Number	Unlabelled CAS Number	Carbon Chain Length
DRE-A15987123MW-50	Perfluorooctanesulfonic acid 13C4 (1,2,3,4-13C4) sodium 50 µg/mL in Methanol:Water	1ml	13C4-PFOS	960315-53-1	4021-47-0	8
DRE-A15987124MW-50	Perfluorooctanesulfonic acid 13C8 sodium 50 µg/mL in Methanol:Water	1ml	13C8-PFOS	2522762-16-7	4021-47-0	8
DRE-A15987127MW-50	1H,1H,2H,2H-Perfluorooctanesulfonic acid 13C2 (1,2-13C2) sodium 50 µg/mL in Methanol:Water	1ml	13C2-6:2 FTS	2708218-89-5	27619-94-9	8
DRE-A15987023MW-50	Perfluorononanesulfonic acid 13C3(1,2,3-13C3) sodium 50 µg/mL in Methanol:Water	1ml	13C3-PFNS	-	98789-57-2	9
DRE-A15986582MW-50	Perfluorodecanesulfonic acid 13C3(1,2,3-13C3) sodium 50 µg/mL in Methanol:Water	1ml	13C3-PFDS	-	2806-15-7	10
DRE-A15986587MW-50	1H,1H,2H,2H-Perfluorodecanesulfonic acid 13C2 (1,2-13C2) sodium 50 µg/mL in Methanol:Water	1ml	13C2-8:2 FTS	2708218-90-8	27619-96-1	10
DRE-A15988996MW-50	Perfluoroundecanesulfonic acid 13C3(1,2,3-13C3) sodium 50 µg/mL in Methanol:Water	1ml	13C3-PFUnDS	-	3016286-99-7	11
DRE-A15986645MW-50	Perfluorodecanesulfonic acid 13C3 (1,2,3-13C3) sodium 50 µg/mL in Methanol:Water	1ml	13C3-PFDoDS	-	79780-39-5	12
DRE-A15986633MW-50	1H,1H,2H,2H-Perfluorodecanesulfonic acid 13C2 D4 (1,2-13C2 D4) sodium 50 µg/mL in Methanol:Water	1ml	13C2-D4-10:2 FTS	-	120226-60-0	12
DRE-A15987996MW-50	Perfluorotridecanesulfonic acid 13C3 (1,2,3-13C3) sodium 50 µg/mL in Methanol:Water	1ml	13C3-PFTrDS	-	174675-49-1	13
Labelled Sulfonamide						
DRE-A15130020MW-50	2-(N-Methylperfluorooctanesulfonamido)acetic acid 13C3 50 µg/mL in Methanol:Water	1ml	13C3-N-MeFOSAA	-	2355-31-9	8
DRE-A13349620MW-50	2-(N-Ethylperfluorooctanesulfonamido)acetic acid 13C4 50 µg/mL in Methanol:Water	1ml	13C4-N-EtFOSAA	-	2991-50-6	8
DRE-C15986909	Perfluorohexanoic acid sodium	100mg	PFHxA	2923-26-4	6	
DRE-A15130010MW-50	2-(N-Methylperfluorooctanesulfonamido)acetic acid D3 (N-methyl D3) 50 µg/mL in Methanol:Water	1ml	D3-N-MeFOSAA	1400690-70-1	2355-31-9	8
DRE-A15115501MW-50	N-Methylperfluorooctanesulfonamide D3 (methyl D3) 50 µg/mL in Methanol:Water	1ml	D3-N-MeFOSA	936109-37-4	31506-32-8	8
DRE-A17004001MW-50	Sulfluramid D5 (ethyl D5) 50 µg/mL in Methanol:Water	1ml	D5-N-EtFOSA	936109-40-9	4151-50-2	8
DRE-A14231571MW-50	N-(2-Hydroxyethyl)-N-methylperfluorooctanesulfonamide D7 50 µg/mL in Methanol:Water	1ml	D7-N-MeFOSE	1265205-95-5	24448-09-7	8
DRE-A13342361MW-50	N-Ethyl-N-(2-hydroxyethyl)perfluorooctanesulfonamide D9 50 µg/mL in Methanol:Water	1ml	D9-N-EtFOSE	1265205-96-6	1691-99-2	8
DRE-A13349601MW-50	2-(N-Ethylperfluorooctanesulfonamido)acetic acid D5 (N-ethyl D5) 50 µg/mL in Methanol:Water	1ml	D5-N-EtFOSAA	1265205-97-7	2991-50-6	8
DRE-A15987112MW-50	Perfluorooctane sulfonamide 18O2 50 µg/mL in Methanol:Water	1ml	18O2-PFOSA	1340480-46-7	754-91-6	8
Labeled Polyfluorinated Phosphate Esters						
DRE-A10655199MW-50	Bis[2-(perfluorohexyl)ethyl] phosphate 13C4 (1,1',2,2'-13C4) sodium 50 µg/mL in Methanol:Water	1ml	13C4-6:2 Di-PAPS	2179232-14-3	57677-95-9	8/8
Mixtures						
DRE-A30000084MW	Korea cosmetics testing Labelled PFAS Mixture 1 µg/mL in Methanol:Water	1ml	-	-	-	-
DRE-A50000795MW	EPA Method 533 Isotope Performance Standard 1-3 µg/mL in Methanol:Water	1ml	-	-	-	-
DRE-A30000089MW	PFAS-C-ES Mixture 2 µg/mL in Methanol:Water	1ml	-	-	-	-
DRE-A50000792MW	EPA Method 1633 EIS Labelled PFAS Mixture 0.25-5 µg/mL in Methanol:Water	1ml	-	-	-	-
DRE-A50000794MW	EPA Method 1633 NIS Labelled PFAS Mixture 0.25-1 µg/mL in Methanol:Water	1ml	-	-	-	-
DRE-A50000793MW	Mass-Labelled PFAS-24ES Extraction Mixture 1 µg/mL in Methanol:Water	1ml	-	-	-	-
DRE-A50000790MW	EPA Method 533 Isotope Dilution Standard 0.5-2 µg/mL in Methanol:Water	1ml	-	-	-	-
DRE-A50000791MW	EPA Method 537.1 Surrogate Primary Dilution Standard 1-4 µg/mL in Methanol:Water	1ml	-	-	-	-
DRE-A30000095MW	EU Drinking Water Directive Labelled PFAS Mixture 2 µg/mL in Methanol:Water	1ml	-	-	-	-
DRE-A30000088MW	Mass-Labelled PFC Injection Standards Solution 2 µg/mL in Methanol:Water	1ml	-	-	-	-
DRE-A50000796MW	EPA Method 537 Internal Standard Primary Dilution Standard 1-4 µg/mL in Methanol:Water	1ml	-	-	-	-

*Concentrations are always given for the free acids, which is also stated on our CoA

Per- and Polyfluorocarboxylic Acids

Per- and Polyfluorocarboxylic Acids (PFCA) are organofluorine analogues of ordinary carboxylic acids. They exhibit hydrophobic properties and are widely used in the synthesis of peptides. Longer chain perfluorocarboxylic acids such as Perfluorooctanoic Acid (PFOA) are used in firefighting foams and are under increased scrutiny due to their bioaccumulative properties.

Product Number	Product Description	Pack Size	Abbreviation	CAS Number	Carbon Chain Length
DRE-A17844400MW-50	Trifluoroacetic acid 50 µg/mL in Methanol:Water	1ml	TFA	76-05-1	2
DRE-C17844415	Trifluoroacetic acid sodium	250mg	TFA	2923-18-4	2
DRE-C17844405	Trifluoroacetic acid 13C2 sodium	10 mg	13C2-TFA	1794767-05-7	2
DRE-CA13802000	Flupropionate	50mg	-	756-09-2	3
DRE-A13802000MW-50	Flupropionate 50 µg/mL in Methanol:Water	1ml	-	756-09-2	3
DRE-CA15986940	Perfluoro-2-methoxyacetic acid (PFMOAA)	25mg	PFMOAA	674-13-5	3
DRE-CA15987240	Perfluoropropionic acid	500mg	PFPrA	422-64-0	3
DRE-A15987240MW-50	Perfluoropropionic acid 50 µg/mL in Methanol:Water	1ml	PFPrA	422-64-0	3
DRE-CA15987242	2H-Perfluoropropionic acid	50mg	-	359-49-9	3
DRE-A15987242MW-50	2H-Perfluoropropionic acid 50 µg/mL in Methanol:Water	1ml	-	359-49-9	3
DRE-A15987241MW-50	Perfluoropropionic acid 13C3 50 µg/mL in Methanol:Water	1ml	13C3-PFPrA	-	3
DRE-A15987253MW-50	Perfluoro-2-propoxypropanoic acid 13C3 (propanoic acid 13C3) 50 µg/mL in Methanol:Water	1ml	13C3-PFPrOPrA	3030247-97-0	3/3
DRE-C15986520	Perfluorobutanoic acid	100mg	PFBA	375-22-4	4
DRE-A15986520MW-50	Perfluorobutanoic acid 50 µg/mL in Methanol:Water	1ml	PFBA	375-22-4	4
DRE-A15986520AL-100	Perfluorobutanoic acid 100 µg/mL in Acetonitrile	1ml	PFBA	375-22-4	4
DRE-A15986523MW-50	Perfluorobutanoic acid 13C3 (2,3,4-13C3) 50 µg/mL in Methanol:Water	1ml	13C3-PFBA	2483735-33-5	4
DRE-A15986521MW-50	Perfluorobutanoic acid 13C4 50 µg/mL in Methanol:Water	1ml	13C4-PFBA	1017281-29-6	4
DRE-A15986612MW-50	Perfluoro-3,6-dioxahexanoic acid 50 µg/mL in Methanol:Water	1ml	NFDHA	151772-58-6	5
DRE-C15987200	Perfluoropentanoic acid	100mg	PFPeA	2706-90-3	5
DRE-A15987200MW-50	Perfluoro-n-pentanoic acid 50 µg/mL in Methanol:Water	1ml	PFPeA	2706-90-3	5
DRE-A15987200MW-100	Perfluoropentanoic acid 100 µg/mL in Methanol:Water	1ml	PFPeA	2706-90-3	5
DRE-A15987201MW-50	Perfluoropentanoic acid 13C5 50 µg/mL in Methanol:Water	1ml	13C5-PFPeA	2283397-79-3	5
DRE-A15986951MW-50	Perfluoro-4-methoxybutanoic acid 13C2 (1,2-13C2) 50 µg/mL in Methanol:Water	1ml	13C2-PFMBA	-	5
DRE-C15986910	Perfluorohexanoic acid	100mg	PFHxA	307-24-4	6
DRE-A15986910AL-100	Perfluorohexanoic acid 100 µg/mL in Acetonitrile	1ml	PFHxA	307-24-4	6
DRE-C15986868	2H,2H-Perfluorohexanoic acid	10mg	-	70887-89-7	6
DRE-C15986912	2H,2H,3H,3H-Perfluorohexanoic acid	50mg	3:3 FTCA	356-02-5	6
DRE-A15986980MW-50	Perfluoro-2-methyl-3-oxahexanoic acid ammonium 50 µg/mL in Methanol:Water	1ml	-	62037-80-3	6
DRE-A15986909MW-50	Perfluorohexanoic acid sodium 50 µg/mL in Methanol:Water	1ml	PFHxA	2923-26-4	6
DRE-A15986914MW-50	Perfluorohexanoic acid 13C2 (1,2-13C2) 50 µg/mL in Methanol:Water	1ml	13C2-PFHxA	960315-47-3	6
DRE-A15986911MW-50	Perfluorohexanoic acid 13C5 (1,2,3,4,6-13C5) 50 µg/mL in Methanol:Water	1ml	13C5-PFHxA	2328024-54-8	6
DRE-A15986618MW-50	3H-Perfluoro-4,8-dioxanonoic acid 50 µg/mL in Methanol:Water	1ml	DONA	919005-14-4	7
DRE-C15986892	7H-Perfluoroheptanoic acid	100mg	HPFHpA	1546-95-8	7
DRE-C15986890	Perfluoroheptanoic acid	100mg	PFHpA	375-85-9	7
DRE-A15986890MW-100	Perfluoro-n-heptanoic acid 100 µg/mL in Methanol:Water	1ml	PFHpA	375-85-9	7
DRE-A15986891MW-50	Perfluoroheptanoic acid 13C4 (1,2,3,4-13C4) 50 µg/mL in Methanol:Water	1ml	13C4-PFHpA	2328024-55-9	7
DRE-A15986875MW-50	Perfluoroheptanoic acid 13C5 (1,2,3,4,5-13C5) 50 µg/mL in Methanol:Water	1ml	13C5-PFHpA	2328024-55-9	7
DRE-C15987147	2H,2H,3H,3H-Perfluorooctanoic acid	50mg	5:3 FTCA	914637-49-3	8
DRE-C15987145	2H,2H-Perfluorooctanoic acid	10mg	FHEA	53826-12-3	8
DRE-C15987162	2H-Perfluoro-2-octenoic acid	50mg	FHUEA	70887-88-6	8
DRE-A15987162AL-100	2H-Perfluoro-2-octenoic acid 100 µg/mL in Acetonitrile	1ml	FHUEA	70887-88-6	8
DRE-CA15986988	Perfluoro-6-methylheptanoic acid sodium	10mg	-	18017-22-6	8
DRE-A15986988MW-50	Perfluoro-6-methylheptanoic acid sodium 50 µg/mL in Methanol:Water	1ml	-	18017-22-6	8
DRE-CA15986984	Perfluoro-7-methyloctanoic acid sodium	10mg	-	-	8

PFAS reference materials

Product Listings

Per- and Polyfluorocarboxylic Acids (PFCA). Continued from previous page

Product Number	Product Description	Pack Size	Abbreviation	CAS Number	Carbon Chain Length
<i>DRE-A15986984MW-50</i>	Perfluoro-7-methyloctanoic acid sodium 50 µg/mL in Methanol:Water	1ml	-	-	8
DRE-C15987150	Perfluorooctanoic acid	100mg	PFOA	335-67-1	8
DRE-A15987150AL-100	Perfluorooctanoic acid 100 µg/mL in Acetonitrile	1ml	PFOA	335-67-1	8
DRE-A15987150MW-100	Perfluorooctanoic acid 100 µg/mL in Methanol:Water	1ml	PFOA	335-67-1	8
<i>DRE-A15987149MW-50</i>	Perfluorooctanoic acid 13C2 (1,2-13C2) sodium 50 µg/mL in Methanol:Water	1ml	13C2-PFOA	864071-08-11	8
<i>DRE-C15987152</i>	Perfluorooctanoic acid ammonium	100mg	PFOA	3825-26-1	8
DRE-A15987150MW-50	Perfluorooctanoic acid 50 µg/mL in Methanol:Water	1ml	PFOA	335-67-1	8
<i>DRE-A15987152ME-100</i>	Perfluorooctanoic acid ammonium 100 µg/mL in Methanol	1ml	PFOA	3825-26-1	8
DRE-A15987154MW-50	Perfluorooctanoic acid 13C4 (1,2,3,4-13C4) 50 µg/mL in Methanol:Water	1ml	13C4-PFOA	960315-48-4	8
DRE-A15987151MW-50	Perfluorooctanoic acid 13C8 50 µg/mL in Methanol:Water	1ml	13C8-PFOA	1350614-84-4	8
DRE-C15987115	Perfluorooctanesulfonamidoacetic acid	10mg	FOSAA	2806-24-8	8/2
DRE-C15987000	Perfluorononanoic acid	100mg	PFNA	375-95-1	9
DRE-A15987000AL-100	Perfluorononanoic acid 100 µg/mL in Acetonitrile	1ml	PFNA	375-95-1	9
DRE-C15987005	2H,2H,3H,3H-Perfluorononanoic acid	50mg	-	27854-30-4	9
DRE-A15987002MW-50	Perfluorononanoic acid 13C5 (1,2,3,4,5-13C5) 50 µg/mL in Methanol:Water	1ml	13C5-PFNA	960315-49-5	9
DRE-A15987001MW-50	Perfluorononanoic acid 13C9 50 µg/mL in Methanol:Water	1ml	13C9-PFNA	2283397-80-6	9
DRE-C15986604	1H,1H,2H,2H-Perfluorodecanoic acid	10mg	7:3 FTCA	812-70-4	10
DRE-C15986598	2H,2H-Perfluorodecanoic acid	10mg	FOEA	27854-31-5	10
DRE-C15986640	2H-Perfluoro-2-decenoic acid	10mg	FOUEA	70887-84-2	10
DRE-C15986608	Perfluoro(3,7-bis(trifluoromethyl))octanoic acid	100mg	PF-3,7-DMOA	172155-07-6	10
DRE-C15986600	Perfluorodecanoic acid	100mg	PFDA	335-76-2	10
DRE-A15986600MW-50	Perfluoro-n-decanoic acid 50 µg/mL in Methanol:Water	1ml	PFDA	335-76-2	10
DRE-A15986596MW-50	Perfluorodecanoic acid 13C2 (1,2-13C2) 50 µg/mL in Methanol:Water	1ml	13C2-PFDA	960315-50-8	10
DRE-A15986599MW-50	Perfluorodecanoic acid 13C6 (1,2,3,4,5,6-13C6) 50 µg/mL in Methanol:Water	1ml	13C6-PFDA	2328024-56-0	10
DRE-C13112600	11-H-Eicosafluoroundecanoic acid	100mg	HPFUnDA	1765-48-6	11
DRE-C15989010	2H,2H,3H,3H-Perfluoroundecanoic acid	50mg	H4PFUnDA	34598-33-9	11
DRE-A15989010ME-50	2H,2H,3H,3H-Perfluoroundecanoic acid 50 µg/mL in Methanol	1ml	H4PFUnDA	34598-33-9	11
DRE-C15989000	Perfluoroundecanoic acid	100mg	PFUnDA	2058-94-8	11
DRE-A15989000MW-50	Perfluoro-n-undecanoic acid 50 µg/mL in Methanol:Water	1ml	PFUnDA	2058-94-8	11
DRE-A15989002MW-50	Perfluoroundecanoic acid 13C2 (1,2-13C2) 50 µg/mL in Methanol:Water	1ml	13C2-PFUnDA	960315-51-9	11
DRE-A15989001MW-50	Perfluoroundecanoic acid 13C7 (1,2,3,4,5,6,7-13C7) 50 µg/mL in Methanol:Water	1ml	13C7-PFUnDA	-	11
DRE-C15986621	2H,2H-Perfluorododecanoic acid	25mg	FDEA	53826-13-4	12
DRE-A15986621MW-100	2H,2H-Perfluorododecanoic acid 100 µg/mL in Methanol:Water	1ml	FDEA	53826-13-4	12
DRE-C15986624	2H-Perfluoro-2-dodecenoic acid	10mg	FDUEA	70887-94-4	12
DRE-A15986624AL-100	2H-Perfluoro-2-dodecenoic acid 100 µg/mL in Acetonitrile	1ml	FDUEA	70887-94-4	12
DRE-C15986620	Perfluorododecanoic acid	50mg	PFDODA	307-55-1	12
DRE-A15986620MW-50	Perfluoro-n-dodecanoic acid 50 µg/mL in Methanol:Water	1ml	PFDODA	307-55-1	12
DRE-A15986628MW-50	Perfluorododecanoic acid 13C2 (1,2-13C2) 50 µg/mL in Methanol:Water	1ml	13C2-PFDODA	960315-52-0	12
DRE-C15988080	Perfluoro-2,5,8-trimethyl-3,6,9-trioxadodecanoic acid	50mg	-	65294-16-8	12
DRE-C15988000	Perfluorotridecanoic acid	100mg	PFTrDA	72629-94-8	13
DRE-A15988000MW-50	Perfluorotridecanoic acid 50 µg/mL in Methanol:Water	1ml	PFTrDA	72629-94-8	13
DRE-A15988003MW-50	Perfluorotridecanoic acid 13C2 (1,2-13C2) 50 µg/mL in Methanol:Water	1ml	13C2-PFTrDA	-	13
DRE-C15987400	Perfluorotetradecanoic acid	50mg	PFTeDA	376-06-7	14
DRE-A15987400MW-50	Perfluorotetradecanoic acid 50 µg/mL in Methanol:Water	1ml	PFTeDA	376-06-7	14
DRE-A15987410MW-50	Perfluorotetradecanoic acid 13C2 (1,2-13C2) 50 µg/mL in Methanol:Water	1ml	13C2-PFTeDA	2708218-82-8	14
DRE-C15986895	Perfluorohexadecanoic acid	50mg	PFHxDA	67905-19-5	16

PFAS reference materials

Product Listings

Per- and Polyfluorocarboxylic Acids (PFCA). Continued from previous page

Product Number	Product Description	Pack Size	Abbreviation	CAS Number	Carbon Chain Length
DRE-A15986895MW-50	Perfluorohexadecanoic acid 50 µg/mL in Methanol:Water	1ml	PFHxDA	67905-19-5	16
DRE-A15986899MW-50	Perfluorohexadecanoic acid 13C2 (1,2-13C2) 50 µg/mL in Methanol:Water	1ml	13C2-PFHxDA	-	16
DRE-C15987080	Perfluorooctadecanoic acid	50mg	PFODA	16517-11-6	18
DRE-A15987080MW-50	Perfluorooctadecanoic acid 50 µg/mL in Methanol:Water	1ml	PFODA	16517-11-6	18



Italic Product Number = Salt version
lgcstandards.com

Per- and Polyfluorosulfonic Acids

Per- and Polyfluorosulfonic Acids (PFSA) are a subgroup of Per- and Polyfluoroalkyl Substances (PFAS) characterised by a sulfonic acid functional group. Due to their exceptional stability and resistance to degradation, PFSA have been widely used in various industrial applications, including stain-resistant fabrics, firefighting foams and non-stick cookware. Their persistence in the environment and human body has raised significant health and environmental concerns, with studies linking PFSA to adverse effects such as cancer, hormone disruption and immune system impacts. Due to these risks, regulatory efforts are increasingly aimed at reducing their use and managing their environmental presence.

Product Number	Product Description	Pack Size	Abbreviation	CAS Number	Carbon Chain Length
DRE-C15986933	Perfluoromethane sulfinate sodium	100mg	Triflinate	2926-29-6	1
DRE-C15986935	Perfluoromethanesulfonic acid	250mg	TFMS	1493-13-6	1
DRE-A15986935MW-50	Perfluoromethanesulfonic acid 50 µg/mL in Methanol:Water	1ml	TFMS	1493-13-6	1
DRE-C17844760	Trifluoromethanesulfonic acid sodium	250mg	TFMS	2926-30-9	2
DRE-CA15986810	Perfluoroethanesulfonic acid	10mg	PFEtS	354-88-1	2
DRE-A15986810MW-50	Perfluoroethanesulfonic acid 50 µg/mL in Methanol:Water	1ml	PFEtS	354-88-1	2
DRE-CA15986815	Perfluoroethanesulfonic acid sodium	10mg	PFEtS	-	2
DRE-CA15987235	Perfluoropropanesulfonic acid sodium	10mg	-	359868-82-9	3
DRE-A15987234MW-50	Perfluoropropanesulfonic acid 50 µg/mL in Methanol:Water	1ml	PFPrS	423-41-6	3
DRE-C15986820	Perfluoro(2-ethoxyethane) sulfonic acid	100mg	PFEESA	113507-82-7	4
DRE-A15986820MW-50	Perfluoro(2-ethoxyethane) sulfonic acid 50 µg/mL in Methanol:Water	1ml	PFEESA	113507-82-7	4
DRE-CA15986515	Perfluorobutanesulfonic acid	100mg	PFBS	375-73-5	4
DRE-A15986515MW-100	Perfluorobutanesulfonic acid 100 µg/mL in Methanol:Water	1ml	PFBS	375-73-5	4
DRE-C15986517	Perfluorobutanesulfonic acid potassium	100mg	PFBSK	29420-49-3	4
DRE-A15986518MW-50	Perfluorobutanesulfonic acid 13C3 (2,3,4-13C3) sodium 50 µg/mL in Methanol:Water	1ml	13C3-PFBS	2708218-84-0	4
DRE-A15986518MW-50	Perfluorobutanesulfonic acid 13C3 (2,3,4-13C3) sodium 50 µg/mL in Methanol:Water	1ml	13C3-PFBS	2708218-84-0	4
DRE-A15986510MW-50	Perfluorobutanesulfonic acid 13C4 sodium 50 µg/mL in Methanol:Water	1ml	13C4-PFBS	-	4
DRE-CA15987190	Perfluoropentanesulfonic acid	25mg	PFPeS	2706-91-4	5
DRE-A15987190MW-100	Perfluoropentanesulfonic acid 100 µg/mL in Methanol:Water	1ml	PFPeS	2706-91-4	5
DRE-A15987205MW-50	Perfluoropentanesulfonic acid sodium 50 µg/mL in Methanol:Water	1ml	PFPeS	630402-22-1	5
DRE-A15987207MW-50	Perfluoropentanesulfonic acid 13C3 (1,2,3-13C3) sodium 50 µg/mL in Methanol:Water	1ml	13C3-PFPeS	-	5
DRE-A15986906MW-50	Perfluorohexanesulfonic acid 13C3 (1,2,3-13C3) sodium 50 µg/mL in Methanol:Water	1ml	13C3-PFHxS	2708218-86-2	6
DRE-A15986902MW-50	Perfluorohexanesulfonic acid 18O2 sodium 50 µg/mL in Methanol:Water	1ml	MPFHxS	-	6
DRE-C15986903	1H,1H,2H,2H-Perfluorohexanesulfonic acid	25mg	4:2 FTS	757124-72-4	6
DRE-A15986903MW-100	1H,1H,2H,2H-Perfluorohexanesulfonic acid 100 µg/mL in Methanol:Water	1ml	4:2 FTS	757124-72-4	6
DRE-A15986629MW-50	1H,1H,2H,2H-Perfluorohexanesulfonic acid 13C2 (1,2-13C2) sodium 50 µg/mL in Methanol:Water	1ml	13C2-4:2 FTS	2708218-88-4	6
DRE-C15986900	Perfluorohexanesulfonic acid	50mg	PFHxS	355-46-4	6
DRE-A15986900MW-50	Perfluorohexanesulfonic acid 50 µg/mL in Methanol:Water	1ml	PFHxS	355-46-4	6
DRE-A15986900AL-100	Perfluorohexanesulfonic acid 100 µg/mL in Acetonitrile	1ml	PFHxS	355-46-4	6
DRE-CA15986905	Perfluorohexanesulfonic acid sodium	10mg	PFHxS	82382-12-5	6
DRE-A15986906MW-50	Perfluorohexanesulfonic acid 13C3 (1,2,3-13C3) sodium 50 µg/mL in Methanol:Water	1ml	13C3-PFHxS	2708218-86-2	6
DRE-CA15986946	Perfluorohexanesulfonic acid (linear+branched)	25mg	-	-	6
DRE-A15986946MW-100	Perfluorohexanesulfonic acid (linear+branched) 100 µg/mL in Methanol:Water	1ml	-	-	6
DRE-CA15986614	7H-Perfluoro-3,6-dioxa-4-methyl-octane-1-sulfonic acid	10mg	Hydro-PS Acid	749836-20-2	7
DRE-CA15986615	Perfluoro-3,6-dioxa-4-methyl-7-octene-1-sulfonic acid	10mg	-	29311-67-9	7
DRE-C15986880	Perfluoroheptanesulfonic acid	50mg	PFHpS	375-92-8	7

PFAS reference materials

Product Listings

Per- and Polyfluorosulfonic Acids. Continued from previous page

Product Number	Product Description	Pack Size	Abbreviation	CAS Number	Carbon Chain Length
<i>DRE-A15986880AL-100</i>	Perfluoroheptanesulfonic acid 100 µg/mL in Acetonitrile	1ml	PFHpS	375-92-8	7
<i>DRE-A15986883MW-50</i>	Perfluoroheptanesulfonic acid 13C3 (1,2,3-13C3) sodium 50 µg/mL in Methanol:Water	1ml	13C3-PFHpS	-	7
<i>DRE-A15986845MW-50</i>	Perfluoro-4-ethylcyclohexylsulfonic acid potassium 50 µg/mL in Methanol:Water	1ml	PFECHS	335-24-0	8
<i>DRE-C15987125</i>	1H,1H,2H,2H-Perfluorooctane sulfonic acid	10mg	6:2 FTS	27619-97-2	8
<i>DRE-A15987125ME-100</i>	1H,1H,2H,2H-Perfluorooctane sulfonic acid 100 µg/mL in Methanol	1ml	6:2 FTS	27619-97-2	8
<i>DRE-A15987126MW-50</i>	1H,1H,2H,2H-Perfluorooctane sulfonic acid sodium 50 µg/mL in Methanol:Water	1ml	6:2 FTS	27619-94-9	8
<i>DRE-A15987127MW-50</i>	1H,1H,2H,2H-Perfluorooctanesulfonic acid 13C2 (1,2-13C2) sodium 50 µg/mL in Methanol:Water	1ml	13C2-6:2 FTS	2708218-89-5	8
<i>DRE-A11459000MW-50</i>	9-Chlorohexadecafluoro-3-oxanonane-1-sulfonic acid potassium 50 µg/mL in Methanol:Water	1ml	9Cl-PF3ONS	73606-19-6	8
<i>DRE-C11041290</i>	Capstone product A	10mg	-	80475-32-7	8
<i>DRE-C11041300</i>	Capstone product B	25mg	-	34455-29-3	8
<i>DRE-A11041320MW-50</i>	Capstone product B 13C3 (dimethyl 13C2, glycine-beta-13C1) 50 µg/mL in Methanol:Water	1ml	-	34455-29-3	8
<i>DRE-A11041300MW-100</i>	Capstone product B 100 µg/mL in Methanol:Water	1ml	-	34455-29-3	8
<i>DRE-CA15987120</i>	Perfluorooctane sulfonic acid	25mg	PFOS	1763-23-1	8
<i>DRE-A15987120MW-100</i>	Perfluorooctane sulfonic acid 100 µg/mL in Methanol:Water	1ml	PFOS	1763-23-1	8
<i>DRE-C15987122</i>	Perfluorooctane sulfonic acid potassium	100mg	PFOS	2795-39-3	8
<i>DRE-A15987122MW-100</i>	Perfluorooctane sulfonic acid potassium 100 µg/mL in Methanol:Water	1ml	PFOS	2795-39-3	8
<i>DRE-C15987130</i>	Perfluorooctane sulfonyl fluoride	100mg	PFOS	307-35-7	8
<i>DRE-CA15987131</i>	Perfluorooctanesulfonic acid (linear+branched)	25mg	PFOS	-	8
<i>DRE-A15987131MW-100</i>	Perfluorooctanesulfonic acid (linear+branched) 100 µg/mL in Methanol:Water	1ml	PFOS	-	8
<i>DRE-A15987123MW-50</i>	Perfluorooctanesulfonic acid 13C4 (1,2,3,4-13C4) sodium 50 µg/mL in Methanol:Water	1ml	13C4-PFOS	960315-53-1	8
<i>DRE-A15987124MW-50</i>	Perfluorooctanesulfonic acid 13C8 sodium 50 µg/mL in Methanol:Water	1ml	13C8-PFOS	2522762-16-7	8
<i>DRE-CA15987022</i>	Perfluorononanesulfonic acid sodium	5mg	PFNS	98789-57-2	9
<i>DRE-A15987022MW-50</i>	Perfluorononanesulfonic acid sodium 50 µg/mL in Methanol:Water	1ml	PFNS	98789-57-2	9
<i>DRE-A15987023MW-50</i>	Perfluorononanesulfonic acid 13C3(1,2,3-13C3) sodium 50 µg/mL in Methanol:Water	1ml	13C3-PFNS	-	9
<i>DRE-A15986582MW-50</i>	Perfluorodecanesulfonic acid 13C3(1,2,3-13C3) sodium 50 µg/mL in Methanol:Water	1ml	13C3-PFDS	-	10
<i>DRE-CA11407100</i>	11-Chloroeicosafluoro-3-oxaundecane-1-sulfonic acid potassium	5mg	11Cl- PF3OUdS	83329-89-9	10
<i>DRE-A11407100MW-50</i>	11-Chloroeicosafluoro-3-oxaundecane-1-sulfonic acid potassium 50 µg/mL in Methanol:Water	1ml	11Cl- PF3OUdS	83329-89-9	10
<i>DRE-C15986585</i>	1H,1H,2H,2H-Perfluorodecanesulfonic acid	25mg	8:2 FTS	39108-34-4	10
<i>DRE-A15986585MW-50</i>	1H,1H,2H,2H-Perfluorodecanesulfonic acid 50 µg/mL in Methanol:Water	1ml	8:2 FTS	39108-34-4	10
<i>DRE-A15986587MW-50</i>	1H,1H,2H,2H-Perfluorodecanesulfonic acid 13C2 (1,2-13C2) sodium 50 µg/mL in Methanol:Water	1ml	13C2-8:2 FTS	2708218-90-8	10
<i>DRE-A15986586MW-50</i>	1H,1H,2H,2H-Perfluorodecanesulfonic acid sodium 50 µg/mL in Methanol:Water	1ml	8:2 FTS	27619-96-1	10
<i>DRE-CA15986580</i>	Perfluorodecanesulfonic acid	5mg	PFDS	335-77-3	10
<i>DRE-A15986580MW-50</i>	Perfluorodecanesulfonic acid 50 µg/mL in Methanol:Water	1ml	PFDS	335-77-3	10
<i>DRE-A15986580AL-100</i>	Perfluorodecanesulfonic acid 100 µg/mL in Acetonitrile	1ml	PFDS	335-77-3	10
<i>DRE-CA15986581</i>	Perfluorodecanesulfonic acid sodium	5mg	PFDS	2806-15-7	10
<i>DRE-A15986581AL-100</i>	Perfluorodecanesulfonic acid sodium 100 µg/mL in Acetonitrile	1ml	PFDS	2806-15-7	10
<i>DRE-A15986581MW-50</i>	Perfluorodecanesulfonic acid sodium 50 µg/mL in Methanol:Water	1ml	PFDS	2806-15-7	10
<i>DRE-A15986582MW-50</i>	Perfluorodecanesulfonic acid 13C3(1,2,3-13C3) sodium 50 µg/mL in Methanol:Water	1ml	13C3-PFDS	-	10
<i>DRE-A15988996MW-50</i>	Perfluoroundecanesulfonic acid 13C3(1,2,3-13C3) sodium 50 µg/mL in Methanol:Water	1ml	13C3-PFUnDS	-	11
<i>DRE-A15988995MW-50</i>	Perfluoroundecanesulfonic acid sodium 50 µg/mL in Methanol:Water	1ml	PFUnDS	441296-91-9	11
<i>DRE-C15986622</i>	1H,1H,2H,2H-Perfluorododecane sulfonic acid	25mg	10:2 FTS	120226-60-0	12

PFAS reference materials

Product Listings

Per- and Polyfluorosulfonic Acids. Continued from previous page

Product Number	Product Description	Pack Size	Abbreviation	CAS Number	Carbon Chain Length
DRE-A15986622MW-100	1H,1H,2H,2H-Perfluorododecane sulfonic acid 100 µg/mL in Methanol:Water	1ml	10:2 FTS	120226-60-0	12
DRE-CA15986632	1H,1H,2H,2H-Perfluorododecanesulfonic acid sodium	10mg	10:2 FTS	108026-35-3	12
DRE-A15986632AC-50	1H,1H,2H,2H-Perfluorododecane sulfonic acid sodium 50 µg/mL in Acetone	1ml	10:2 FTS	108026-35-3	12
DRE-A15986627MW-50	Perfluorododecanesulfonic acid sodium 50 µg/mL in Methanol:Water	1ml	PFD _o DS	1260224-54-1	12
DRE-A15986645MW-50	Perfluorododecanesulfonic acid 13C3 (1,2,3-13C3) sodium 50 µg/mL in Methanol:Water	1ml	13C3-PFD _o DS	-	12
DRE-A15986633MW-50	1H,1H,2H,2H-Perfluorododecanesulfonic acid 13C2 D4 (1,2-13C2 D4) sodium 50 µg/mL in Methanol:Water	1ml	13C2-D4-10:2 FTS	-	12
DRE-A15987995MW-50	Perfluorotridecanesulfonic acid sodium 50 µg/mL in Methanol:Water	1ml	PFT _r DS	791563-89-8	13
DRE-A15987996MW-50	Perfluorotridecanesulfonic acid 13C3 (1,2,3-13C3) sodium 50 µg/mL in Methanol:Water	1ml	13C3-PFT _r DS	-	13



Fluorotelomer Alcohols

Volatile precursors to Perfluorinated Carboxylic Acids (PFCAs), Fluorotelomer Alcohols (FTOHs) are fluorotelomers with an alcohol functional group. FTOHs are environmental contaminants characterised by high volatility, prone to long-range atmospheric transport. FTOHs have been identified in consumer products, including stain repellents and windscreen wiper fluid.

Product Number	Product Description	Pack Size	Abbreviation	CAS Number	Carbon Chain Length
DRE-CA15987237	1H,1H-Perfluoropropanol	100mg	-	422-05-9	3
DRE-CA14194190	1,1,1,1,3,3,3-Hexafluoro-2-propanol	1g	-	920-66-1	3
DRE-C15986540	1H,1H-Perfluorobutanol	100mg	-	375-01-9	4
DRE-C15986913	1H,1H-Perfluorohexanol	50mg	-	423-46-1	6
DRE-C15986915	1H,1H,2H,2H-Perfluoro-1-hexanol	100mg	4:2 FTOH	2043-47-2	6
DRE-C15986886	1H,1H,7H-Perfluoroheptanol	1g	-	335-99-9	7
DRE-CA15986887	1H,1H,2H,2H,3H,3H-Perfluoroheptanol	100mg	-	83310-97-8	7
DRE-CA15987225	1H,1H,1H,2H-Perfluoro-2-heptanol	25mg	5:2 sFTOH	914637-05-1	7
DRE-C15987160	1H,1H,2H,2H-Perfluoro-1-octanol	100mg	6:2 FTOH	647-42-7	8
DRE-C15987161	1H,1H,8H-Perfluorooctanol	250mg	-	10331-08-5	8
DRE-C15986893	1H,1H,1H,2H-Perfluoro-2-nonanol	10mg	7:2 sFTOH	24015-83-6	9
DRE-CA15986922	1H,1H,2H,2H,3H,3H-Perfluorononanol	50mg	-	80806-68-4	9
DRE-C15986601	1H,1H,2H,2H-Perfluoro-1-decanol	100mg	8:2 FTOH	678-39-7	10
DRE-CA15989090	1H,1H-Perfluoroundecanol	50mg	-	307-46-0	11
DRE-CA15989100	1H,1H,11H-Perfluoro-1-undecanol	50mg	9:2 sFTOH	307-70-0	11
DRE-C16986625	1H,1H,2H,2H-Perfluoro-1-dodecanol	100mg	10:2 FTOH	865-86-1	12
DRE-CA15987415	1H,1H,2H,2H-Perfluorotetradecanol	50mg	-	39239-77-5	14



Perfluoroalkyl Sulfonamides

Perfluoroalkyl Sulfonamides (PFASAs) are a subclass of Per- and Polyfluoroalkyl Substances (PFAS) characterised by their strong carbon-fluorine bonds and the presence of a sulfonamide group. These compounds are known for their exceptional chemical stability, resistance to heat, and water and oil repellency. Due to these properties, PFASAs have been widely used in various industrial applications, including as surfactants, lubricants, and in the production of non-stick coatings and textiles. However, their environmental persistence and potential health impacts have raised significant concerns, leading to increased regulatory scrutiny and efforts to find safer alternatives.

Product Number	Product Description	Pack Size	Abbreviation	CAS Number	Carbon Chain Length
DRE-C15128000	N-(2-Hydroxyethyl)-N-methyl-perfluorobutane sulfonamide	25mg	-	34454-97-2	4
DRE-C10653522	N,N-Bis(2-hydroxyethyl)perfluorobutane sulfonamide	100mg	-	34455-00-0	4
DRE-C12723258	N-[3-(Dimethylamino)propyl]perfluoro1- butanesulfonamide	10mg	-	68555-77-1	4
DRE-C15115450	N-Methylperfluorobutanesulfonamide	50mg	MeFBSA	68298-12-4	4
DRE-C15986519	Perfluorobutanesulfonamide(N-methyl)acetate	10mg	MeFBSAA	159381-10-9	4
DRE-C15986512	Perfluorobutane sulfonamide	50mg	PFBSA	30334-69-1	4
DRE-C12723276	N-[3-(Dimethylamino)propyl]perfluoro1- hexanesulfonamide	10mg	-	50598-28-2	6
DRE-C15115470	N-Methylperfluorohexane sulfonamide	50mg	N-Me-FHxSA	68259-15-4	6
DRE-C15986898	Perfluorohexane sulfonamide	50mg	FHxSA	41997-13-1	6
DRE-A15986898MW-50	Perfluorohexanesulfonamide 50 µg/mL in Methanol:Water	1ml	PFHxSA	41997-13-1	6
DRE-A13349620MW-50	2-(N-Ethylperfluorooctanesulfonamido)acetic acid 13C4 50 µg/mL in Methanol:Water	1ml	13C4-N-Et-FOSAA	-	8
DRE-C15987115	2-(N-Ethylperfluorooctanesulfonamido)acetic acid D5 (N-ethyl D5) 50 µg/mL in Methanol:Water	1ml	D5-N-EtFOSAA	1265205-97-7	8
DRE-CA15130000	2-(N-Methylperfluorooctanesulfonamido)acetic acid	10mg	N-MeFOSAA	2355-31-9	8
DRE-A15130000MW-50	2-(N-Methylperfluorooctanesulfonamido)acetic acid 50 µg/ mL in Methanol:Water	1ml	N-MeFOSAA	2355-31-9	8
DRE-A15130000AL-100	2-(N-Methylperfluorooctanesulfonamido)acetic acid 100 µg/mL in Acetonitrile	1ml	N-MeFOSAA	2355-31-9	8
DRE-CA13349600	2-(N-Ethylperfluorooctanesulfonamido)acetic acid	25mg	N-EtFOSAA	2991-50-6	8
DRE-A13349600MW-50	2-(N-Ethylperfluorooctanesulfonamido)acetic acid 50 µg/ mL in Methanol:Water	1ml	N-EtFOSAA	2991-50-6	8
DRE-A13349600AL-100	2-(N-Ethylperfluorooctanesulfonamido)acetic acid 100 µg/ mL in Acetonitrile	1ml	N-EtFOSAA	2991-50-6	8
DRE-A15130020MW-50	2-(N-Methylperfluorooctanesulfonamido)acetic acid 13C3 50 µg/mL in Methanol:Water	1ml	13C3-N-Me-FOSAA	-	8
DRE-A15130010MW-50	2-(N-Methylperfluorooctanesulfonamido)acetic acid D3 (N-methyl D3) 50 µg/mL in Methanol:Water	1ml	D3-N-MeFOSAA	1400690-70-1	8
DRE-C14231570	N-(2-Hydroxyethyl)-N-methylperfluorooctanesulfonamide	25mg	N-MeFOSE	24448-09-7	8
DRE-A14231570ME-100	N-(2-Hydroxyethyl)-N-methylperfluorooctanesulfonamide 100 µg/mL in Methanol	1ml	N-MeFOSE	24448-09-7	8
DRE-A14231571MW-50	N-(2-Hydroxyethyl)-N-methylperfluorooctanesulfonamide D7 50 µg/mL in Methanol:Water	1ml	D7-N-MeFOSE	1265205-95-5	8
DRE-C12723278	N-[3-(Dimethylamino)propyl]perfluoro1- octanesulfonamide	10mg	-	13417-01-1	8
DRE-C13342360	N-Ethyl-N-(2-hydroxyethyl)perfluorooctanesulfonamide	50mg	N-EtFOSE	1691-99-2	8
DRE-A13342360ME-100	N-Ethyl-N-(2-hydroxyethyl)perfluorooctanesulfonamide 100 µg/mL in Methanol	1ml	N-EtFOSE	1691-99-2	8
DRE-A13342361MW-50	N-Ethyl-N-(2-hydroxyethyl)perfluorooctanesulfonamide D9 50 µg/mL in Methanol:Water	1ml	D9-N-EtFOSE	1265205-96-6	8
DRE-C15115500	N-Methylperfluorooctanesulfonamide	50mg	N-MeFOSA	31506-32-8	8
DRE-A15115500MW-100	N-Methylperfluorooctanesulfonamide 100 µg/mL in Methanol:Water	1ml	N-MeFOSA	31506-32-8	8
DRE-A15115501MW-50	N-Methylperfluorooctanesulfonamide D3 (methyl D3) 50 µg/mL in Methanol:Water	1ml	D3-N-MeFOSA	936109-37-4	8
DRE-C15987110	Perfluorooctane sulfonamide	100mg	PFOSA	754-91-6	8
DRE-A13349601MW-50	Perfluorooctane sulfonamide 18O2 50 µg/mL in Methanol:Water	1ml	18O2-PFOSA	1340480-46-7	8
DRE-C17004000	Sulfuramide	100mg	EtFOSA	4151-50-2	8
DRE-A17004001MW-50	Sulfuramid D5 (ethyl D5) 50 µg/mL in Methanol:Water	1ml	D5-N-EtFOSA	936109-40-9	8
DRE-A15987112MW-50	Perfluorooctanesulfonamidoacetic acid	10mg	-	-	8/2
DRE-C12723280	N-[3-(Dimethylamino)propyl]-1H,1H,2H,2H-perfluoro-1-octanesulfonamide	10mg	-	34455-22-6	13

Italic Product Number = Salt version

Polyfluoroalkyl Acrylates

Polyfluoroalkyl Acrylates are a class of Per- and Polyfluoroalkyl Substances (PFAS) used primarily as polymerisation agents in the production of water- and stain-resistant coatings, textiles and paper products. Despite their beneficial applications, polyfluoroalkyl acrylates raise environmental and health concerns due to their persistence and potential to degrade into more stable and bioaccumulative PFAS, such as Perfluorooctanoic Acid (PFOA). Consequently, there is growing regulatory scrutiny and research focused on safer alternatives and the environmental impact of these substances.

Product Number	Product Description	Pack Size	Abbreviation	CAS Number	Carbon Chain Length
DRE-CA15987227	1H,1H,5H-Perfluoropentyl methacrylate	100mg	-	355-93-1	5
DRE-CA15987169	1H,1H,2H,2H-Perfluorooctyl methacrylate	250mg	6:2 FTMAC	2144-53-8	8
DRE-C15987170	1H,1H,2H,2H-Perfluorooctyl acrylate	100mg	-	17527-29-6	8
DRE-CA15986595	1H,1H,2H,2H-Perfluorodecyl methacrylate	100mg	8:2 FTMAC	1996-88-9	10
DRE-C15986602	1H,1H,2H,2H-Perfluorodecyl acrylate	100mg	8:2FTAc	27905-45-9	10
DRE-CA15986658	1H,1H,2H,2H-Perfluorododecyl methacrylate	50mg	-	2144-54-9	12
DRE-C15986630	1H,1H,2H,2H-Perfluorododecyl acrylate	50mg	10:2FTAc	17741-60-5	12

Polyfluorinated Phosphate Esters

Polyfluorinated Phosphate Esters are a class of PFAS compounds that have useful properties, including high stability, low surface tension and resistance to heat and chemicals. Therefore, they are often used as surfactants, lubricants and in fire-retardant applications. However, despite their useful applications, they have raised concerns in the scientific community due to their persistence in the environment and potential health risks.

Product Number	Product Description	Pack Size	Abbreviation	CAS Number	Carbon Chain Length
DRE-A10655199MW-50	Bis[2-(perfluorohexyl)ethyl] phosphate 13C4 (1,1',2,2'-13C4) sodium 50 µg/mL in Methanol:Water	1ml	13C4-6:2 Di-PAPS	2179232-14-3	8/8
DRE-A10655210AL-100	Bis(2-perfluorooctylsulfonyl-N-ethylaminoethyl) phosphate 100 µg/mL in Acetonitrile	1ml	diSAmPAP	2965-52-8	8/8
DRE-C10655185	Bis[2-(perfluorohexyl)ethyl] phosphate	10mg	6:2 Di-PAPS	57677-95-9	8/8
DRE-C15312100	Mono[2-(perfluorooctyl)ethyl] phosphate	5mg	8:2 PAPS	57678-03-2	10
DRE-CA17894475	Tris(1H,1H,2H,2H-perfluorooctyl)phosphate	10mg	6:2 Tri-PAPS	165325-62-2	8/8/8
DRE-C15987168	(1H,1H,2H,2H-Perfluorooctyl-1H,1H,2H,2H-perfluorodecyl)phosphate	5mg	6:2/8:2 Di-PAPS	943913-15-3	8/10
DRE-C10655195	Bis[2-(perfluorooctyl)ethyl] phosphate	10mg	8:2 Di-PAPS	678-41-1	10/10
DRE-C17894473	Tris(1H,1H,2H,2H-perfluorodecyl)phosphate	10mg	-	149790-22-7	10/10/10
DRE-C15986664	1H,1H,2H,2H-Perfluorododecyl phosphate bis(cyclohexylammonium)	10mg	-	-	12

Per- and Polyfluoroalkyl Ether Carboxylic Acids

Per- and Polyfluoroalkyl Ether Carboxylic Acids are a type of perfluorinated compound used in various industrial applications such as surfactants, lubricants and coatings. They are particularly sought after as they are exceptionally stable, water and oil repellent and resistant to heat and chemicals. However, like with other PFAS compounds, they are under scrutiny for their persistence in the environment and potential health effects, making them a crucial group of compounds to test for.

Product Number	Product Description	Pack Size	Abbreviation	CAS Number	Carbon Chain Length
DRE-A15986940MW-50	Perfluoro-2-methoxyacetic acid 50 µg/mL in Methanol:Water	1ml	PFMCAA	674-13-5	3
DRE-CA15986970	Perfluoro-3-methoxypropanoic acid	50mg	PFMOPrA	377-73-1	4
DRE-A15986970MW-50	Perfluoro-3-methoxypropanoic acid 50 µg/mL in Methanol:Water	1ml	PFMOPrA	377-73-1	4
DRE-C15986612	Perfluoro-3,6-dioxaheptanoic acid	100mg	NFDHA	151772-58-6	5
DRE-C15986950	Perfluoro-4-methoxybutanoic acid	25mg	PFMOBA	863090-89-5	5
DRE-A15986950MW-50	Perfluoro-4-methoxybutanoic acid 50 µg/mL in Methanol:Water	1ml	PFMOBA	863090-89-5	5
DRE-CA15987250	Perfluoro-2-propoxypropanoic acid	50mg	PFPrOPrA	13252-13-6	6
DRE-A15987250MW-100	Perfluoro-2-propoxypropanoic acid 100 µg/mL in Methanol:Water	1ml	PFPrOPrA	13252-13-6	6
DRE-C15986625	Perfluoro-3,6-dioxaoctanoic acid	10mg	-	80153-82-8	6
DRE-C15986618	3H-Perfluoro-4,8-dioxanonanoic acid	10mg	DONA	919005-14-4	7
DRE-C15986605	Perfluoro-2,5-dimethyl-3,6-dioxanonanoic acid	100mg	HFPO-TA	13252-14-7	9



Other PFAS

Product Number	Product Description	Pack Size	Abbreviation	CAS Number	Carbon Chain Length
DRE-C10655193	Bis(perfluoromethanesulfonyl)imide lithium	100mg	-	90076-65-6	1/1
DRE-CA15115455	Perfluorobutanoic acid-methyl ester	100mg	-	356-24-1	4
DRE-C15987500	Perfluorotributylamine	100mg	PFTBA	311-89-7	4/4/4
DRE-A15987500MB-1000	Perfluorotributylamine (PFTBA) 1000 µg/mL in Methyl-tert-butyl ether	1ml	PFTBA	311-89-7	4/4/4
DRE-CA15986555	1H,1H,2H-Perfluorocyclopentane	500mg	-	15290-77-4	5
DRE-A17846310AL-100	1,1,2-Trifluoro-2-(heptafluoropropoxy)ethene 100 µg/mL in Acetonitrile	1ml	-	1623-05-8	5
DRE-CA15115550	Perfluoropentanoic acid-methyl ester	50mg	-	13038-26-1	5
DRE-C15986990	Perfluoro-2-methyl-3-pentanone	100mg	-	756-13-8	6
DRE-CA15115480	Perfluorohexanoic acid-methyl ester	50mg	-	424-18-0	6
DRE-C10655190	Bis(perfluorohexyl)phosphinic acid	25mg	6:6 PFPi	40143-77-9	6/6
DRE-C15986553	(Perfluorocyclohexyl)methanol	100mg	-	28788-68-3	7
DRE-CA15986865	Perfluoroheptane	100mg	-	335-57-9	7
DRE-CA15986527	Perfluorobutanoic acid-ethyl ester	25mg	-	356-27-4	8
DRE-C15115250	Perfluorooctanoic acid methyl ester	100mg	-	376-27-2	8
DRE-C15987155	Perfluorooctanoic acid-ethyl ester	100mg	-	3108-24-5	8
DRE-C15987172	(1H,1H,2H,2H-Perfluorooctyl)triethoxysilane	100mg	-	51851-37-7	8
DRE-CA15987100	Perfluorooctane	100mg	-	307-34-6	8
DRE-CA15987175	(1H,1H,2H,2H-Perfluorooctyl)trimethoxysilane	50mg	-	85857-16-5	8
DRE-C15987101	1H,1H,2H,2H-Perfluorooctanephosphonic acid	50mg	PFPA	252237-40-4	8
DRE-CA15987139	1H,1H,1H,2H-Perfluorooctanethiol	50mg	-	34451-26-8	8
DRE-C15987178	Perfluorooctanephosphonic acid	10mg	-	40143-78-0	8
DRE-CA15986999	Perfluorononane	250mg	-	375-96-2	9
DRE-C15986560	Perfluorodecanephosphonic acid	10mg	PFDPa	52299-26-0	10
DRE-C15986603	1H,1H,2H,2H-Perfluorodecyl acetate	100mg	8:2 FTOAc	37858-04-1	10
DRE-CA15986606	(1H,1H,2H,2H-Perfluorodecyl)triethoxysilane	100mg	-	101947-16-4	10
DRE-C15986559	Perfluorodecane	100mg	-	307-45-9	10
DRE-C15988980	1H-Perfluoroundecane	100mg	-	66256-53-9	11
DRE-C15986638	Perfluorododecane	100mg	-	307-59-5	12
DRE-C15987980	Perfluorotridecane	500mg	-	376-03-4	13
DRE-C15987390	Perfluorotetradecane	100mg	-	307-62-0	14
DRE-CA15986929	1H,1H,2H,2H-Perfluoro-1-iodotetradecane	50mg	-	30046-31-2	14
DRE-C15987184	Perfluoropentadecane	100mg	-	133010	15
DRE-C15986884	Perfluorohexadecane	100mg	-	355-49-7	16

Matrix Materials

Due to their persistent nature, PFAS compounds bioaccumulate in various matrices, including water, soil, sediment, air and biological matrices. The interaction between PFAS compounds and these matrices can complicate analyses. Therefore, the use of matrix materials plays a crucial role in understanding these effects. LGC offers a comprehensive range of matrix materials for PFAS testing, thanks to our partnerships with leading reference material producers.

Part Number	Part Description
IRMM-427	Pike Perch - PFASs in fish tissue
IRMM-428	PFASs in Water
NIST-2781	Domestic Sludge
NIST-2585	Organic Contaminants in House Dust
NIST-1957	Organic Contaminants in Non-Fortified Human Serum (Freeze-Dried)
NIST-1958	Organic Contaminants in Fortified Human Serum (Freeze-Dried)
NIST-1950	Metabolites in Human Plasma
NIST-8690	PFAS in Aqueous Film-Forming Foams (AFFF) Formulation III
NIST-8691	PFAS in Aqueous Film-Forming Foams (AFFF) Formulation III
NIST-8692	PFAS in Aqueous Film-Forming Foams (AFFF) Formulation III
NIST-8693	PFAS in Aqueous Film-Forming Foams (AFFF) Formulation III
NIST-2586	Trace Elements in Soil Containing Lead From Paint
BAM-B003	Per- and polyfluorinated alkyl substances (PFAS) in textiles



Italic Product Number = Salt version

[lgcstandards.com](https://www.lgcstandards.com)

PFAS Proficiency Testing

AXIO PROFICIENCY TESTING

At LGC, reference materials are just the start of the quality assurance process. Dr. Ehrenstorfer's sister brand, LGC AXIO Proficiency Testing (PT), offers a wide range of PFAS Proficiency Testing (PT) to ensure that your testing remains accurate, whether you work with water, soil, food or beverages. Furthermore, our best-in-class data reporting platform, PORTAL, provides your laboratory with the tools to monitor and effectively analyse your performance over time, allowing you to identify potential performance issues well in advance. All of our PFAS PT samples are designed in line with industry regulations, such as the EU Drinking Water Directive, plus the UK Drinking Water Inspectorate's guidelines.

Part Number	Part Name	Analytes	Matrix	Format	Dispatch date*
PT-AQ-67	PFAS in potable water (EU Drinking Water Directive)	20 PFAS	Potable water	2L + spike	11-May-25
PT-AQ-70	PFAS in potable water (DWI Guidance)	48 PFAS	Potable water	2L + spike	08-Jun-26 30-Nov-26
PT-AQ-73	Fluorotelomer Alcohols (FTOHs) in potable water	6:2 Fluorotelomer Alcohol; 8:2 Fluorotelomer Alcohol	Potable water	1L + spike	08-Jun-26
PT-CN-32	PFAS in soil	Perfluorononanoic acid (PFNA); Perfluorooctanoic acid (PFOA); Perfluorohexane sulfonic acid (PFHxS); Perfluorooctane sulfonic acid (PFOS); Sum of 4 PFAS	Soil	30g	07-Apr-26
PT-MT-765	PFAS in meat or fish	At least 10 PFAS from a list (see scheme description)	Fish or meat	30g	16-Mar-26
PT-FC-910	PFAS in egg powder	At least 10 PFAS from a list (see scheme description)	Egg powder	30g	20-Oct-25 19-Oct-26
PT-BA-11	PFAS in beer	Perfluorobutanoic acid (PFBA); Perfluorooctanoic acid (PFOA); Perfluorononanoic acid (PFNA); Perfluorodecanoic acid (PFDA); Perfluorohexanoic acid (PFHxA); Perfluoroheptanoic acid (PFHpA); Perfluorohexane sulfonic acid (PFHxS); Perfluoroheptane sulfonic acid (PFHpS); Perfluorooctane sulfonic acid (PFOS); GenX (HFPODA); NaDONA; 11CI-PF3ONS; 9CI-PF3ONS	Beer	50ml	16-Mar-26

**PT samples are dispatched to all customers at the same time. To get involved in the next PT round, please order your sample at least a week before the dispatch date. Results must generally be submitted within five weeks following the dispatch date.*




To explore AXIO's full range of food, beverage, environment schemes and more, scan or [click here](#)

For full listing of offices and distributors,
please contact us: lgcstandards.com



Contact us



Sales offices

LGC

Science
for a Safer World

lgcstandards.com