

Your **Power** for Health


greiner bio-one



VACUETTE® VISIO PLUS Needle

Bringing Vision to
Blood Collection

VACUETTE® VISIO PLUS Needle

The **VACUETTE®** VISIO PLUS Needle with optical venipuncture control can make an important contribution to simplifying blood collection routine for phlebotomists.

The needle has a translucent plastic hub, which functions as a view window.

If venipuncture is successful, the view window is coloured red due to the blood flow.

The **VACUETTE®** VISIO PLUS Needle guarantees a safe blood collection procedure for the patient, whilst at the same time ensuring best possible comfort. Optical venipuncture control via the view window makes it much easier to find the vein lumen.

Re-positioning in the vein which may cause patient discomfort, can be avoided with the **VACUETTE®** VISIO PLUS venipuncture aid.

Perfected grinding technology allows for sharp needle tip and optimal tissue penetration.

To be able to differentiate from the **VACUETTE®** standard needles, the **VACUETTE®** VISIO PLUS blood collection needles have transparent needle containers.

VACUETTE® VISIO PLUS – the Greiner Bio-One innovation.

Item n°	Description	Size	Packaging Box / Carton
450040	VACUETTE® VISIO PLUS Needle green, sterile, latex-free	21G x 1 ½" / 0.8 x 38 mm	100 pcs./2000 pcs.
450041	VACUETTE® VISIO PLUS Needle black, sterile, latex-free	22G x 1 ½" / 0.7 x 38 mm	100 pcs./2000 pcs.
450042	VACUETTE® VISIO PLUS Needle green, sterile, latex-free	21G x 1" / 0.8 x 25 mm	100 pcs./2000 pcs.
450043	VACUETTE® VISIO PLUS Needle black, sterile, latex-free	22G x 1" / 0.7 x 25 mm	100 pcs./2000 pcs.
450239	VACUETTE® QUICKSHIELD Complete PLUS single-packed, sterile Safety Tube Holder & VACUETTE® VISIO PLUS Needle Combination Product	21G x 1 ½" / 0.8 x 38 mm	30 pcs./300 pcs.
450228	VACUETTE® QUICKSHIELD Complete PLUS single-packed, sterile Safety Tube Holder & VACUETTE® VISIO PLUS Needle Combination Product	22G x 1 ½" / 0.7 x 38 mm	30 pcs./300 pcs.
450235	VACUETTE® QUICKSHIELD Complete PLUS single-packed, sterile Safety Tube Holder & VACUETTE® VISIO PLUS Needle Combination Product	21G x 1" / 0.8 x 25 mm	30 pcs./300 pcs.
450226	VACUETTE® QUICKSHIELD Complete PLUS single-packed, sterile Safety Tube Holder & VACUETTE® VISIO PLUS Needle Combination Product	22G x 1" / 0.7 x 25 mm	30 pcs./300 pcs.

Perfected grinding technology allows for sharp needle tip and optimal tissue penetration.

VACUETTE® QUICKSHIELD Complete PLUS

Simply Practical.

The **VACUETTE®** QUICKSHIELD Safety Tube Holder is especially suitable for routine blood collection. This product can provide the user with the simplest handling and reliable infection protection with no change to usual collection technique.

The pre-assembled product, **VACUETTE®** VISIO PLUS needle with **VACUETTE®** QUICKSHIELD Safety Tube Holder, has the advantage of the following features:

- The holder comes with pre-attached needle, sterile and single-packed.
- Rotational shield: Re-position the shield for the best possible position during blood collection.
- Transparent view window for optical venipuncture control.
- It is not necessary to thread in the needle.
- A practical, ready-to-use solution that saves time.
- Handling is optimised.

Safe disposal due to protected needle.

Your **Power** for Health


greiner bio-one



For further information, please visit our website www.gbo.com/preanalytics or contact us:

Austria (Main Office)

Greiner Bio-One GmbH
Phone +43 7583 6791-0
Fax +43 7583 6318
E-Mail office@at.gbo.com

Germany

Greiner Bio-One GmbH/Preanalytics
Phone +49 7022 948-0
Fax +49 7022 948-514
E-Mail office@de.gbo.com

Netherlands

Greiner Bio-One B.V.
Phone +31 172 4209 00
Fax +31 172 4438 01
E-Mail info@nl.gbo.com

United Kingdom

Greiner Bio-One Ltd.
Phone +44 1453 8252 55
Fax +44 1453 8262 66
E-Mail info@uk.gbo.com

Brazil

Greiner Bio-One Brasil
Phone +55 19 3468-9600
Fax +55 19 3468-9621
E-Mail office@br.gbo.com

Hungary

Greiner Bio-One Hungary Kft.
Phone +36 96 213 088
Fax +36 96 213 198
E-Mail office@hu.gbo.com

Spain

VACUETTE Espana S.A.
Phone +34 91 652 77 07
Fax +34 91 652 33 35
E-Mail info@vacuette.es

USA

Greiner Bio-One North America Inc.
Phone +1 704 261-7800
Fax +1 704 261-7899
E-Mail info@us.gbo.com

China

Greiner Bio-One Suns Co., Ltd.
Phone +86 10 83 55 19 91
Fax +86 10 63 56 69 00
E-Mail office@cn.gbo.com

India

Greiner Bio-One INDIA Pvt., Ltd.
Phone +91 120 456 8787
Fax +91 120 456 8788
E-Mail info@gboindia.com

Switzerland

Greiner Bio-One VACUETTE
Schweiz GmbH
Phone +41 7 12 28 55 22
Fax +41 7 12 28 55 21
E-Mail office@ch.gbo.com

France

Greiner Bio-One SAS
Phone +33 1 69 86 25 25
Fax +33 1 69 86 25 35
E-Mail office@fr.gbo.com

Japan

Greiner Bio-One Co. Ltd.
Phone +81 3350 58875
Fax +81 3505 8974
E-Mail info@jp.gbo.com

Thailand

Greiner Bio-One Thailand Ltd
Phone +66 38 4656 33
Fax +66 38 4656 36
E-Mail office@th.gbo.com



980184 rev.02, 02.2016 e