


Trace Elements Formulation

Cat. No.	25-021	25-022	25-023
Description	(A) 1000x	(B) 1000x	(C) 1000x
Units	mg/L	mg/L	mg/L

Components			
CuSO ₄ • 5H ₂ O	1.60	--	--
ZnSO ₄ • 7H ₂ O	863.00	--	--
Selenite • 2Na	17.30	--	--
Ferric citrate	1155.10	--	--
MnSO ₄ • H ₂ O	--	0.17	--
Na ₂ SiO ₃ • 9H ₂ O	--	140.00	--
Molybdic acid, ammonium salt	--	1.24	--
NH ₄ VO ₃	--	0.65	--
NiSO ₄ • 6H ₂ O	--	0.13	--
SnCl ₂ (anhydrous)	--	0.12	--
AlCl ₃ • 6H ₂ O	--	--	1.20
AgNO ₃	--	--	0.17
Ba(C ₂ H ₃ O ₂) ₂	--	--	2.55
KBr	--	--	0.12
CdCl ₂	--	--	2.28
CoCl ₂ • 6H ₂ O	--	--	2.38
CrCl ₃ (anhydrous)	--	--	0.32
NaF	--	--	4.20
GeO ₂	--	--	0.53
KI	--	--	0.17
RbCl	--	--	1.21
ZrOCl ₂ • 8H ₂ O	--	--	3.22



www.corning.com/lifesciences/solutions

At Corning, cells are in our culture. In our continuous efforts to improve efficiencies and develop new tools and technologies for life science researchers, we have scientists working in Corning R&D labs across the globe, doing what you do every day. From seeding starter cultures to expanding cells for assays, our technical experts understand your challenges and your increased need for more reliable cells and cellular material.

It is this expertise, plus a 160-year history of Corning innovation and manufacturing excellence, that puts us in a unique position to offer a beginning-to-end portfolio of high-quality, reliable cell culture consumables.

For additional product or technical information, please visit www.cellgro.com or call 1.800.235.5476.

Mediatech, Inc.
A Corning Subsidiary
9345 Discovery Boulevard
Manassas, VA 20109
t 800.235.5476
t 703.471.5955
f 703.467.9851
www.cellgro.com

CORNING

FALCON

AXYGEN

GOSSELIN

PYREX

Corning cellgro products are manufactured by Mediatech, Inc.
For a listing of Corning trademarks, visit us at www.cellgro.com/about-us/trademarks.html
All other trademarks in this document are the property of their respective owners.