


Overall dimensions




Please note images are for representation purposes only

Packed weight: 36.3kg
Product weight: 29.1kg

Packed dimensions:
A: 465mm
B: 610mm
C: 755mm

Product dimensions:
A: 245mm
B: 430mm
C: 640mm

Working area



Working area:
D: 118mm
E: 154mm
F: 145mm
(10mm below top)

LT ecocool™ 100 minimum immersion volume depth: 85mm

LT ecocool™ 150 minimum immersion volume depth: 130mm

Please note images are for representation purposes only

Vessel capacity without lid

Duran 21-801: 50ml	9 ⁺
Duran 21-801: 100ml	4
Duran 21-801: 150ml	3
Duran 21-801: 250ml	2
Duran 21-801: 500ml	1
Duran 21-801: 750ml	1
Duran 21-801: 1000ml	1
Duran YOUTILITY: 125ml	4
Duran YOUTILITY: 250ml	2
Duran YOUTILITY: 500ml	2
Duran YOUTILITY: 1000ml	1

Vessel capacity with tank lid

Erlenmeyer flask: 100ml	4
Erlenmeyer flask: 250ml	2
Erlenmeyer flask: 500ml	1
Erlenmeyer flask: 1000ml	1*

* Bridge plate released

+ May need height support in the tank depending on fill level

