



designed for scientists



## HABITAT cell sw 0.5 single-wall, for 0.5 l

/// Data Sheet

The 0.5-liter vessel package for the bioreactor HABITAT research contains a 0.5-liter single-walled glass vessel for cell culture applications as well as the appropriate stirring drive, including the following sensors: temperature-, pH-, DO-, level- and foam sensors. Additionally, the necessary connection accessories are included. Together with the separately available control tower package (e.g. HABITAT cell) you receive all necessary components for a successful cultivation.

Scope of delivery

- Stand 0.5
- HA.gv.sw.0.5 Glass vessel, single-wall, for 0.5 l



designed for scientists

- HA.mt.s.2 Motor, small, for 0.5, 1 and 2 l
- HA.hp.s.0.5 Harvest pipe, straight, for 0.5 l
- HA.fd Feeding port (2 pcs.)
- HA.sp.m.0.5 Micro sparger, for 0.5 l
- HA.cn Condenser
- HA.ip.pi.0.5 3-pitched blade impeller, for 0.5 l
- HA.s.tm.1 Temperature sensor, for 1 l
- HA.s.ph.0.5 Hamilton pH sensor, for 0.5 l
- HA.s.do.0.5 Hamilton DO analogue sensor, for 0.5 l
- HA.s.fo Foam sensor
- HA.s.lv.1 Fill level sensor, for 0.5 to 1 l
- HA.hbl.0.5 Heating blanket, for 0.5 l
- HA.cab.sw Cable & tube set, for single-wall vessels
- HA.sf.250 Sample flask, 250 ml (4 pcs.)

# IKA

designed for scientists