



designed for scientists



HABITAT cell sw 2 single-wall, for 2 l

/// Data Sheet

The 2-liter vessel package for the bioreactor HABITAT research contains a 2-liter single-walled glass vessel for cell culture applications as well as the appropriate stirring drive, including the following sensors: temperature-, pH-, DO-, level- and foam sensors. Additionally, the necessary connection accessories are included. Together with the separately available control tower package (e.g. HABITAT cell) you receive all necessary components for a successful cultivation.

Scope of delivery

- Stand 2
- HA.gv.sw.2 Glass vessel, single-wall, for 2 l



designed for scientists

- HA.mt.s.2 Motor, small, for 0.5, 1 and 2 l
- HA.hp.s.2 Harvest pipe, straight, for 2 l
- HA.ino Port for inoculation
- HA.sp.m.2 Micro sparger, for 2 l
- HA.cn Condenser
- HA.ip.pi.2 3-pitched blade impeller, for 2 l
- HA.s.tm.2 Temperature sensor, for 2 l
- HA.s.ph.2 Hamilton pH sensor, for 2 l
- HA.s.do.2 Hamilton DO analogue sensor, for 2 l
- HA.s.fo Foam sensor
- HA.s.lv.2 Fill level sensor, for 2 l
- HA.hbl.2 Heating blanket, for 2 l
- HA.cab.sw Cable & tube set, for single-wall vessels
- HA.sf.250 Sample flask, 250 ml (4 pcs.)

IKA

designed for scientists