

BAX® System

Q7 Workstation

With a BAX® System Q7 workstation, customers who own a 7500 Fast* Real-Time PCR System from Applied Biosystems can expand the functionality of their instrument by running BAX® System applications on the same instrument. The BAX® System computer is connected to a 7500 Fast instrument through a labeled data switch. When you wish to run BAX® System applications, simply set the switch to BAX; for 7500 Fast native applications, set the switch to AB.



BAX® System Q7 Workstation

Part # ASY2013

Includes computer with BAX® System software, monitor, A/B switch, USB cables, PCR tube racks and PCR tube holder trays

The BAX® System Q7 Workstation

The BAX® System Q7 workstation consists of a computer, monitor, A/B switch and peripherals that ship with BAX® System software already installed. Customers who already own both the AB 7500 Fast* and BAX® System instruments need only the Q7 workstation package to expand the functionality and throughput of their 7500 Fast instrument.

BAX® System Q7 Equipment Package

Part # ASY2014 (120V) Dry Block Heaters

Part # ASY2015 (220V) Dry Block Heaters

Part # ASY2089 (120V) (220V) Automated Thermal Block

Includes memory stick, heating and cooling blocks with inserts, capping/decapping tools, cluster tubes with caps and racks, adjustable mechanical pipettes, repeating pipette, multi-channel pipette, pipette tips, nitrile gloves and user documentation

The BAX® System Q7 Equipment Package

Customers who do not already own a BAX® System can purchase the BAX® System Q7 workstation plus the equipment package, which includes the additional equipment necessary for BAX® System assays, along with sufficient supplies for the first 96 tests.

Find Out More

For more information about how your business can benefit from a BAX® System Q7 workstation in your lab, contact your account manager or call **1-800-863-6842**.

* Due to the FDA certification status of the AB 7500 Fast Dx instrument, the Q7 workstation should not be used with this model.



Find support documents, instructional videos, and more at **www.hygiena.com**