



Budget series

Underbench Refrigerator Model: BIO115FRB

The **BioCOLD Budget** series of refrigeration has been built to meet the requirement for general laboratory storage. The Budget range does not have the options list available however it is still manufactured to the same high standards as all the Biocold products.



Standard Features	Options
<ul style="list-style-type: none"> Reversible doors Replaceable door gasket CFC and HCFC free refrigeration system 	N/A

Technical Specification			
Dimensions		Temperature	
External (H x W x D), (mm)	850 x 475 x 445	Set Point (adjustable) (°C)	+4
Weight (kg)	28	Operating range (°C)	+2 to +8
Capacity		Electrical Supply	
Refrigerator Volume (L)	95	Voltage (V)	220-240
Shelves (Glass)	2	Frequency (Hz)	50
Defrost Type	Manual	Fuse (A)	13



www.biocold.co.uk