



**Budget series**

## Underbench Refrigerator Model: BIO125FRB

The **BioCOLD Budget** series of refrigeration has been built to meet the requirement for general laboratory storage. The Budget range does not have the options list available however it is still manufactured to the same high standards as all the Biocold products.



Standard Features	Options
<ul style="list-style-type: none"> <li>Reversible doors</li> <li>Replaceable door gasket</li> <li>CFC and HCFC free refrigeration system</li> </ul>	N/A

Technical Specification			
<b>Dimensions</b>		<b>Temperature</b>	
External (H x W x D), (mm)	847 x 546 x 556	Set Point (adjustable) (° C)	+4
Weight (kg)	30	Operating range (° C)	+2 to +8
<b>Capacity</b>		<b>Electrical Supply</b>	
Refrigerator Volume (L)	125	Voltage (V)	220-240
<b>Shelves (Glass)</b>	2	Frequency (Hz)	50
<b>Defrost Type</b>	Manual	Fuse (A)	13



[www.biocold.co.uk](http://www.biocold.co.uk)