



- Extensive Range
- 90+ Years of History
- Extended Aftercare Support
- UKAS Laboratory

#### Vaccine Carriers and Bags

## Labcold Vaccine Carrier RPDF0012F

This Labcold portable vaccine carrier is designed to solve the problem of cold chain compliance on the move as it can be powered by both mains and vehicular 12 volt power sources.

There is an easy to read digital temperature display and because it is accurate to one decimal place, this Labcold vaccine carrier can be calibrated so you can be confident that the temperature displayed is as experienced by the contents of the carrier.

Using a thermoelectric, solid state cooling system there are no dangerous liquids that could escape or moving parts that could get damaged. The carrier comes complete with a convenient shoulder strap and the interior features a divider and cushioned lids to protect the contents in transit.

Purpose designed for any professional who needs to give vaccines in the field, this carrier can be used in conjunction with Labcold data loggers for extra temperature security.

Please be aware this product is not designed to be run for extended periods of time nor to hold temperature when not connected to power source.

Please be aware this product is not suitable for the transport of blood or blood products.

### Specifications

Product Code	RPDF0012F
Gross Capacity	12 litres
Style	Portable
UK Warranty	1 Year
External Dimensions (mm) HxWxD	295 x 430 x 295
Temperature Range	2°C to 10°C
Net Weight	4.6kg
Temperature Display	Yes



## Specifications (cont.)

Digital Temperature Control	Yes
Carry Strap	Yes
Moulded feet for stability in transit	Yes
Volt Operation	12 Volt DC and 240 AC Volt
Energy Consumption	1.61 kw/24hr
Ambient Operating Temperature	16°C to 32°C

## Features

- Temp range: +2° to +10°C
- Safe gas free peltier refrigeration
- Digital temperature display
- Tough plastic exterior
- Supplied with mains and in car power cables
- Lightweight and easily transportable
- Easy digital temperature adjustment
- Insulation pad

## Models

	Model	Gross Capacity
	Labcold Vaccine Carrier RPDF0019F	19 litres
	Labcold Vaccine Carrier RPDF0012F	12 litres