



PRODUCT SPECIFICATION

OVERSHOES ANTI-SLIP NON-STERILE

TuffGuard **6**

PRODUCT INFORMATION

PRODUCT DESCRIPTION

Disposable Overshoe

PACKAGING

Primary packaged in polythene packaging

SHELF LIFE

3 years

PRODUCT FEATURES

White
Anti-slip
Material - Polypropylene Laminate
Weight - 60gsm
Length - 400mm/tolerance (+/-5%)
Width - 180mm/tolerance (+/-5%)
Stretched elastic opening - 320mm/
tolerance (+/-5%)

QUALITY ASSURANCE

Product supplied under BS EN ISO 9001
registered systems

STORAGE CONDITIONS

Store in a dry, well-ventilated area in
ambient conditions

ORDERING/PACKING INFORMATION

Non-Sterile Code	Sterile Code	Description	Size	Pack Qty	Case Qty
SNWH03016		Overshoe White Anti-Slip Non-Sterile	180 mm x 400 mm	100	10 Packs of 100