



## SLS Lab Pro Eco Drying Cabinets

### Technical Specifications

	SLS5028	SLS5026	SLS5024	SLS5022
<b>Fan Assisted</b>	Yes	Yes	No	No
<b>External Dimensions (W x D x H mm)</b>	1168 x 650 x 1870	635 x 650 x 1870	1000 x 550 x 807	765 x 450 x 732
<b>Internal Dimensions (W x D x H mm)</b>	1030 x 590 x 1470	495 x 590 x 1470	860 x 480 x 525	625 x 375 x 450
<b>Capacity</b>	885L	426L	227L	113L
<b>Shelves</b>	4	4	3	3
<b>Power Rating</b>	2500 watts	1500 watts	750 watts	500 watts



SLS5028

SLS5026

SLS5024

SLS5022