



## Product Specification Sheet

### **TS/15-A-NACL95**

Sterile blue sponge dosed with 0.9% Sodium Chloride (NaCl) in 100 ml plastic container, includes glove

Product Code	TS/15-A-NACL95
Description	Sterile blue sponge dosed with 0.9% Sodium Chloride (NaCl) in 100 ml plastic container, includes glove
Intended Use	Sampling of meat carcasses for microbiological analysis
Components	Blue Cellulose Sponge, Clear plastic container, white screw lid, Plastic glove.
Contents	De-ionised Water – 1 Litre, Sodium Chloride –9g.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	12 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Some Fastidious organisms may struggle/fail to recover without extra supplement from the carcass/external sources.
Pack Size	95 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility. ii) Organism recovery (S.typhimurium)

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 2 Date: 24.04.2020  
ProductSpec\_TS15-A-NACL95