



Product Specification Sheet

TS/15-B:NaCl

Sterile blue sponge in easy open stomacher pouch, individually packed & dosed with 0.9% Sodium Chloride (NaCl)

Product Code	TS/15-B:NaCl
Description	Sterile blue sponge in easy open stomacher pouch, individually packed & dosed with 0.9% Sodium Chloride (NaCl)
Intended Use	Sampling of meat carcasses for microbiological analysis
Components	Blue Cellulose Sponge, easy open stomacher pouch.
Contents	De-ionised Water - 1 Litre, NaCl - 9g.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	14 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Some Fastidious organisms may struggle/fail to recover without extra supplement from the carcass/external sources.
Pack Size	200 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility. ii) Organism recovery (<i>S.typhimurium</i>)

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 4 Date: 03.09.2024

ProductSpec_TS15-B:NaCl