



Product Specification Sheet

TS/15-BPW

Sterile blue sponge in easy open stomacher pouch, individually packed & dosed with Buffered Peptone Water

Product Code	TS/15-BPW
Description	Sterile blue sponge in easy open stomacher pouch, individually packed & dosed with Buffered Peptone Water
Intended Use	Sponge for pathogen monitoring in poultry houses
Components	Blue Cellulose Sponge, easy open stomacher pouch.
Contents	De-ionised Water, Bacteriological Peptone, Phosphate Buffers.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	14 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Buffered Peptone water (BPW) is a generalised media hence will recover all organisms. Further selective steps are required
Pack Size	200 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility. ii) Organism recovery (<i>S.typhimurium</i>)

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 5 Date: 03.09.2024
ProductSpec_TS15-BPW