



Product Specification Sheet

TS/15-B

Sterile blue sponge in easy open stomacher pouch, individually packed & dosed with 10ml of Neutralising Buffer

Product Code	TS/15-B
Description	Sterile blue sponge in easy open stomacher pouch, individually packed & dosed with 10ml of Neutralising Buffer
Intended Use	Sampling for environmental and hygiene monitoring of surfaces and cleaning validation
Components	Blue Cellulose Sponge, easy open stomacher pouch.
Contents	De-ionised Water, Tween 80, Lecithin, Sodium Thiosulphate, Sodium Thioglycollate, Buffers to maintain pH - 7.2 ± 0.2
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	14 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	1. Neutralisation is limited to surface residues. 2. Due to product characteristics a small precipitate may be observed which does not affect performance
Pack Size	200 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility. ii) Organism recovery

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 11 Date: 03.09.2024

ProductSpec_TS15-B