



Product Specification Sheet

TS/15-CP

12 inch Foam tipped swab with a blue flexible polypropylene shaft, in polypropylene outertube.

Product Code	TS/15-CP
Description	12 inch Foam tipped swab with a blue flexible polypropylene shaft, in polypropylene outertube.
Intended Use	Environmental and drain sampling.
Components	Large blue plastic cap, 30cm (12") blue flexible plastic shaft, blue Foam Sponge Head, 30cm (12") polypropylene tube, Plug.
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Do not use in solvents
Pack Size	100 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 4 Date: 24.04.2020

ProductSpec_TS15-CP