
Enviroscreen

EnviroStik (Sponge with Handle)



- All in one **sampling solution**
- Stomachable & **Biocide Free** sponge
- **Quick,** touch free consistent sampling
- Peel open - **no scissors** needed
- Resealable **secure sample transport**
- **Easy Break Handle** for lab processing
- Reliable **collection & release**
- **Neutralise** disinfectant residues
- Vital for **HACCP** avoiding costly Recalls
- **Uncover** hidden contamination
- Strengthens **brand protection**

Ready to use & pre-moistened for convenience.
Easy open & re-sealable Stomachable pouch.
Ideal for environmental & hygiene monitoring of awkward areas, in/around equipment or critical surfaces.

ISO 9001:2015 Compliant
ISO 18593:2018 Compliant

12 Months Shelf Life
Sterilised by
Gamma Irradiation

Technical Service Consultants Ltd

T : +44 (0)1706 620600
F : +44 (0)1706 620445
E : sales@tscswabs.co.uk
W : www.tscswabs.co.uk

EnviroStik Sponge

EnviroStik offers end-users an ideal monitoring solution that is a convenient & reliable way to confirm the presence of living microorganisms on any critical surfaces, making it ideal for general hygiene sampling & environmental monitoring. The new design increases useability allowing concealed or awkward spaces to be easily monitored beyond just the top surface with improved control and precision. Touching or contaminating the sample is virtually impossible with EnviroStik as the breakpoint handle allows its quick removal eliminating any accidental contamination keeping counts accurate.

In many cases a surface will contain traces of disinfectants and/or sanitisers. Left alone and over time, these cleaning fluid traces will not kill the microorganisms, but will just prevent any growth (bacteriostatic effect).

TSC Neutralising Buffer will neutralise most common disinfectants including peroxides and derivatives, phenols, quaternary ammonium compounds (QAC's), amphoterics, chlorines and associated halogens, providing end-users with total testing assurance.

Product Features


- New non-rivet design allows trouble free stomaching
- Pre-moist sponge matrix offers enhanced sample collection & release
- Re-sealable stomacher pouch for sample integrity & easy processing
- Large white writing panel to record sample details
- Eliminates multiple sample transfers
- Reduces errors & contamination risk
- True sample counts for accurate results & assured product release




Download Certificates of Conformity from **www.tscswabs.co.uk**

Sampling Areas & Applications

- Ideal for sampling larger surfaces in a wide range of environments including food production areas, clinical & veterinary infection control and pharmaceutical & cosmetics facilities
- Recommended for food surfaces where disinfectants & sanitisers are used, production equipment, food contact surfaces & surrounding areas or use in cleaning validation
- Exposes potential microorganism contamination by removing chemical residues whilst preserving the sample for analysis
- Perfect for transporting samples for external processing
- Protecting your brand and preventing costly recalls





Code	Description	Qty	Recommended Use			
			Environmental	Food & Beverage	Clinical & Veterinary Infection Control	Pharmaceutical & Cosmetics
TS/15-SH (Straight handle)	> Resealable stomacher for easy sample processing > High visibility Blue Sponge > Pre-moistened with Neutralising Buffer > Ideal for Surface Sampling	125	✓	✓	✓	✓
TS/15-SA (Saline)	> Resealable stomacher for easy sample processing > High visibility Blue Sponge > Pre-moistened with Saline > Ideal for direct food contact or carcasses	125	✓	✓	✓	✓
TS/15-T40	> Sterile template > Allows easy compliance to recommended Standards for Sampling (eg ISO)	40	✓	✓	✓	✓

Don't forget - Technical Service Consultants Ltd cover all your swabbing needs - we have an extensive range of sponge swabs at www.tscswabs.co.uk

Technical Service Consultants Ltd

Technical Service Consultants Ltd
Microbiology House, Fir Street, Heywood, Lancashire OL10 1NW, UK
T : +44 (0)1706 620600 F : +44 (0)1706 620445
E : sales@tscswabs.co.uk W : www.tscswabs.co.uk