Technical Service Consultants Ltd
Microbiology House, Fir Street, Heywood, Lancashire OL10 1NW, UK
T : +44 (0)1706 620600 F : +44 (0)1706 620445
E : sales@tscswabs.co.uk W : www.tscswabs.co.uk

Product Specification Sheet

TS/15-BPS
EnviroStik - Sponge with Straight Handle dosed in Buffered peptone Water in easy open stomacher pouch

Product Code	TS/15-BPS
Description	EnviroStik - Sponge with Straight Handle dosed in Buffered peptone Water in easy open stomacher pouch
Intended Use	Sponge on Stick for environmental and hygiene monitoring including general pathogen testing.
Components	Blue Cellulose Sponge, Blue plastic stick with built-in breakpoint, easy open stomacher pouch.
Contents	De-ionised Water, Bacteriological Peptone, Phosphate Buffers.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	14 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Buffered Peptone water (BPW) is a generalised media hence will recover all organisms therefore cleaning after use is essential. Further selective steps are required.
Pack Size	125 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility ii) Product Specification iii) Organism recovery

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 0 Date: 22.10.2024
ProductSpec_TS15-BPS